



Aug 22, 2021

Microfluidics 3 - Mold Fabrication: 3D Printing

Serhat Sevlı¹, C. Yunus Sahan¹¹Nehir Biyoteknoloji Ltd. www.nehirbt.com

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bxj2pkqe

NehirBT

Tech. support email: bilgi@nehirbt.com.tr

Serhat S

ABSTRACT

Su8 is a photoresist resin used in MEMS technology for thick structures. PDMS microfluidic chips fabrication is suitable for SU8+Si wafer molds. This protocol describes the Direct UV laser lithography application performed in our laboratories.

DOI

dx.doi.org/10.17504/protocols.io.bxj2pkqe

PROTOCOL CITATION

Serhat Sevlı, C. Yunus Sahan 2021. Microfluidics 3 - Mold Fabrication: 3D Printing. **protocols.io**<https://dx.doi.org/10.17504/protocols.io.bxj2pkqe>

KEYWORDS

Microfluidics, Spin coating, SU8, Si wafer

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 20, 2021

LAST MODIFIED

Aug 22, 2021

PROTOCOL INTEGER ID

52570

MATERIALS TEXT

3D Printing Resin; OOkuma, FP60 functional prototyping resin,
3D Printing Resin; ResinWorks, Master Mold Resin for PDMS
Cleaning; Isopropanol alcohol and distilled water

SAFETY WARNINGS

* All the related steps must be done in a class-D cleanroom (minimum) and under yellow or red light.

* Please learn well about the hazards of 3D printing resin chemicals.

BEFORE STARTING

This protocol is derived from

3D Printer Setup 30m

- 1 Molds for PDMS microfluidic chips can be fabricated using 3D Printing technology. SLA(/DLP) printers are mostly preferred for microfluidics due to its high resolution, relative low cost, and being suitable for direct microchannel formation. ^{30m}

NehirBT uses a SLA/DLP 3D printing technology for transferring drawings by CAD software to channels directly as mold.

The important parameters arranged are;

A	B
Specs	Info
Total size	50x110x150
Light source	UV LED/DLP
Point size	30-50 microns
Velocity	30 sec per layer
Resolution	30-50 microns
Aspect ratio	1:100
3D structures	unlimited

3D SLA/DLP Printer
Creative CadWorks MiiCraft 100

Microfluidic Channel Design 5h

- 2 The designs are generated using a CAD software in STL format and with millimeter as a unit. ^{5h}

The flow simulation is performed using a multiphysics finite element analysis software, so the design optimization is performed by simulation support and correcting the errors.

MiiCraft software is used for converting the design file to MII format which is suitable for MIICRAFT100 instrument.

MiiCraft software controls the MIICRAFT100 instrument by arranging the duration of exposure per layer and layer thickness.

3D Printing 2d

- 3 The resin choice is critical for the 3D Printing of PDMS molds since chemicals inside resin materials inhibit the curation of PDMS polymer. Only there are a few suitable resins and NehirBT uses; ^{5h}
1- OOkuma; FP60 functional prototyping resin,
2- ResinWorks; Master Mold Resin for PDMS

Resin is put inside the reservoir of the 3D printing device.

The printing is initiated by the setup of

A	B
Specs	Info
Layer thickness	30-50 microns
Layer curing time	1-3 sec
Gap adjustment	0,1 mm
Base layer curing time	1-20 sec
Base layer layer number	1-10

Post Fabrication Recipe 10m

3d

- 4 After fabrication of the PDMS mold, it must be thermally and optically cured.
Optical curing is exposure to UV LED (385-405 nm) for several hours and thermal treatment at 95°C for several hours is needed.



Not enough post-fabrication treatment causes remaining inhibitory chemicals which causes non-cured PDMS.

Hot Plate
Electromag LB.EM.M4060

UV LED Illuminator
custom made 385-405nm

Washing of 3D Printed Mold 30m

30m

- 5 The prepared mold is washed with isopropanol after completing the treatments and then dried with pressurized nitrogen gas.

The mold is now ready for PDMS microchannel manufacturing. Put the mold in a petri dish, and label it with all appropriate information. Pouring PDMS material on mold in a petri dish is performed in class2 laminar flow hood.