



JUL 31, 2023

🌐 ITS1 Amplicon Prep

aglazer¹

¹TGen North



aglazer

ABSTRACT

This protocol was provided by IDT in the the xGen ITS1 Amplicon Panel manual.

OPEN  ACCESS



ATTACHMENTS

[xgen-16s-amplicon-panel-v2-and-xgen-its1-amplicon-panel-protocol.pdf](#)

Protocol Citation: aglazer
2023. ITS1 Amplicon Prep.
protocols.io
<https://protocols.io/view/its1-amplicon-prep-cxyfxptn>

MANUSCRIPT CITATION:
xGen 16S Amplicon Panel v2
xGen ITS1 Amplicon Panel.
(n.d.). Retrieved July 30,
2023, from
[https://sfvideo.blob.core.wind
ows.net/sitefinity/docs/default
t-source/protocol/xgen-16s-
amplicon-panel-v2-and-xgen-
its1-amplicon-panel-
protocol.pdf?](https://sfvideo.blob.core.windows.net/sitefinity/docs/default-source/protocol/xgen-16s-amplicon-panel-v2-and-xgen-its1-amplicon-panel-protocol.pdf?sfvrsn=acace007_16)
sfvrsn=acace007_16

License: This is an open
access protocol distributed
under the terms of
the [Creative Commons
Attribution License](#), which
permits unrestricted use,
distribution, and reproduction
in any medium, provided the
original author and source
are credited

Protocol status: In
development
We are still developing and
optimizing this protocol

Created: Jul 31, 2023


Last Modified: Jul 31, 2023

PROTOCOL integer ID:
85735

Prepare panel-specific Multiplex PCR Reaction Mix

- 1 Before mixing reagents, calculate the total volume of the Master Mix based on the number of reactions required, with 10% overage for pipetting.
- 2 Make the Multiplex PCR Reaction Mix according to the table below and **pipette-mix** the Master Mix thoroughly.

Note

Vortex components G1 (ITS1 Amplicon Panel Primer) and G2 and pulse-spin tubes to collect contents before use. Gently shake Enzyme G3 at room temperature for  00:05:00 minutes or until all solutes appear to be in solution. Place back on ice for the remainder of the protocol.

Note

Keep prepared Master Mix and all components **ON ICE** until ready to use.

Reagents	Volume per sample (uL)
Reagent G1 - ITS1 Amplicon Panel Primer	2
Reagent G2	3
Enzyme G3	15
Total Volume	20

Multiplex PCR Master Mix Recipe

- 3 Use a multichannel to pipette out  20 µL the PCR Mastermix into each well of a mid-skirt plate as needed and keep ON ICE until template addition.

Perform Multiplex PCR



- 4 Pre-program the thermal cycler for multiplex PCR according to the table below and **allow the block to reach 98°C before loading samples**

Note

Confirm that lid heating is turned **ON** and is set to 105°C

A	B	C
	98C	30 sec
	98C	10 sec
4 cycles	63C	5 min
	65C	1 min
14 cycles	98C	10 sec
	64C	1 min
	65C	1 min
	4C	Hold

Multiplex PCR Thermal Cycler Program

- 5 Use a multichannel to add  10 µL of template to the  20 µL Multiplex PCR Reaction Mix and pipette-mix thoroughly.
- 6 Place the mid-skirt plate containing template and Master Mix in the thermal cycler and run the


Note

Near the completion of the thermal cycler run, prepare the Indexing Reaction Mix in the post-PCR area. **Assemble this reaction mix on ice and keep cold until adding it to samples in the indexing step.** All components except indexes may be added to the Master Mix when running multiple samples in parallel.

Perform Post-Multiplex PCR Cleanup






10m 30s

- 7 Make sure that beads and samples are at room temperature. Spin down the sample plate and briefly vortex beads to homogenize before use.

- 8 Add  30 μL of magnetic beads (ratio: 1.0x) to each 30 μL sample after thermal cycling and mix by **vortexing**. Pulse-spin the samples in a microcentrifuge tube to collect contents.


Note

Make sure no bead-sample suspension droplets are left on the sides of the tube. If droplets are visible, pulse-spin the samples in a microcentrifuge to collect contents.


- 9 Incubate the samples for  00:05:00 minutes at room temperature, off the magnet. After incubation, place the sample tubes on a magnetic rack until the solution clears and a pellet is formed (about  00:05:00 minutes).
- 10 While the sample is on the magnet, remove and discard the supernatant without disturbing the pellet (about  5 μL may be left behind) using a clean pipette tip. Leave the plate on the magnet.
- 11 Add  180 μL of **freshly prepared** 80% ethanol solution to the pellet while it is still on the magnet. Use care not to disturb the pellet. Incubate for  00:00:30 seconds, then carefully remove the ethanol solution.

10m

30s

- 12 Repeat the previous step, for a second wash with the ethanol solution.
- 13 Pulse-spin the samples in a microcentrifuge, place back onto the magnet, and remove any residual ethanol solution from the bottom of the tube with a small volume tip.
- 14 Resuspend each bead pellet in  17.4 µL post-PCR TE Buffer.

Note

If an "off-bead" PCR is preferred, place tubes back on the magnet and transfer  17.4 µL of eluate to a fresh plate.

Note

Keep samples at room temperature. At no time should "with bead" samples be stored on ice, as this could affect binding to magnetic beads.

Prepare Indexing PCR Master Mix

- 15 Before mixing reagents, calculate the total volume of master mix based on the number of reactions required, with 10% overage for pipetting. It is recommended to prepare at least 10 reactions at any one time to maintain a volume of Enzyme I3 that can be accurately pipetted.
- 16 Prepare the Indexing Master Mix **ON ICE** according to the following table

Note

All components **except indexes** may be added to the Master Mix when running multiple samples in parallel.

Component	Volume per Sample (uL)
-----------	------------------------

Component	Volume per Sample (uL)
Reagent I1	3.3
Enzyme I2	0.5
Enzyme I3	0.1
Enzyme I4	25
Total Volume	28.9

Indexing Master Mix Recipe

- 17 Set aside the Indexing MMX **on ice** until template addition.

Perform Indexing PCR

- 18 Program the thermal cycler with the Indexing PCR program described in the table below and **allow the block to reach 37°C before loading samples.**

Note





Confirm lid heating is turned **ON** and set to reach 105°C.

Note

The PCR cycle number can be increased for samples that may give low yields.


A	B	C
	37C	20 min
	98C	30 sec
7 cycles	98C	10 sec
	60C	30 sec
	66C	1 min
	4C	Hold

Indexing PCR Program

- 19 Add  3.7 μL of the pre-mixed xGen Amplicon UDI primer pair to each sample ( 17.4 μL).
- 20 Add  28.9 μL of the cold indexing PCR reaction mix, stored on ice, to each sample. Mix thoroughly by **pipetting** (total volume will be  50 μL).
- 21 Spin down the sample plate, place in the thermal cycler, and run the Indexing PCR program.

Perform Post-Indexing PCR Cleanup

12m 30s

- 22 Make sure the PEG NaCl solution is at room temperature. Spin down the sample plate and briefly vortex the PEG NaCl solution to homogenize before use.
- 23 Add  42.5 μL of PEG NaCl solution (ratio: 0.85x) to each 50 μL sample after thermal cycling and mix by **vortexing**.

Note


Make sure no bead-sample suspension droplets are left on the sides of the tube. If droplets are visible, pulse-spin the samples in a microcentrifuge to collect contents.


Note



If performing an "off-bead PCR," use 42.5 μL (ratio: 0.85x) of fresh magnetic beads.

- 24 Incubate the samples for  00:05:00 minutes at room temperature, off the magnet. Pulse-spin the samples in a microcentrifuge. Afterwards, place the sample tubes on a magnetic rack

10m



until the solution clears and a pellet is formed (about  00:05:00 minutes).


25 While the sample is on the magnet, remove and discard the supernatant without disturbing the pellet (about  5 μL may be left behind) using a clean pipette tip. Leave the plate on the magnet.

26 Add  180 μL of **freshly prepared** 80% ethanol solution to the pellet while it is still on the magnet. Use care not to disturb the pellet. Incubate for  00:00:30 seconds, then carefully remove the ethanol solution. 30s

27 Repeat the previous step, for a second wash with the ethanol solution.

28 Pulse-spin the samples in a microcentrifuge, place back onto the magnet, and remove any residual ethanol solution from the bottom of the tube with a small volume tip.

29 Immediately add  20 μL of post-PCR TE buffer and resuspend the pellet by pipetting up and down until homogenous. Incubate for  00:02:00 minutes off the magnet. 2m

30 Place the samples back on the magnet and transfer the clean  20 μL library eluates to a fresh plate. Make sure that the eluate does not contain magnetic beads (indicated by brown color in eluate). If magnetic beads are present, place the plate on the magnet, wait for pellets to re-form, then transfer eluate again.

31 **SAFE STOP:** Store freshly prepared libraries at -20°C .