



Aug 13, 2020

# Enzyme linked immunosorbent assay for investigating the binding of protein-LG (SpLG) to immunoglobulins.

Angel A Justiz-Vaillant<sup>1</sup>, Norma McFarlane-Anderson<sup>2</sup><sup>1</sup>University of the West Indies St. Augustine; <sup>2</sup>University of West Indies. Mona Campus**1** Works for me [dx.doi.org/10.17504/protocols.io.bjpmk6](https://dx.doi.org/10.17504/protocols.io.bjpmk6)[University of the West Indies](#) [angel.vaillant@sta.uwi.edu](mailto:angel.vaillant@sta.uwi.edu)Angel Justiz-Vaillant  
University of the West Indies St. Augustine

## ABSTRACT

This SpLG ELISA can be used to detect specific antibodies in various animal species including human, mouse, rat, dog, rabbit, chicken, monkey, pig and hamster [1].

1. Kihlberg BM, Sjöbring U, Kastern W, Björck L. Protein LG: a hybrid molecule with unique immunoglobulin binding properties. *J Biol Chem.* 1992;267(35):25583-25588.

## DOI

[dx.doi.org/10.17504/protocols.io.bjpmk6](https://dx.doi.org/10.17504/protocols.io.bjpmk6)

## PROTOCOL CITATION

Angel A Justiz-Vaillant, Norma McFarlane-Anderson 2020. Enzyme linked immunosorbent assay for investigating the binding of protein-LG (SpLG) to immunoglobulins.. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjpmk6>

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 13, 2020

## LAST MODIFIED

Aug 13, 2020

## PROTOCOL INTEGER ID

40397

- 1 This ELISA is used to study the interaction of recombinant protein LG (SpLG) with different immunoglobulin preparations.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 µg/µl per well of SpLG in carbonate-bicarbonate buffer pH 9.6.
- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.

- 4 50 µl of animal serum (1 mg/ml) is added and incubated for 1 h at room temperature and the microplate is rewashed 4X with PBS-Tween.
- 5 Then 50 µl of peroxidase-labeled SpLG conjugate diluted 1:5000 in PBS-non-fat milk is added to each well and incubated for 1 h at RT. The plate is washed 4X with PBS-Tween.
- 6 50 µl of 4 mg/ml o-phenylenediamine solution (OPD) is added and the plate is incubated 15 minutes at RT in the dark.
- 7 The reaction is stopped with 50 µl of 3M H<sub>2</sub>SO<sub>4</sub> solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 492 nm.
- 9 A cut-off point should be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher the affinity of SpLG to immunoglobulin G.