



Feb 25, 2021

WU sc-prep Protocol for Solid Tumors v2.1

Reyka Jayasinghe¹, Li Ding¹, Feng Chen¹¹Washington University, Saint Louis

1

Works for me

dx.doi.org/10.17504/protocols.io.bsnqnddw

Reyka Jayasinghe

SUBMIT TO PLOS ONE

ABSTRACT

Single-cell Isolation for solid tumors Ding Lab 2021

DOI

dx.doi.org/10.17504/protocols.io.bsnqnddw

PROTOCOL CITATION

Reyka Jayasinghe, Li Ding, Feng Chen 2021. WU sc-prep Protocol for Solid Tumors v2.1. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bsnqnddw>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 22, 2021

LAST MODIFIED

Feb 25, 2021

PROTOCOL INTEGER ID

47536

Reagents

- 1 -[Human tumor dissociation kit](#) (Miltenyi Biotec #130-095-929)
 - PBS
 - PBS, 0.5% BSA (keep at 4°C)
 - (optional) Red Blood Cell (RBC) Lysis Solution ACK buffer Fisher (A1049201)
 - Falcon 15mL/50mL Conical Centrifuge Tubes (Corning Life Sciences DL, Fisher Scientific)
 - Large filter: Fisherbrand Cell Strainers. Nylon mesh 40um (Fischer Scientific, Cat # 22-363-547)
 - Smaller filter: PluriSelect Mini-Strainer 40 uM (Cat #43-10040-60)
 - Petri dishes
 - LoBind Protein or Genomic Microcentrifuge Tubes, Eppendorf®, Volume=2 ml Fisher 13-698-792
 - All steps should be performed under sterile conditions
 - For tissue fragments >1000mg, scale up reagent volumes (volumes listed are for 1 g tissue)

Octo Dissociation

- 2 Mince tissue into small pieces with razor blade in petri dish.

- 3 Resuspend tissue in enzyme mix (2.5 mL DMEM to 100 uL of Enzyme H, 50 uL Enzyme R, and 12.5 uL enzyme A) in C-Tube
- 4 Incubate on gentleMACS™ Octo Dissociator - 37h_TDK_1 setting (incubation can be shortened if tissue is very soft ~20-25 min). Check sample every 10-15 minutes.
- 5 After dissociation, filter cells through 40um filter and spin down in 15ml tube at 450g for 5 min at 4 degrees
- 6 If there is any red visible in the cell pellet proceed to ACK RBC lysis step below (Refer to bottom of protocol for examples).
RBC ACK Lysis Solution
-Slowly resuspend cell pellet in 500 uL-5mL (5ml for larger pellets) of RBC ACK Lysis Solution
-Incubate on ice for 30 sec to 5 min (longer incubation for more RBCs)
-Quench reaction with 10mL of BSA/PBS buffer
-Centrifuge at 450g 5min at 4 degrees
- 7 Do 1-3 1mL washes of 0.5% BSA in PBS on ice (spin at 450g 5min at 4 degrees) to remove dead cells and debris. Do washes in 2mL tube. After each wash aspirate off BSA/PBS and leave ~50-100 uL of BSA/PBS buffer. Resuspend cells using 1000 uL wide-bore pipette tip and count 10 uL of sample using a hemacytometer and trypan blue.

After first wash take 10 uL of sample and visualize under microscope (for smaller amounts take 2 uL of sample and 8 uL of wash buffer). If after 1 wash there are too many dead cells – consider using dead cell removal kit or debris removal kit:

- Dead Cell Removal: <https://www.miltenyibiotec.com/US-en/products/macs-sample-preparation/removal-reagents/dead-cell-removal-kit.html>
- Debris Removal: <https://www.miltenyibiotec.com/US-en/products/macs-sample-preparation/removal-reagents/debris-removal-solution.html>

Note: For Debris removal – use slow speed 450-500 g for all spinning steps.

If cell suspension looks pretty clean after first wash, then second wash should be clean to load.

If you have approximately 400-2000 cells/uL then continue on to GEMs if doing single cell. Load 20,000 cells.