

Jun 30, 2020

$^{14}\text{CO}_2$ -based assay for measuring Rubisco activity & activation state

In 1 collection


Cristina Rodrigues Gabriel Sales¹, Anabela Silva², Elizabete Carmo-Silva¹¹Lancaster Environment Centre, Lancaster University, Library Avenue, Lancaster, LA1 4YQ, UK;²Biosystems & Integrative Sciences Institute (BioISI), Science Faculty of Lisbon University, Lisbon, 1749-016, Portugal

1 Works for me dx.doi.org/10.17504/protocols.io.bf8cjrsw

Cristina Rodrigues Gabriel Sales

ABSTRACT

The Rubisco activity $^{14}\text{CO}_2$ -based assay measures the incorporation of $^{14}\text{CO}_2$ into the acid-stable product 3-phosphoglycerate (3-PGA). This protocol is based on Parry et al. (1997).



THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Sales CRG, Silva AB, Carmo-Silva E. 2020. Measuring Rubisco activity: challenges and opportunities of NADH-linked microtiter plate-based and ^{14}C -based assays. Journal of Experimental Botany, <https://doi.org/10.1093/jxb/eraa289>

DOI

dx.doi.org/10.17504/protocols.io.bf8cjrsw

PROTOCOL CITATION

Cristina Rodrigues Gabriel Sales, Anabela Silva, Elizabete Carmo-Silva 2020. $^{14}\text{CO}_2$ -based assay for measuring Rubisco activity & activation state. **protocols.io**
[dx.doi.org/10.17504/protocols.io.bf8cjrsw](https://doi.org/10.17504/protocols.io.bf8cjrsw)

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Sales CRG, Silva AB, Carmo-Silva E. 2020. Measuring Rubisco activity: challenges and opportunities of NADH-linked microtiter plate-based and ^{14}C -based assays. Journal of Experimental Botany, <https://doi.org/10.1093/jxb/eraa289>

COLLECTIONS ⓘ



Protocols from Sales et al. (2020) Rubisco activity: challenges and opportunities of NADH-linked microtiter plate-based and ^{14}C -based assays



Protocols from Sales et al. (2020) Rubisco activity: challenges and opportunities of NADH-linked microtiter plate-based and ^{14}C -based assays

KEYWORDS

Enzyme activity assay, Rubisco, Radiometric assay, Crop improvement

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution](https://creativecommons.org/licenses/by/4.0/)

[License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 11, 2020

LAST MODIFIED

Jun 30, 2020

PROTOCOL INTEGER ID

36836

PARENT PROTOCOLS

Part of collection

[Protocols from Sales et al. \(2020\) Rubisco activity: challenges and opportunities of NADH-linked microtiter plate-based and ¹⁴C-based assays](#)

[Protocols from Sales et al. \(2020\) Rubisco activity: challenges and opportunities of NADH-linked microtiter plate-based and ¹⁴C-based assays](#)

GUIDELINES

1. Check the "Materials" tab for a list of all the chemicals used in this protocol.
2. In the "Steps" tab, there is a brief description of the materials and equipment necessary for the protocol execution.
3. In the "Steps" tab, there is information on preparation of solutions, procedures for determining Rubisco initial and total activities, and notes to take into consideration to ensure reliable results.
4. The references cited are at the end of the "Materials" tab.

MATERIALS

NAME	CATALOG #	VENDOR
Bicine	B3876	Sigma Aldrich
Magnesium chloride hexahydrate (MgCl₂·6H₂O)	M2393	Sigma Aldrich
Ethylenediaminetetraacetic acid disodium salt dihydrate (EDTA)	E1644	Sigma Aldrich
Benzamidine	B6506	Sigma Aldrich
ε-Aminocaproic acid	A2504	Sigma Aldrich
Sodium hydroxide (NaOH)	S5881	Sigma Aldrich
2-Mercaptoethanol	M6250	Sigma Aldrich
DL-Dithiothreitol (DTT)	43819	Sigma Aldrich
Phenylmethanesulfonyl fluoride (PMSF)	P7626	Sigma Aldrich
Protease inhibitor cocktail	P9599	Sigma Aldrich
D-Ribulose 1,5-bisphosphate sodium salt hydrate (RuBP)	83895	Sigma Aldrich
Sodium bicarbonate [¹⁴C] (NaH¹⁴CO₃)	NEC086H005MC	Perkin Elmer
Potassium hydroxide (KOH)	P5958	Sigma Aldrich
Formic acid	F0507	Honeywell
Gold star quanta scintillation cocktail	QSQ1	Meridian
Ethanol absolute 99.8 %	10437341	Fisher Scientific

MATERIALS TEXT



Carmo-Silva E, Andralojc PJ, Scales JC, Driever SM, Mead A, Lawson T, Raines CA, Parry MAJ (2017). Phenotyping of field-grown wheat in the UK highlights contribution of light response of photosynthesis and flag leaf longevity to grain yield. *Journal of Experimental Botany* 68: 3473-3486.
<https://doi.org/10.1093/jxb/erx169>



Kane HJ, Wilkin JM, Portis AR, Andrews TJ (1998). Potent inhibition of ribulose-bisphosphate carboxylase by an oxidized impurity in ribulose-1,5-bisphosphate. *Plant Physiology* 117: 1059-1069.
<http://10.1104/pp.117.3.1059>



Parry MAJ, Andralojc PJ, Parmar S, Keys AJ, Habash D, Paul MJ, Alred R, Quick WP, Servaites JC (1997). Regulation of Rubisco by inhibitors in the light. *Plant, Cell & Environment* 20: 528-534.
<https://doi.org/10.1046/j.1365-3040.1997.d01-85.x>



Sharwood RE, Sonawane BV, Ghannoum O, Whitney SM (2016). Improved analysis of C4 and C3 photosynthesis via refined in vitro assays of their carbon fixation biochemistry. *Journal of Experimental Botany* 67: 3137-3148.
<https://doi.org/10.1093/jxb/erw154>



Wong C-H (1980). Practical enzymatic syntheses of ribulose 1,5-bisphosphate and ribose 5-phosphate. *Journal of the American Chemical Society* 102: 7938-7939.
<https://doi.org/10.1021/ja00547a023>

SAFETY WARNINGS

Work with radiation should follow local safety procedures.

Before using the protocol always check the Safety Data Sheet (SDS) for each chemical.

BEFORE STARTING

MATERIAL & EQUIPMENTS (for list of chemicals check "Materials" tab)

- Leaf sample frozen in -80°C
- Centrifuge for microtubes (speed 14000 *g*, 4 °C; VWR, Mega Star 600R)
- Dry block heating system (Grant Instruments QBD4)
- Vortex
- Chronometer
- Fume hood
- Pipette set
- Mortar and pestle
- Glass vials for liquid scintillation counting (Perkin Elmer 6000167, 7 mL)
- 1.5 mL microtubes

REAGENTS & SOLUTIONS

1

REAGENTS & SOLUTIONS TO PREPARE BEFOREHAND

1.1 Basic extraction buffer (1x)

[M]50 Milimolar (mM) Bicine-NaOH pH8.2

[M]20 Milimolar (mM) MgCl₂·6H₂O

[M]1 Milimolar (mM) EDTA

[M]2 Milimolar (mM) Benzamidine

[M]5 Milimolar (mM) ε-Aminocaproic acid

- Dissolve in ultrapure H₂O; adjust pH to 8.2 with NaOH; degas the solution bubbling with nitrogen (5 min/100 mL), then add:

[M]50 Milimolar (mM) 2-Mercaptoethanol

- Adjust for the final volume; it can be dispensed in aliquots (e.g. 50 mL Falcon tubes).

♻ -20 °C (storage)

1.2 [M]1 Molarity (M) DTT

- Dissolve in ultrapure H₂O. ♻ 4 °C (storage)

1.3 [M]100 Milimolar (mM) PMSF

- Dissolve in ethanol 99%. ♻ 4 °C (storage)

1.4 Plant protease inhibitor cocktail

♻ -20 °C (storage)

1.5 Basic assay buffer (2x)

[M]200 Milimolar (mM) Bicine-NaOH

[M]40 Milimolar (mM) MgCl₂·6H₂O

- Dissolve in ultrapure H₂O; adjust pH to 8.2 with NaOH; adjust for the final volume; degas the solution bubbling with nitrogen (5 min/100 mL). It can be dispensed in aliquots (e.g. 50 mL Falcon tubes).

♻ -20 °C (storage)

1.6 [M]**30 Milimolar (mM)** RuBP

⌚ **-20 °C (storage)**



High purity RuBP ($\geq 99\%$) is required to avoid interference in measurable activity due to the presence of RuBP-analogs that inhibit carboxylation (Kane et al., 1998; Sharwood et al., 2016). It is available commercially or it can be produced enzymatically from AMP-5' monohydrate and ATP disodium salt (Wong, 1980).

1.7 [M]**0.1 Molarity (M)** $\text{NaH}^{14}\text{CO}_3$ (0.5 Ci/mol)

⌚ **-20 °C (storage)**

1.8 [M]**0.3 Molarity (M)** KOH

⌚ **Room temperature (storage)**

1.9 [M]**10 Molarity (M)** Formic acid

⌚ **Room temperature (storage)**

1.10 Gold Star Quanta scintillation cocktail

⌚ **Room temperature (storage)**

2 SOLUTIONS TO PREPARE JUST BEFORE USE

- Prepared with reagents/solutions described in step 1.

2.1 Complete extraction buffer

1x Basic extraction buffer (from step 1.1)

[M]**10 Milimolar (mM)** DTT (from step 1.2)

[M]**1 Milimolar (mM)** PMSF (from step 1.3)

[M]**1 % (v/v)** Plant protease inhibitor cocktail (from step 1.4)

- Prepare the volume considering the number of extractions to be performed throughout the day plus two extras (to have a little excess). Mix all together. ⌚ **On ice**



The volume of extraction buffer will depend on the size of the leaf sample and the protein content, therefore it is species dependent and should be tested beforehand. It is important to ensure that the Rubisco concentration in the assays does not compromise the sensitivity of the assays (e.g. too much would consume all the substrate too quickly).

2.2 Complete assay buffer (volume per vial or reaction)

■ **250 µl** 2x Basic assay buffer (from step 1.5) - final concentration in the solution = 1x

■ **50 µl** 100 mM NaH¹⁴CO₃ (from step 1.7) - final concentration in the solution = 10 mM

■ **165 µl** ultrapure H₂O

- Prepare the volume considering that the activities are measured in duplicates (2 technical replicates).
- If the goal is to assay Rubisco initial and total activities, for each sample extract, it is necessary to have the volume for 4 vials (2 initials + 2 totals).
- In addition, account for additional volume for at least two more vials for Blanks (background counts, in the absence of RuBP), plus one extra (to have a little excess).



Example for 10 extractions: 10x (2xInitial + 2xTotal activity assays) + 2xBlanks + extra = volume for 43 vials.



■ **2 Milimolar (mM)** KH₂PO₄ can be added in the complete assay buffer for field samples to maximise the measurable Rubisco activity (Carmo-Silva et al., 2017).

PROCEDURE

3 START

- Thaw the frozen solutions that will be used in the day.
- Turn on the heat block and set to the temperature to be used for Rubisco activity measurements.



The temperature to be used for the Rubisco activity measurement depends on the experiment goals. Typical measurement temperatures are **25 °C** (standard) and **30 °C**, depending on the species. Assays can be performed at a range of temperatures, however special care should be taken at high temperatures (e.g. **50 °C**), as this can lead to evaporation of the assay mix; rates will be faster, i.e. the assay might become less sensitive; and may cause difficulty in handling hot vials.

- Turn on the centrifuge and set to **4 °C**.
- Collect samples from **-80 °C** into liquid nitrogen.
- Prepare the complete extraction buffer (step 2.1) and the complete assay buffer (step 2.2) and keep it **On ice**

4 EXTRACTIONS & RUBISCO ASSAYS

- 4.1
- Immediately before starting the extraction, place 4 vials in the heat block (for 2 initial and 2 total activity assays).
 - Add **465 µl** of complete assay buffer (from step 2.2.) to each of the 4 vials.
 - Add **10 µl** of 30 mM RuBP (from step 1.6) to the 2 vials for initials.

4.2 Extraction

- Add the complete extraction buffer to an ice-cold mortar.
- Take a sample from the liquid nitrogen container and add to the mortar.
- Grind the sample thoroughly for 🕒00:00:30 to maximum of 🕒00:01:00 .
- Collect the homogenate into an ice-cold 1.5 mL microtube and centrifuge 🌀14000 x g, 4°C 00:01:00 .



To prevent Rubisco deactivation (or even denaturation) the extraction should not take more than 1 min and it should be done in a ice-cold mortar, keeping the sample cold at all times. In our hands, with the extraction buffer described (containing protease inhibitors, mercaptoethanol and DTT, which keeps the enzyme reduced) 1 min centrifugation does not impact Rubisco activity. However, this should be tested for each species and extraction buffer used.

- When centrifugation stops, take the extract supernatant into another ice-cold 1.5 mL microtube.
- Proceed with the Rubisco assays straight away.
- Add 🧴25 µl of sample extract consecutively to each of 4 vials, with 15 s intervals (Total1, Total2, Initial1, Initial2). All additions are completed within 1.5 min of finishing centrifugation.



The Initial activity assays start with extract addition, while the Total activity assays start with addition of RuBP after 3 min of extract incubation with CO₂ and Mg²⁺ to allow for Rubisco carbamylation.



This protocol can be adapted for measuring Rubisco activity with purified enzyme. In this case, Rubisco is frequently pre-activated and initial activity assays are performed.

- Quench the assays after 30s by adding 🧴100 µl 10 M formic acid (from step 1.9).
- A possible time-line for the assays would be:

Step	Solution to be added	Initial activity		Total activity	
		Rep 1	Rep 2	Rep 1	Rep 2
1	Add 465 µL assay buffer	Before extraction	Before extraction	Before extraction	Before extraction
2	Add 10 µL of 30 mM RuBP	Before extraction	Before extraction	NA	NA
Proceed to the extraction. Follow the next steps when the extract supernatant is ready					
Start the chronometer					
3	Add 25 µL of extract	00:45	01:00	00:15	00:30
4	Quench with 100 µL 10M formic acid	01:15	01:30	NA	NA
5	Add 10 µL of 30 mM RuBP	NA	NA	03:15	03:30
6	Quench with 100 µL 10 M formic acid	NA	NA	03:45	04:00

- In the interval between quenching initial activities and initiating the reaction for total activities, it is possible to prepare the vials for the next assays (steps 1 and 2 in this table).
- Repeat for all the extractions of the day.
- Blanks can be prepared at the start or end of the day, by adding to each vial **465 µl** of complete assay buffer (from step 2.2.) and **100 µl** 10M formic acid (from step 1.9).
- In addition, it is useful to prepare at least two background checks per experiment by adding to each vial **465 µl** of complete assay buffer (from step 2.2.) , **25 µl** sample extract, and 3 minutes later **100 µl** 10M formic acid (from step 1.9) to test for background levels of carboxylation due to RuBP that may be present in the leaf extracts. In our hands, this tends to be negligible.



We typically do a simple test to verify the total amount of radioactivity present in the complete assay buffer, prepared just before starting the extractions. This serves to verify that the amount of radioactivity in the solution is reliable and comparable across days of assays. For this,
Add **490 µl** 0.3 M KOH (from step 1.8) plus **10 µl** complete assay buffer (from step 2.2). Add **3.6 mL** Gold Star Quanta scintillation cocktail (from step 1.10). Close the vial as soon as possible and mix well.

- Dry all the vials at **100 °C** in the heat block (it takes approximately **01:00:00**).
- Let vials cool, then add **400 µl** ultrapure H₂O to each vial to re-hydrate acid stable compounds.
Wait **00:05:00**
- Add **3.6 mL** of Gold Star Quanta scintillation cocktail (from step 1.10). Close the vials and vortex/mix well.
- Determine ¹⁴C radioactivity using a scintillation counter.

CALCULATIONS

5 Assumptions:

- NaH¹⁴CO₃: 0.5 Ci/mol CO₂ = 0.5 µCi/µmol CO₂
- 1 µCi = 2220000 disintegration per minute (dpm);

Example:

Blank: 89.03 dpm
Vial 1: 23378.90 dpm

1. correct vial DPM value by subtracting background counts (Blank):

$$23378.9 - 89.03 = 23289.87 \text{ dpm}$$

2. convert dpm to µCi: $23289.87 / 2220000 = 0.010491 \text{ µCi}$

3. convert to CO₂ concentration: $0.010491 / 0.5 = 0.0210 \text{ µmol CO}_2$

Rubisco activity is then converted to the unit of interest by accounting for the reaction time, the volume of sample extract used and the corresponding sample leaf area (µmol CO₂ m⁻² s⁻¹) or protein content (µmol CO₂ min⁻¹ mg⁻¹).

From the Rubisco activity calculations above for initial (*V_i*) and total activity (*V_t*), the Rubisco activation state (AS, %) can be calculated:


$$AS = 100 \times V_i/V_t$$