



OCT 02, 2023

## Biochemical assays and evaluation

Iparra<sup>1</sup>

<sup>1</sup>University of California San Diego



Leonardo A Parra-Rivas  
University of California, San Diego

### ABSTRACT

Preparation of brain and Neuro2A lysates

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.5jyl8pey7g2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl8pey7g2w/v1)

**Protocol Citation:** Iparra 2023. Biochemical assays and evaluation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8pey7g2w/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Oct 02, 2023

**Last Modified:** Oct 02, 2023

## Preparation of Brain and Neuro2A Lysates

1

Whole mouse brains were homogenized with a Dounce tissue grinder (20 times) in neuronal protein extraction reagent (N-PER) (Thermo Scientific) containing protease/phosphatase inhibitors (Cell Signaling #5872).

2

Triton X-100 (Sigma #X100-100ML) was added at a final concentration of 1%, and the samples were incubated with rotation for 1 hour at 4 °C.

3

Samples were centrifuged at  $10,000 \times g$  for 10 min at 4 °C, and the supernatant was collected.


4

For Neuro2A lysates, cells were washed with 1X PBS 3 times and incubated for 5 min on ice in the presence of N-PER reagent supplemented with protease inhibitors.

5

Samples were centrifuged at  $10,000 \times g$  for 10 min at 4 °C to remove cellular debris.

6



Total protein

concentration was measured in brain and Neuro2A lysates (DCTM Protein Assay Kit II, Biorad), and samples were used in subsequent experiments.