



APR 05, 2023

OPEN  ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.rm7vzb9x4vx1/v1

Protocol Citation: Grayson Huston 2023. [Modified] DNeasy PowerSoil Pro Kit_Increased Sediment Volume & Optional Concentration. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzb9x4vx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Other
Protocol successful at detecting fish sedDNA collected in a stream during a fish migration
Protocol unsuccessful at detecting fish sedDNA from Maine lakes

Created: Mar 10, 2023

Last Modified: Apr 05, 2023

PROTOCOL integer ID:
78502

[Modified] DNeasy PowerSoil Pro Kit_Increased Sediment Volume & Optional Concentration

Grayson Huston¹

¹Maine Center for Genetics in the Environment, University of Maine



Grayson Huston

ABSTRACT

Protocol (both increased sediment amount up to 2.0g as well as concentrating DNA post-extraction) unsuccessful at detecting fish sedDNA from lakes in Maine, USA.


Both protocols successful at detecting fish sedDNA collected in streams during anadromous fish sea-run migrations

Modified PowerSoil Pro extraction - sample preparation & c.. 31m

1 **CENTRIFUGE** sediment samples briefly to separate pore water

DISCARD pore water to retain only sediment samples

2 **SPIN** the PowerBead Pro Tube briefly to ensure that the beads have settled at the bottom

ADD up to  2.0 g of wet sediment to the PowerBead Pro Tube

Note

Qiagen recommends  .25 g

In this test,  0.25 g,  0.5 g,  1.0 g, and  2.0 g of sediment were tested


ADD  800 μ L Solution CD1



VORTEX briefly to mix

3 **SECURE** PowerBead Pro Tubes horizontally to a 1.5mL-2.0mL Vortex Adapter 20m

VORTEX for  00:10:00

ROTATE tubes so caps are oriented in opposite direction

VORTEX for another  00:10:00

4 **CENTRIFUGE** PowerBead Pro Tube at  15000 x g for  00:01:00 1m

TRANSFER all supernatant to a clean 2 mL Microcentrifuge Tube

Modified PowerSoil Pro extraction - inhibitor removal 31m

5 **ADD**  200 μ L of Solution CD2

VORTEX briefly to mix


6 **CENTRIFUGE** at  15000 x g for  00:01:00

1m


AVOIDING the pellet, transfer **all** supernatant to a clean 2 mL Microcentrifuge Tube

Modified PowerSoil Pro extraction - bind DNA

31m

7 **ADD**  600 µL of Solution CD3

VORTEX briefly to mix

8 **LOAD**  650 µL of the lysate onto a MB Spin Column

1m

CENTRIFUGE at  15000 x g for  00:01:00

DISCARD the liquid flow-through

9 **REPEAT** step 8 to ensure all of the lysate has passed through the MB Spin Column

CAREFULLY place the MB Spin Column into a clean 2 mL Collection Tube

Modified PowerSoil Pro extraction - wash spin column


31m

10 **ADD**  500 µL of Solution EA to the MB Spin column

1m

CENTRIFUGE at  15000 x g for  00:01:00

DISCARD the liquid flow-through and place the MB Spin Column into same 2 mL Collection Tube

11 **ADD**  500 µL of Solution C5 to the MB Spin Column

1m

CENTRIFUGE at  15000 x g for  00:01:00

DISCARD the liquid flow-through and place the MB Spin Column into a **new** 2 mL Collection Tube


12 **CENTRIFUGE** at  16000 x g for  00:02:00

2m

CAREFULLY place the MB Spin Column into a new 1.5mL Elution Tube


Modified PowerSoil Pro extraction - elute the DNA

31m

13 **ADD**  100 µL of Solution C6 to the center of the white membrane in the MB Spin Column

2m

INCUBATE at  Room temperature for  00:01:00

CENTRIFUGE at  15000 x g for  00:01:00

14 **PIPETTE** the liquid flow-through and re-add it to the center of the white membrane in the MB Spin Column

2m

INCUBATE at  Room temperature for  00:01:00

CENTRIFUGE at  15000 x g for  00:01:00

DISCARD the MB Spin Column

DNA is now ready for downstream applications

(Optional) DNA Concentration with PALL Nanosep 30K Cent...

2m

15 **ENSURE** that the sample reservoir is firmly placed into the filtrate receiver

PIPETTE  50-100 µL of DNA extract into the sample reservoir

CAP the Nanosep device

16 **CENTRIFUGE** at  5000 x g for  00:02:00

2m

RECOVER concentrated sample from the sample reservoir with a micropipette

TRANSFER concentrated sample to a new 1.5 mL Microcentrifuge Tube

Concentrated DNA is now ready for downstream applications

Note

If the sample appears to have "spun dry", recover the sample by pipetting ~ 20uL of elution buffer onto the membrane and recovering with a micropipette