



Oct 17, 2021

Evagreen Dye Fluorescence Verification

Chia-Hsien Shih¹¹Chung Shan Medical University

1

[dx.doi.org/10.17504/protocols.io.by5ppy5n](https://doi.org/10.17504/protocols.io.by5ppy5n)

CSMU_Taiwan

Chia-Hsien Shih
Chung Shan Medical University

This protocol is to establish the relationship between Evagreen Dye fluorescent intensity and DNA concentration.

DOI

[dx.doi.org/10.17504/protocols.io.by5ppy5n](https://doi.org/10.17504/protocols.io.by5ppy5n)

Chia-Hsien Shih 2021. Evagreen Dye Fluorescence Verification. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.by5ppy5n>



protocol ,

Oct 17, 2021

Oct 17, 2021

54159

Preparation

- 1 Take **30 μ L** of 10 μ M ssDNA (T4 ligation primer) into a new eppendorf
- 2 Dilute 10 μ M ssDNA into 1 μ M :
Take **3 μ L** of 10 μ M ssDNA and **27 μ L** RNase-free water into a new eppendorf
- 3 Dilute 1 μ M ssDNA into 100nM :
Take **3 μ L** of 1 μ M ssDNA and **27 μ L** RNase-free water into a new eppendorf

- 4 Dilute 100nM ssDNA into 10nM :
Take **3 μ L** of 100nM ssDNA and **27 μ L** RNase-free water into a new eppendorf
- 5 Dilute 10nM ssDNA into 1nM :
Take **3 μ L** of 10nM ssDNA and **27 μ L** RNase-free water into a new eppendorf
- 6 Dilute 1nM ssDNA into 100pM :
Take **3 μ L** of 1nM ssDNA and **27 μ L** RNase-free water into a new eppendorf
- 7 Dilute 100pM ssDNA into 10pM :
Take **3 μ L** of 100pM ssDNA and **27 μ L** RNase-free water into a new eppendorf
- 8 Add **1.5 μ L** of 20X evagreen dye into each eppendorf.

Measure

- 9 Respectively load **20 μ L** of 10 μ M, 1 μ M, 100nM, 10nM, 1nM, 100pM and 10pM ssDNA into each well
- 10 Measure the fluorescence excitation and emission intensity