

Oct 22, 2025

## 2D TEM CLEM (Correlative Light Microscopy and Electron Microscopy)

DOI

[dx.doi.org/10.17504/protocols.io.261gend2jg47/v1](https://dx.doi.org/10.17504/protocols.io.261gend2jg47/v1)

Pietro De Camilli<sup>1,2,3,4,5,6</sup>, Yumei Wu<sup>1,2,3,4</sup>

<sup>1</sup>Department of Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>2</sup>Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>3</sup>Howard Hughes Medical Institute, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>4</sup>Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>5</sup>Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>6</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland, 20815, USA



Yumei Wu

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.261gend2jg47/v1>

**Protocol Citation:** Pietro De Camilli, Yumei Wu 2025. 2D TEM CLEM (Correlative Light Microscopy and Electron Microscopy). protocols.io <https://dx.doi.org/10.17504/protocols.io.261gend2jg47/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** March 03, 2022

**Last Modified:** October 22, 2025

**Protocol Integer ID:** 59054



**Keywords:** Correlative Light Microscopy and Electron Microscopy (CLEM), Transmission Electron Microscopy (TEM), conventional chemical fixation, ASAPCRN, 2d tem clem, 2d transmission electron microscopy, electron microscopy, procedure of correlative light microscopy, correlative light microscopy, general procedure of correlative light microscopy, clem, tem, imaging

## Abstract

This protocol details the general procedure of Correlative Light Microscopy and Electron Microscopy (CLEM) with conventional chemical fixation and 2D Transmission Electron Microscopy (TEM) imaging.

## Attachments



[373-827.docx](#)

15KB

## Materials

### Materials

- 4% PFA
- 0.25% glutaraldehyde
- 2.5% glutaraldehyde
- imaging buffer
- 0.1 M sodium cacodylate
- 2% OsO<sub>4</sub>
- 1.5% K<sub>4</sub>Fe(CN)<sub>6</sub> (Sigma-Aldrich)
- MiliQ water
- 2% aqueous uranyl acetate
- 50% etoh
- 75% etoh
- 95% etoh
- 100% etoh
- 50% Epon
- pure Epon 812
- microtome


















poly-d-lysine-coated 35 mm Dish | No. 1.5 Gridded Coverslip (D) | 14 mm Glass Diameter **MatTek Corporation** Catalog #P35G-1.5-14-CGRD-D

## Troubleshooting



## General preparation:


6h 5m

- 1 Culture cells on poly-d-lysine-coated 35 mm MatTek dish (P35G-1.5-14-CGRD) and transfect the cells with the plasmids of the interest.
- 2 Pre-fix the cells in  37 °C -warmed 4% PFA and 0.25% glutaraldehyde in the imaging buffer.
- 3 Wash in the same imaging buffer. 
- 3.1 Wash  00:05:00 in the same imaging buffer. (1/3) 5m
- 3.2 Wash  00:05:00 in the same imaging buffer. (2/3) 5m
- 3.3 Wash  00:05:00 in the same imaging buffer. (3/3) 5m
- 4 Image with a fluorescence microscopy, find the regions of the interest and their coordinates on the dish, then take the Z-stack fluorescence images and the grids map/coordinate by using phase contrast.
- 5 Fix the cells in 2.5% glutaraldehyde and  0.1 Molarity (M) sodium cacodylate buffer for  01:00:00 . 1h
- 6 Wash in  0.1 Molarity (M) sodium cacodylate buffer. 
- 6.1 Wash in  0.1 Molarity (M) sodium cacodylate buffer for  00:05:00 . (1/4) 5m
- 6.2 Wash in  0.1 Molarity (M) sodium cacodylate buffer for  00:05:00 . (2/4) 5m
- 6.3 Wash in  0.1 Molarity (M) sodium cacodylate buffer for  00:05:00 . (3/4) 5m




- 6.4 Wash in IM 0.1 Molarity (M) sodium cacodylate buffer for 🕒 00:05:00 . (4/4) 5m
- 7 Fix in 2% OsO<sub>4</sub> and 1.5% K<sub>4</sub>Fe(CN)<sub>6</sub> (Sigma-Aldrich) in IM 0.1 Molarity (M) sodium cacodylate buffer for 🕒 01:00:00 . 1h
- 8 Wash in MiliQ water. 🔥
- 8.1 Wash in MiliQ water for 🕒 00:05:00 . (1/4) 5m
- 8.2 Wash in MiliQ water for 🕒 00:05:00 . (2/4) 5m
- 8.3 Wash in MiliQ water for 🕒 00:05:00 . (3/4) 5m
- 8.4 Wash in MiliQ water for 🕒 00:05:00 . (4/4) 5m
- 9 En bloc stained with 2% aqueous uranyl acetate 🕒 Overnight at 🌡 4 °C . 5m 🔄
- 10 Wash in MiliQ water. 🔥
- 10.1 Wash in MiliQ water for 🕒 00:05:00 . (1/4) 5m
- 10.2 Wash in MiliQ water for 🕒 00:05:00 . (2/4) 5m
- 10.3 Wash in MiliQ water for 🕒 00:05:00 . (3/4) 5m
- 10.4 Wash in MiliQ water for 🕒 00:05:00 . (4/4) 5m




11 Dehydration in 50% etoh for  00:05:00 .

5m


12 Dehydration in 75% etoh for  00:05:00 .

5m


13 Dehydration in 95% etoh for  00:05:00 .

5m


14 Dehydration in 100% etoh.

14.1 Dehydration in 100% etoh for  00:10:00 . (1/3)

10m

14.2 Dehydration in 100% etoh for  00:10:00 . (2/3)

10m

14.3 Dehydration in 100% etoh for  00:10:00 . (3/3)


10m

15 Infiltration in 50% Epon for 0.5-1 hour.

16 Infiltration in pure Epon 812.

16.1 Infiltration in pure Epon 812 for  01:00:00 . (1/2)

1h

16.2 Infiltration in pure Epon 812 for  01:00:00 . (2/2)

1h

17 Add two drops of pure Epon into the center of the glass bottom of the culture dish.



18 Polymerize at  60 °C for 1-2 days.



- 19 Cell of interest were relocated based on the pre-recorded coordinates under a stereomicroscope.
- 20 Trim and cut ultrathin sections (50-60 nm) with microtome.
- 21 Observe the ultrathin sections in a electron microscope.

