


Aug 01, 2024

# Protocol for the production and storage of *Diplodia sapinea* pycnidiospores

DOI

[dx.doi.org/10.17504/protocols.io.kqdg3271pv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3271pv25/v1)



Macroscopic and microscopic appearance of a VMM cultivation plate after incubation under constant light (5000–6000 lx) for 21 d on VMM. Pycnidia release spore-containing fluid droplets.

Anne Geertje Oostlander<sup>1</sup>, Laura Brodde<sup>2</sup>, Miriam von Bargaen<sup>1</sup>, Rasmus Enderle<sup>3</sup>, Marco Leiterholt<sup>1</sup>, Dagmar Trautmann<sup>3</sup>, Malin Elfstrand<sup>4</sup>, Jan Stenlid<sup>4</sup>, André Fleißner<sup>1</sup>

<sup>1</sup>Institute of Genetics, Technische Universität Braunschweig, Braunschweig, Germany;

<sup>2</sup>Swedish University of Agricultural Sciences, Uppsala, Sweden;

<sup>3</sup>Julius Kühn Institute (JKI)—Federal Research Centre for Cultivated Plants, Braunschweig, Germany;

<sup>4</sup>Department of Forest Mycology and Plant Pathology, Swedish University of Agricultural Sciences, Uppsala, Sweden



Anne Geertje Oostlander

Technische Universität Braunschweig

OPEN ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.kqdg3271pv25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3271pv25/v1)

**Protocol Citation:** Anne Geertje Oostlander, Laura Brodde, Miriam von Bargaen, Rasmus Enderle, Marco Leiterholt, Dagmar Trautmann, Malin Elfstrand, Jan Stenlid, André Fleißner 2024. Protocol for the production and storage of *Diplodia sapinea* pycnidiospores.

**protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg3271pv25/v1>

## Manuscript citation:

Oostlander AG, Brodde L, Bargaen M von, Leiterholt M, Trautmann D, Enderle R et al. A reliable and simple method for the production of viable pycnidiospores of the Pine pathogen *Diplodia sapinea* and a spore-based infection assay on Scots Pine. Plant Dis 2023; 107(11):3370–7. doi.org/10.1094/PDIS-01-23-0107-RE

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 25, 2024



**Last Modified:** August 01, 2024

**Protocol Integer ID:** 104068

**Keywords:** Diplodia sapinea, Diplodia tip blight, in vitro sporulation, pycnidiospores, Sphaeropsis sapinea

## Abstract

An efficient, standardized protocol for the production and storage of highly viable pycnidiospores of *Diplodia sapinea*.

## Image Attribution

Oostlander AG, Brodde L, Bargen M von, Leiterholt M, Trautmann D, Enderle R et al. A reliable and simple method for the production of viable pycnidiospores of the Pine pathogen *Diplodia sapinea* and a spore-based infection assay on Scots Pine. Plant Dis 2023; 107(11):3370–7.

## Materials

### **Material, media and solutions:**

- Petri dishes (5.5 cm)
- *D. sapinea* (e.g. ex-type strain CBS 138184)
- Light shelf with cold white daylight (Osram Lumilux 18W/865), intensity from 5000 to 6500 lux
- Cleanbench
- Centrifuge

### **Minimal medium (VMM):**

as described in Vogel, 1964

20 ml Vogel's solution

20 g sucrose

15 g agar

add 1 l water autoclave

### **Trace element solution:**

50 g citric acid

50 g zinc sulfate (J.T. Baker)

10 g ammonium iron(II) sulfate

2,5 g copper sulfate

0,5 g manganese sulfate

0,5 g boric acid

0,5 g sodium molybdate

add 1 l water

Chloroform (1 ml) is added as a preservative, store at room temperature.

### **Vogel's solution:**

125 g sodium citrate

250 g potassium dihydrogen phosphate

100 g ammonium nitrate

10 g magnesium sulfate

5 g calcium chloride

5 ml trace element solution

2,5 ml biotin solution

add 1 l water

Chloroform (2 ml) is added as a preservative, store at room temperature.

### **Biotin solution:**

0,1 g biotin

add 1 l water

Store at -20°C



**0,01 % (v/v) Tween:**

0,1 ml Tween

add 1 l water

aliquod and autoclave

**40% (v/v) glycerol**

200 ml glycerol

add 1 l water

aliquod and autoclave



- 1 Inoculate a 5.5 cm Petri dish containing 10 to 15 ml of solid VMM with *D. sapinea* with an agar plug with mycelia from a preculture on VMM.
- 2 Incubate the plate for 21 d at 26 to 29°C with constant light (cold daylight at an intensity of 5000 to 650 lux). Pycnidia will form on the surface of the culture and open to release a droplet of liquid containing spores.
- 3 Harvest the spores from the plate by adding approximately 2 ml of 0.01% (v/v) Tween to the surface and rinse the surface of the colony several times by pipetting. Repeat the procedure for a higher yield. Typically, between  $1 \times 10^6$  and  $2 \times 10^6$  spores can be harvested per plate. Spores can be separated from the liquid by centrifugation for 10 s at 5000 rpm and then resuspended in your desired volume.
- 4 For storage, mix 500  $\mu$ l of the spore suspension and 500  $\mu$ l 40% glycerol in a 2 ml screw top tube or cryovial. Freeze the glycerol stock tube at -80°C.

## Protocol references

Vogel HJ. A convenient growth medium for *Neurospora crassa* (Medium N). Microbial genetics bulletin 1956; (13):42–3.