



Aug 27, 2020

BG11 hypersaline agar plates

Bosak Lab Protocols¹¹Massachusetts Institute of Technology**1** Works for me dx.doi.org/10.17504/protocols.io.bkchkst6

Bosak Lab

Emilie J Skoog
Massachusetts Institute of Technology

DOI

dx.doi.org/10.17504/protocols.io.bkchkst6

PROTOCOL CITATION

Bosak Lab Protocols 2020. BG11 hypersaline agar plates. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bkchkst6>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 26, 2020

LAST MODIFIED


Aug 27, 2020

PROTOCOL INTEGER ID

41065

- 1 Add the following to a beaker and dissolve using a stir bar on a magnetic stir plate.

BG11 hypersaline agar plate medium	amount	unit
initial Milli-Q water	968	ml
NaCl	49.8	g
NaNO3	1.5	g
Na2CO3	0.02	g
KCl	1.3867	g
MgCl2*6H2O	8.569	g
MgSO4*7H2O	6.4998	g
CaCl2*2H2O	2.042	g
Stock A	10	ml
Stock B	10	ml
Stock C	10	ml
Stock 5 (trace metals)	1	ml
Vitamins	1	ml
Agar	15	g

- 2 Add Milli-Q water to bring total volume to  1 L
- 3 Autoclave media 40m
- 4 While waiting for autoclave cycle to finish, set up desired number of sterile petri dishes or other culturing plates in a sterile biohood.
- 5 Remove media from autoclave and let cool until just warm enough to touch.
- 6 Pour media into each petri dish and let solidify.*

*TIP: Before agar solidifies, use a sterile pipette tip to move any bubbles that may form to the edges of the petri dish.