



Version 3 ▾

Sep 13, 2022

ONT Q20+ (V12) Adapter Ligation for Fungal DNA Barcoding V.3

In 2 collections

Stephen Douglas Russell¹, Stephen Douglas Russell¹¹The Hoosier Mushroom Society

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.dm6gpb5zdlzp/v3

The Hoosier Mushroom Society

Stephen Douglas Russell

ABSTRACT

This process will take your A-tailed library and add the nanopore adapters. Simply combine several chemicals for a single reaction and do a bead cleanup.

Tested with:

Flowcells: Flongle 9.4.1 or MinION 10.4.1

Ligation Kit: V12 - LSK112

Time required: ~45 minutes

DOI

dx.doi.org/10.17504/protocols.io.dm6gpb5zdlzp/v3

PROTOCOL CITATION

Stephen Douglas Russell, Stephen Douglas Russell 2022. ONT Q20+ (V12) Adapter Ligation for Fungal DNA Barcoding. **protocols.io**
<https://protocols.io/view/ont-q20-v12-adapter-ligation-for-fungal-dna-barcode-cghrtt56>
Version created by Stephen Douglas Russell

COLLECTIONS ⓘ

ONT DNA Barcoding Fungal Amplicons w/ MinION & Flongle **ONT DNA Barcoding Fungal Amplicons w/ MinION & Flongle**

KEYWORDS

nanopore, fungi, flongle, fmol, library preparation

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 13, 2022

LAST MODIFIED

Sep 13, 2022

PROTOCOL INTEGER ID

69905

PARENT PROTOCOLS

Part of collection

[ONT DNA Barcoding Fungal Amplicons w/ MinION & Flongle](#)

[ONT DNA Barcoding Fungal Amplicons w/ MinION & Flongle](#)

MATERIALS TEXT

Reagents

[Ligation Sequencing Kit \(Q20 \) Oxford Nanopore](#)

Technologies Catalog #SQK-LSK112 In 5 steps

: \$694.43 per 6 reactions

[NEBNext Quick Ligation Module New England](#)

Biolabs Catalog #E6056S Step 1

: \$361.00 per

20 reactions

*note: This kit has two components. We use one. NEB checking on whether the single one is available for purchase. Samples of this kit should be available to start.

[HighPrep™ PCR Clean-up System MagBio Genomics](#)

Inc. Catalog #AC-60005

: \$117.88 per 50

mL. \$0.047 per rxn.

Total per Flongle run (1/2 rxns): \$66.95

Total per MinION run: \$133.89

Total per 96 samples: \$13.38

Total per sample (Flongle: 480 samples): \$0.139

Consumables

Eppendorf DNA LoBind 1.5mL tubes

10uL pipette tips

100-200uL pipette tips

Equipment

PCR tube rack

Vortex mixer

Mini centrifuge

PCR cleanup magnet

10uL Pipette

100uL Pipette

Hula mixer (Ebay): \$200.00 (optional)

Quantus or Qubit Fluorometer (optional)

Adapter Ligation

- 1 Spin down the Adapter Mix H (AMX H) and Quick T4 Ligase, and place on ice.

AMX H -

[Ligation Sequencing Kit \(Q20\) Oxford Nanopore](#)

Technologies Catalog #SQK-LSK112

[NEBNext Quick Ligation Module New England](#)

Quick T4 Ligase - **Biolabs Catalog #E6056S**

- 2 Thaw Ligation Buffer (LNB) at room temperature, spin down and mix by pipetting. Due to viscosity, vortexing this buffer is ineffective. Place on ice immediately after thawing and mixing.

LNB -

[Ligation Sequencing Kit \(Q20\) Oxford Nanopore](#)

Technologies Catalog #SQK-LSK112

- 3 Thaw the Elution Buffer (EB) at room temperature, mix by vortexing, spin down and place on ice.

EB - Lig

[Ligation Sequencing Kit \(Q20\) Oxford Nanopore](#)

Technologies Catalog #SQK-LSK112

- 4 Thaw one tube of Short Fragment Buffer (SFB) at room temperature, mix by vortexing, spin down and place on ice.

SFB -



[Ligation Sequencing Kit \(Q20\) Oxford Nanopore](#)

Technologies Catalog #SQK-LSK112





- 5 In a 1.5 ml Eppendorf DNA LoBind tube, mix in the following order:

Between each addition, pipette mix 10-20 times.

Reagent	9.4.1 Flongle Volume	R10.4.1 MinION Volume
DNA sample from the previous step	30 µl	60 µl
Ligation Buffer (LNB)	12.5 µl	25 µl
NEBNext Quick T4 DNA Ligase	5 µl	10 µl
Adapter Mix H (AMX H)	2.5 µl	5 µl
Total	50 µl	100 µl

- 6  00:00:05 Spin down with a mini centrifuge for  00:00:05 .

10s

- 7 Incubate the reaction for ⌚00:10:00 at room temperature. 10m
- 8 Resuspend magnetic bead stock by vortexing.
- 9 Add 20 µL (Flongle) or 40 µL (MinION) of resuspended beads to the reaction and mix by flicking the tube.
- 10 Incubate on a Hula mixer (rotator mixer) for ⌚00:05:00 at room temperature. 5m
- 11 Spin down the sample for ⌚00:00:05 and pellet on a magnet for ⌚00:02:00 . Keep the tube on the magnet, and pipette off the supernatant. 2m 5s
- 12 Wash the beads by adding 150 µL (Flongle) or 250 µL (MinION) of Short Fragment Buffer (SFB). Flick the beads to resuspend, spin down for ⌚00:00:05 , then return the tube to the magnetic rack for ⌚00:02:00 and allow the beads to pellet. Remove the supernatant using a pipette and discard. 2m 5s
- Note: flicking the tube does not seem to fully resuspend the beads. Just flick 10 times or so.
- SFB -
[🔗Ligation Sequencing Kit \(Q20\) Oxford Nanopore Technologies Catalog #SQK-LSK112](#)
- 13 Repeat the previous step. ↺
- 14 Spin down for ⌚00:00:05 and place the tube back on the magnet. Pipette off any residual supernatant. Allow to dry for ~30 seconds, but do not dry the pellet to the point of cracking. 5s
- 15 Remove the tube from the magnetic rack and resuspend the pellet in 15 µl Elution Buffer (EB). Incubate for ⌚00:10:00 at room temperature. 10m

16 Pellet the beads on a magnet until the eluate is clear and colorless, for at least 00:01:00 .

1m

17 Remove and retain 15 µl of eluate containing the DNA library into a clean 1.5 ml Eppendorf DNA LoBind tube.

Store on ice until you are ready to load in your flowcell.

Quantification

18 If you have access to a Quantus or Qubit fluorimeter, now is a good time to quantify 2uL of DNA in your sample.

It is recommend loading 5 fmol to 10 fmol of this final prepared library onto your flow cells. Loading more than 20 fmol of DNA can reduce the rate of duplex read capture. Dilute the library in Elution Buffer if required.

<https://www.promega.com/resources/tools/biomath/>

For 900bp length DNA (what our ITS1F-4 rxns appear to average), we are looking for

Flongle R9.4.1: 5 fmol - 20 fmol = .003ug - .012ug of DNA.

MinION R9.4.1: 5 - 50 fmol = .003ug - .029ug of DNA.

MinION R10.4.1: 25 - 75 fmol = .015ug - .044ug of DNA.

For a 22 ng/uL sample (Quantus quantification):

$$\frac{22\text{ng}}{1\text{uL}} \times \frac{1\text{ug}}{1000\text{ng}} = 0.022\text{ug/uL} \times 13\text{uL (elution buffer; or x15 if you did not quantify using 2uL)} = 0.286\text{ ug DNA in}$$

sample of 13uL elution buffer.

*Note: the 0.286 in the calculations below will change based on your individual DNA amount.

**Also note: The ONT protocol suggests using additional EB in order to make concentration adjustments. As there is not a lot of the reagent in the standard packets, and it is not possible to buy more individually, I have been using molecular water for this step with no ill effect.

Flongle 9.4.1

How much additional molecular water to have 5.5uL needed for the next step give us correct amount of DNA?

$$0.286\text{ug} / \text{xuL} = 0.010\text{ug (17 fmol DNA)} \quad \text{x} = 28.6\text{uL} \quad \text{x} 5.5\text{uL} = 157.3\text{uL} - 13\text{uL (or 15 if you did not quantifly)} = 144.3\text{uL}$$

So at 0.022ug/uL quantification, add an additional 144.3uL of molecular water to have right concentration to use 5.5uL for the next step with Flongle.

11ng/uL sample comes out to adding an additional 65.7 uL of molecular water.
31ng/uL sample comes out to adding an additional 208.7 uL of molecular water.

I would use 150 uL of extra molecular water if you are not able to quantify your sample.

MinION R9.4.1

$0.33\text{ug} / \text{xuL} = 0.025\text{ug}$ (42 fmol DNA) $\times = 13.2\text{uL}$ $\times 11\text{uL} = 145.2\text{uL} - 15\text{uL} = 130\text{uL}$ elution buffer addition.

MinION R10.4.1

$0.33\text{ug} / \text{xuL} = 0.04\text{ug}$ (67 fmol DNA) $\times = 8.25\text{uL}$ $\times 11\text{uL} = 90.75\text{uL} - 15\text{uL} = 75\text{uL}$ elution buffer addition.