



Sep 06, 2022

# Preparing of genomic DNA from in vitro cultured cells

In 1 collection

Hanqin Li<sup>1</sup>, Yogendra Verma<sup>1</sup>, Dirk Hockemeyer<sup>1</sup>, Frank Soldner<sup>2</sup><sup>1</sup>University of California, Berkeley; <sup>2</sup>Albert Einstein College of Medicine

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.b4ntqven](https://dx.doi.org/10.17504/protocols.io.b4ntqven)

Devin E Snyder

## ABSTRACT

This protocol describes a standard procedure for preparing crude cell lysate which can be further analyzed by PCR.

## Protocol overview:

- A. Preparing crude cell lysate directly from hPSCs culture
- B. Preparing crude cell lysate from dissociated cells
- C. Preparing cell lysate from FACS-sorted cells (small cell number)

## General notes

1. Throughout this protocol, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.

## DOI

[dx.doi.org/10.17504/protocols.io.b4ntqven](https://dx.doi.org/10.17504/protocols.io.b4ntqven)

## PROTOCOL CITATION

Hanqin Li, Yogendra Verma, Dirk Hockemeyer, Frank Soldner 2022. Preparing of genomic DNA from in vitro cultured cells. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b4ntqven>



## FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's  
Grant ID: ASAP-000486

## COLLECTIONS ⓘ

**Genotyping by next generation sequencing**

## KEYWORDS

ASAPCRN

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Feb 03, 2022

## LAST MODIFIED

Sep 06, 2022

## PROTOCOL INTEGER ID

57779

## PARENT PROTOCOLS

Part of collection

[Genotyping by next generation sequencing](#)

## MATERIALS TEXT

| Item  | Vendor           | Catalog # |
|---|------------------|-----------|
| DMEM/F12                                    | Thermo<br>Fisher | 11320082  |
| DPBS<br>w/o Calcium and<br>magnesium (DPBS) | Corning          | MT21031CV |
| Newborn Calf Serum                          | Sigma            | N4762     |
| Penicillin & Streptomycin<br>(100X)         | Thermo<br>Fisher | 15140163  |
| 0.25% Trypsin with EDTA<br>(Trypsin)        | Thermo<br>Fisher | 25200114  |
| Proteinase K                                | Sigma            | P6556     |
| Exact N Amp Blood PCR<br>Kit                | Sigma            | XNAT2-1KT |

## A. Preparing crude cell lysate directly from hPSCs culture

15m

- 1 Add 500 µl crude lysis buffer (Proteinase K added) to a fully confluent well of a 6 -well plate (~1.5 million cells)

### 1.1

#### Crude lysis buffer

| A                             | B         |
|-------------------------------|-----------|
| KCl                           | 50 mM     |
| MgCl <sub>2</sub>             | 2 mM      |
| NP-40                         | 0.45%     |
| Tween-20                      | 0.45%     |
| Tris                          | 10 mM     |
| Proteinase K (add before use) | 250 µg/ml |

- 2 Incubate at **37 °C** **00:10:00**

10m

- 3 Collect the crude lysate into a 1.7 ml Eppendorf tube.

- 4 Incubate at **50 °C** **Overnight**

- 5 Transfer 50 µl to a 200 µl microcentrifuge tube

- 6 Heat inactivate at **95 °C** for **00:05:00** in a thermocycler.

5m

- 7 Chill **On ice**

8 The crude cell lysate is ready for PCR

B. Preparing crude cell lysate from dissociated cells 10m

9 Wash cells with DPBS once

10 For each well in a 6-well plate, use 0.5 ml Trypsin

11 Incubate  **37 °C**  **00:05:00** 5m

12 Inactivate Trypsin using 2 ml wash medium

12.1 Wash medium

| A                                | B      |
|----------------------------------|--------|
| DMEM/F12                         | 470 ml |
| Newborn Calf Serum               | 25 ml  |
| Penicillin & Streptomycin (100X) | 5 ml   |

Final volume; 500 ml

13 Centrifuge at  **200-300 x g, Room temperature, 00:05:00** 5m

14 Remove supernatant

15 For a fully confluent well of a 6-well plate, mix the cell pellet with 500 µl crude lysis buffer (Proteinase K added). Reduce the amount of crude lysis buffer if the culture is not fully confluent.

16 Transfer the crude lysate into a 1.7ml Eppendorf tube.

17 Incubate at  **50 °C**  **Overnight**

18 Transfer 50 µl to a 200 µl microcentrifuge tube

19 Heat inactivate at  **95 °C**  **00:05:00** in a thermocycler

5m

20 Chill  **On ice**

21 The crude cell lysate is ready for PCR

#### C. Preparing cell lysate from FACS-sorted cells (small cell number)

22 For small cell numbers, we usually use a direct cell to DNA extraction kit (Exact N Amp Blood PCR Kit [Sigma, # XNAT2-1KT]) according to the manufacturers' instructions.