



MAY 08, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.5jyl8j5zrg2w/v1

Protocol Citation: Dan Dou, C. Alexander Boecker, Erika L.F. Holzbaur 2023. Cell lysis and immunoblotting for protein and phospho-protein quantification. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8j5zrg2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Oct 13, 2022

Last Modified: May 08, 2023

PROTOCOL integer ID:
 71304

Keywords: Lysate, Immunoblot, Western, Phospho, LRRK2


Cell lysis and immunoblotting for protein and phospho-protein quantification

Dan Dou^{1,2}, C. Alexander Boecker³, Erika L.F. Holzbaur^{1,2}

¹Department of Physiology, Perelman School of Medicine, University of Pennsylvania, Philadelphia, PA 19104, USA;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, USA;

³Department of Neurology, University Medical Center Goettingen, 37077 Goettingen, Germany



Dan Dou

ABSTRACT

Here, we describe the procedure by which human iPSC-derived neurons or mouse embryonic fibroblasts (MEFs) were lysed and probed for levels of proteins of interest using Western blot.

ATTACHMENTS




[552-1148.pdf](#)

MATERIALS

Reagents

■ RIPA buffer

A	B
Tris-HCl	50 mM
NaCl	150 mM
Triton X-100	0.1%
Deoxycholate	0.5%
SDS	0.1%

-  HALT phosphatase and protease inhibitor cocktail (100x) Thermo Fisher Scientific Catalog #78442
-  Microcystin-LR Microcystis aeruginosa - CAS 101043-37-2 - Calbiochem Merck MilliporeSigma (Sigma-Aldrich) Catalog #475815m
-  Pierce BCA Protein Assay Kit Thermo Fisher Scientific Catalog #23225

■ 4x Protein Loading Buffer

A	B
Tris-HCl, pH 6.8	125 mM
Glycerol	50%
SDS	4%
Orange G	0.2%

■ Acrylamide

■ 4x Running buffer




A	B
Trizma base	48 g
Glycine	230.4 g
NaN ₃	20 mL
ddH ₂ O	Diluted to 4 L

■ Running buffer

A	B
4x running buffer	250 mL
ddH ₂ O	750 mL
10% SDS	10 mL

■ Transfer buffer

A	B
4x running buffer	125 mL
ddH2O	875 mL
10% SDS	500 µL
For RABs add 20% Methanol	


-  Immobilon®-FL PVDF Membrane Merck MilliporeSigma (Sigma-Aldrich) Catalog #ipfl00010
-  Chameleon® Duo Pre-stained Protein Ladder LI-COR Catalog #928-60000
-  Revert™ 700 Total Protein Stain for Western Blot Normalization (250 ml) LI-COR Catalog #926-11021

■ Revert Wash Solution

A	B
Acetic acid	6.7%
Methanol	30%
in ddH2O	

■ Revert Reversal Solution

A	B
NaOH	0.1 M
Methanol	30%
in ddH2O	

-  EveryBlot Blocking Buffer 500 ml Bio-Rad Laboratories Catalog #12010020
- Primary antibodies (see Materials and Methods for specific antibodies used)
- Secondary antibodies (see Materials and Methods for specific antibodies used)

Equipment

- ODYSSEY CLx Imaging System (LI-COR)

Equipment	
Mini-PROTEAN Tetra Vertical Electrophoresis Cell	NAME
Electrophoresis Cell	TYPE
Bio-Rad	BRAND
1658004	SKU
https://www.bio-rad.com/en-in/product/mini-protean-tetra-vertical-electrophoresis-cell?ID=N3F2UD4VY	LINK

Equipment	
Mini Trans-Blot Electrophoretic Transfer Cell	NAME
Electrophoretic Transfer Cell	TYPE
Bio-Rad	BRAND
1703930	SKU
https://www.bio-rad.com/en-in/sku/1703930-mini-trans-blot-electrophoretic-transfer-cell?ID=1703930	LINK

SAFETY WARNINGS




- Microcystin-LR is an extremely potent hepatotoxin and should be handled with great care.
- Acrylamide is a neurotoxin and should be handled with care.
- Methanol-containing reagents should be handled carefully, as methanol can penetrate single-layer laboratory gloves.

Preparation of cell lysates




- 1 Quickly wash cells twice with ice-cold PBS. After the second wash, tilt the dish and completely aspirate all residual PBS.



- 2 Immediately add ice-cold lysis buffer, ensuring that the entire surface is covered by lysis buffer. Place cells  On ice.



Note


The amount of lysis buffer to use depends on cell confluency / cell number, cell type, and cell culture dish. In most cases, using  100 μL –  150 μL lysis buffer per well of a 6-well plate should result in a protein concentration >  1 $\mu\text{g}/\mu\text{L}$.

Note

Halt protease and phosphatase inhibitor cocktail and microcystin-LR should be added fresh on the day of use.



- 3 Scrape cells off the dish using a cell lifter.

3.1


Transfer the cell lysate to an Eppendorf tube  On ice .

3.2

20m

Leave lysates  On ice for  00:20:00 to allow for efficient lysis.

Note

Cell lysates can be snap frozen in liquid nitrogen and stored at  -80 °C for future use.

4



Centrifuge for  00:10:00 at  17000 x g and  4 °C .

10m



4.1



Discard pellet and use clarified supernatant to determine protein concentration by BCA assay following the manufacturer's instructions, performing all measurements in triplicates.



5



Add  100 µL β-mercaptoethanol to  900 µL of 4x Protein Loading Buffer and mix well.

5.1



Add complete 4x Protein Loading Buffer to cell lysates, mix well, and boil for  00:05:00 at  95 °C .

5m

Note

Do not store Protein Loading Buffer with BME for more than two weeks.


SDS-polyacrylamide gel electrophoresis

- 6 Load samples onto 8% (for LRRK2 protein) to 15% (for PPM1H and Rab proteins) acrylamide gels alongside Chameleon Duo pre-stained protein ladder (LI-COR).





Note

Carefully rinse wells with running buffer before loading cell lysates.

- 7 Start electrophoresis at 80 V for  00:20:00, then increase to 120 V and electrophorese until orange dye runs out.

20m

Protein transfer

- 8 Activate Immobilon-FL PVDF membrane by submerging in methanol for  00:00:30 -  00:01:00.

1m 30s

- 8.1 Wash in ddH₂O and equilibrate in transfer buffer.



- 8.2 Soak sponges in methanol, wash in ddH₂O and equilibrate in transfer buffer.



- 8.3 Equilibrate filter paper in transfer buffer.


- 8.4 Assemble blotting sandwich.

8.5

Carefully remove any air bubbles between layers using a roller.


9 Fill transfer tank with ice-cold transfer buffer.



9.1 Place transfer system  On ice .

9.2





Transfer proteins from gel onto PVDV membrane at 100 V for  01:15:00 .



1h 15m

Total protein stain, membrane blocking, and antibody incub.

1h 41m



10 After wet-tank transfer, let membrane dry completely for at least  01:00:00 at  Room temperature .

1h

11 Rehydrate membrane for  00:01:00 in 100% methanol, then wash  00:05:00 in 1x TBS.

6m



12 Incubate for  00:05:00 in Revert total protein stain (LI-COR) while gently shaking at  Room temperature .

5m



12.1



Wash twice with Revert wash solution, then rinse in ddH₂O and image membrane on ODYSSEY CLx imaging system.

13

Remove Revert total protein stain by incubating membrane in Revert reversal solution for

10m

⌚ 00:10:00 at 🌡 Room temperature while gently shaking.



13.1

Rinse membrane with ddH₂O.

14

Block membrane for ⌚ 00:05:00 in Everyblot Blocking Buffer (Bio-Rad) at

5m

🌡 Room temperature .

15

Dilute primary antibodies in Everyblot Blocking Buffer and incubate at 🌡 4 °C ⌚ Overnight .

5m



16

Wash membrane in 1x TBS + 0.1% Tween-20 (TBS-T) at 🌡 Room temperature (4 washes,

5m

⌚ 00:05:00 each).



17

Dilute secondary antibodies 1:20,000 in Everyblot Blocking Buffer and 0.02% SDS.

**Note**

Incubate membrane in secondary antibodies for ⌚ 01:00:00 at 🌡 Room temperature .

18

Wash membrane in TBS-T at 🌡 Room temperature (4 washes, ⌚ 00:05:00 each).

5m



19



Rinse membrane with TBS (no detergent), then image on ODYSSEY CLx imaging system. Quantify signal intensity using Image Studio Software.