



Jun 02, 2022

# Yeast cells live fluorescence imaging

Florian Wilfling<sup>1</sup>, Fabian Fiedler<sup>1</sup>, Anna Bieber<sup>2</sup>, Cristina Capitanio<sup>2</sup><sup>1</sup>Max Planck Institute of Biophysics, Frankfurt; <sup>2</sup>Max Planck Institute of Biochemistry

1

[dx.doi.org/10.17504/protocols.io.n92ldzqjnv5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldzqjnv5b/v1)**Cristina Capitanio**  
Max Planck Institute of Biochemistry

Protocol for sample preparation and live fluorescence imaging of yeast cells.

DOI

[dx.doi.org/10.17504/protocols.io.n92ldzqjnv5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ldzqjnv5b/v1)

Florian Wilfling, Fabian Fiedler, Anna Bieber, Cristina Capitanio 2022. Yeast cells live fluorescence imaging. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.n92ldzqjnv5b/v1>



#yeast, #fluorescence, #glass slides coating

 protocol ,

May 24, 2022

Jun 02, 2022

63135

#### MEDIA

- LoFlo medium (Yeast Nitrogen Base w/o Amino Acids and w/o Folic Acid and Riboflavin. Formedium)

#### REAGENTS

- Concanavalin A solution: 1 mg ml<sup>-1</sup> Concanavalin A dissolved in purified water. (sterile filtered 0.2 µm)

#### INSTRUMENTS

- Nikon Ti2 Eclipse microscope, with Olympus Apo TIRF 100x 1.49 oil objective and Hamamatsu ORCA-Flash 4.0 LT+ Digital CMOS camera

#### SOFTWARE

- CellCounter plugin, ImageJ

### Yeast cell culture in low fluorescence medium 12h

1 Grow yeast cells overnight in Low Fluorescence (LoFlo) medium. 12h

1.1 In the morning dilute cells in LoFlo to OD 0.1.

1.2 Grow until mid-log phase (0.5–0.8 OD).

### Pre-treat microscope slides with Concanavalin A solution 1h

2 Add 10 µl Concanavalin A to the glass surface.

2.1 Let the slides dry completely before use.

### Prepare the sample for imaging 20m

3 Once the cell culture has reached 0.5–0.8 OD, spin down 2 ml of the original culture. Resuspend the pellet in 30 µl LoFlo to obtain a concentrated sample.

3.1 Apply 3.5 µl of the concentrated yeast culture on the Concanavalin A-treated glass slide.

3.2 Cover the sample with a glass coverslip.

Imaging procedure: 1h

4 First, choose a group of cells to use to set up adequate imaging settings. Use lasers corresponding to the right fluorophores and regulate exposure time based on the brightness of the features of interest.

4.1 Set up stack acquisition. For yeast cells set a total stack size of 6-8 µm. We used a step size of 300 nm.

4.2 In order to have an unbiased approach, use Transmitted Light to select a new a field of view for data acquisition.  
The field of view should contain 25-100 cells in the same focal plane.

4.3 Acquire stacks.

Image analysis

5 Count cells and analyse signal with Fiji Cell counter Plug-in.