



Jun 02, 2022

LRRK1 Immunoprecipitation kinase assay

Asad Malik¹, Athanasios Karapetsas¹, Dario R Alessi¹¹Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Street, Dundee DD1 5EH, UK

1

dx.doi.org/10.17504/protocols.io.kqdg3p2wpl25/v1

asap

Dario Alessi

We describe a non-radioactive assay that we deploy for analysing LRRK1 protein kinase activity *in vitro* using Rab7A as a substrate. This assay can be used to measure the intrinsic activity of LRRK1 immunoprecipitated from cells as well as to assess the impact of mutations on LRRK1 activity.

[382-844.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.kqdg3p2wpl25/v1

Asad Malik, Athanasios Karapetsas, Dario R Alessi 2022. LRRK1 Immunoprecipitation kinase assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg3p2wpl25/v1>



transfection of HEK293, LRRK1 activity, Immunoprecipitation , ASAPCRN

protocol ,

Mar 21, 2022

Jun 02, 2022

Mar 21, 2022 renuka.s

Mar 25, 2022 Dario Alessi




Figure 1: Optimisation of recombinant Rab7A phosphorylation *in vitro* by immunoprecipitated GFP-LRRK1. HEK293 cells were transfected with GFP-LRRK1 WT. The cells were harvested 24 h post transfection in lysis buffer containing either 1% Triton, 1% Digitonin or 0.3% CHAPS and supplemented with 1 μ M GTP- γ -S as indicated. 1mg of cell extract was subjected to a GFP immunoprecipitation, followed by a 45-minute kinase reaction at 30°C in the presence of 5 μ g recombinant Rab7A and excess Mg-32P-ATP. 75% of kinase reactions were separated on SDS-PAGE. Following electrophoresis, gels were fixed (50% (v/v) methanol, 10% (v/v) acetic acid), stained in Coomassie brilliant blue, dried and exposed to a phospho-imaging screen for assessing radioactive 32P incorporation.




Figure 2: *In vitro* phosphorylation of recombinant Rab7A by immunoprecipitated GFP-LRRK1 WT. HEK293 cells were transfected with GFP-LRRK1 WT or GFP-LRRK1[D1409A] (kinase inactive mutant). The cells were lysed 24 h post transfection and the 1mg of cell extract were subjected to a GFP immunoprecipitation, followed by a 45 minute kinase reaction at 30°C in the presence of 5 µg recombinant Rab7A and excess Mg-ATP. 25% of GFP-LRRK1 WT and GFP-LRRK1[D1409A] kinase reactions were subjected to immunoblot analysis with the indicated antibodies and the membranes were developed using the Odyssey CLx scan Western Blot imaging system.

Reagents:

For transfection of GFP-tagged LRRK1 in HEK293 cells:

cDNA for expression of human GFP-tagged LRRK1 in mammalian cells; Polyethylenimine "Max" (MW 4,000) (Polysciences, Inc., cat no 24885): **1 mg/mL** stock in de-ionised H₂O, filtered.

Lysis buffer:

A	B
HEPES pH 7.5*	50 mM
(v/v) CHAPS hydrate*	0.30%
Na ₃ VO ₄ **	1 mM
NaF	50 mM
β-glycerophosphate	10 mM
sodium pyrophosphate	5 mM
sucrose	0.27 M
cOmplete™, EDTA-free Protease Inhibitor Cocktail (Roche, 11836170001)**	
GTP-γ-S**	1mM
Microcystin-LR (Enzo Life Sciences, ALX-350-012)**	1 µg/ml

*Prior optimisation has shown that LRRK1 best retains its kinase activity when cells are harvested in a lysis buffer containing 0.3% CHAPS with 50mM HEPES pH 7.5 (Figure 1).

**To be added fresh before use.

[Polyethylenimine Hydrochloride Linear \(MW 4000\) Polysciences](#)

Inc Catalog #24885

Bradford assay kit

[Pierce™ Coomassie Plus \(Bradford\) Assay Kit Thermo](#)

Fisher Catalog #23236

Resin for LRRK1 immunoprecipitation: aGFP16-aGFP2-His6 NHS-activated Sepharose beads (available from MRC Reagents and Services: <https://mrcpureagents.dundee.ac.uk/>) for GFP-tagged LRRK1

IP wash buffers: Lysis buffer supplemented with **[M]300 millimolar (mM)** NaCl;

[M]50 millimolar (mM) HEPES **pH 7.5**.

Kinase assay buffer:

A	B
HEPES pH 7.5	50 mM
MgCl ₂	10 mM
ATP	1 mM

Recombinant Rab7A protein (available from MRC Reagents and Services: <https://mrcpureagents.dundee.ac.uk/>)

4X Loading buffer:

[NUPAGE LDS sample buffer \(4x\) Thermo Fisher](#)

Scientific Catalog #NP0007

or

4X SDS loading buffer:

A	B
Tris-HCl, pH6.8	250 mM
SDS	8% (w/v)
Glycerol	40% (v/v)
Bromophenol blue	0.02% (w/v)

Equipment:

Refrigerated Centrifuge
Centrifuge

Eppendorf EP-5417R [↗](#)

Plate reader for Protein quantification (BioTek Epoch, or equivalent)
Thermo mixer (Eppendorf ThermoMixer, or equivalent)

Digital Dry Baths/Block Heaters
Dry bath/heat block

Thermo Scientific™ 88870005 [↗](#)

Transient transfection of HEK293 cells for analysis of over-expressed LRRK1 activity in vitro:

1



Transfect HEK293 cells at around 60-70% confluency. For a 10cm dish, add **10 µg** DNA (GFP-tagged LRRK1 or GFP-empty vector) and **30 µL** of **1 mg/mL** PEI solution to **1 mL** of Opti-MEM™ Reduced Serum Medium and vortex for 20/30 seconds.

Note: We recommend including a GFP-empty vector transfection as well as a GFP-tagged LRRK1 D1409A (kinase dead) transfection to control for specificity of LRRK1 immunoprecipitation and activity.

2



20m

Incubate at **Room temperature** for **00:20:00** to allow the DNA/PEI complexes to form.

3



Add the transfection mix to the culture medium in each dish and incubate cells at **37 °C** after transfection.

4

Lyse cells 20-24 hours after transfection.

Preparation and quantification of cell lysates:

5



Quickly rinse cells in the tissue culture dish by carefully pouring **Room temperature** culture media without Foetal bovine serum (FBS) into the dish.

6

Pour off media from the culture dish and completely aspirate any residual media. Immediately add freshly prepared ice-cold lysis buffer, ensuring that the entire surface is covered by lysis buffer.

Note: The amount of lysis buffer to use will depend on cell type. As a guideline, use **400 mL** of lysis buffer for a 10 cm dish for HEK293 cells.

7

Immediately transfer the cell dishes to ice.

8


Scrape the cells on the dish using a cell lifter (Sigma-Aldrich CLS3008, or equivalent) to ensure


all cells are detached from the dish.

9 Using a pipette, transfer cell lysate to an Eppendorf tube  **On ice** .

10 Leave samples  **On ice** for  **00:20:00** to allow for efficient lysis. 20m


11  10m

Clarify lysates by centrifugation at  **20800 x g, 4°C, 00:10:00** .


12 Transfer the supernatants into new Eppendorf tubes and discard the pellet. Keep the tubes  **On ice** .

Note: Cell lysates can be snap frozen in liquid nitrogen and stored at -80°C for future use. When assessing kinase activity of LRRK1 immunoprecipitated from cells, we do not recommend more than one freeze/thaw cycle.



13 Determine the protein concentration of cell lysates by Bradford assay according to the manufacturer's instructions, performing measurements in triplicate.

Note: Ensure the concentration of the samples is in the linear range for the Bradford assay. If it isn't, prepare appropriate dilutions in water of each lysate. Generally, protein concentrations of near confluent cells lysed as described above should result in protein concentrations of at least  **2 ug/ul** .

Immunoprecipitation of over-expressed LRRK1 from HEK293 cells:

14 

Note: When comparing multiple GFP-tagged variants of LRRK1, we recommend assessing the levels of LRRK1 in the lysates prior to immunoprecipitation by subjecting ~10 ug cell extract to immunoblotting, normalizing total LRRK1/Tubulin levels and adjusting how much cell lysate is to be used to immunoprecipitate LRRK1 based on this quantification, to ensure that the amount of enzyme between reactions is as close as possible.

Add  **20 mL** of aGFP16-aGFP2-His6 NHS-activated Sepharose beads (washed 3 times in PBS and resuspended in PBS to make a 1:1 slurry) to  **1 mg** of cell extract.

Note: The immunoprecipitation conditions (amount of resin and amount of cell extract) might need optimisation.

15



3h

Incubate at  **4 °C** for  **03:00:00**, under mild agitation.



16



2m

Collect the resin by centrifugation at  **2500 x g, 4°C, 00:02:00**. Discard supernatant.

17

Resuspend resin in  **500 µL** of lysis buffer supplemented with  **300 millimolar (mM)** NaCl.

18

Repeat steps 16 and 17 twice.

19



2m

Collect the resin by centrifugation at  **2500 x g, 4°C, 00:02:00**. Discard supernatant.

20

Resuspend resin in  **500 µL** of  **50 millimolar (mM)** HEPES  **pH 7.5**.

21 Repeat step 16 and 20.

22



2m

Collect the resin by centrifugation at **2500 x g, 4°C, 00:02:00**. Discard supernatant.

23

Resuspend the resin in **50 millimolar (mM)** HEPES **pH 7.5** (1:1 ratio).

24

Aliquot the resin into Eppendorf tubes kept **On ice** (one Eppendorf tube per reaction, **10 µL** resin each).

In vitro kinase assay:

25

The amount of kinase and that of substrate, as well as the duration of the reaction should be optimised to ensure that the measured activity of LRRK1 is in the linear range. We recommend performing a pilot study to establish the optimal parameters to use for the kinase assay under the experimental conditions of choice.

Prepare a “master mix” containing **50 millimolar (mM)** HEPES **pH 7.5**, **50 millimolar (mM)** KCl, **10 millimolar (mM)** MgCl₂, **1 millimolar (mM)** ATP, and recombinant Rab protein.

26

Start the kinase reaction by adding the master mix to the immunoprecipitated kinase and transferring the Eppendorf tubes to the thermo mixer set at **30 °C**, **1000 rpm**.

27

Stop the kinase reaction by adding 4X LDS loading buffer to the reaction mix to a final concentration of 2X.



10m

28



Incubate the mixture at **70 °C** on a heat block for **00:10:00** to elute LRRK1 from the resin.

29

Collect the eluent by centrifugation through a 0.22-µm-pore-size Spinex column.

30

Supplement the samples with 2-Mercaptoethanol to 1% (v/v).

31



5m

Incubate the samples for **00:05:00** at **70 °C** on a heat block before proceeding to quantitative immunoblotting analysis.

Analysis of kinase reaction products by quantitative immunoblotting analysis

32

The reaction products can be analysed by quantitative immunoblotting analysis (as described in XXXX). **Table 1** lists the primary antibodies that we recommend using, which include antibodies to detect Rab7A phosphorylation at Serine-72.

A	B	C	D	E
AntibodyTarget	Company	Cat. number	Host species	Dilution
Rab7A (Total)	Sigma	R8779	Mouse	1 ug/ml
alpha-tubulin	Cell Signaling Technology	3873	Mouse	1:5,000
pS72 Rab7A	Abcam Inc.	MJF-38, Clone 1	Rabbit	1 ug/ml
LRRK1 (total) (C-terminus)	MRC-PPU Reagents and Services, University of Dundee	S405C	Sheep	1 ug/ml