




Upload image

Apr 30, 2020

# Ab Lysis Buffer

Colin Daniel<sup>1</sup><sup>1</sup>Oregon Health and Science University**1** Works for me dx.doi.org/10.17504/protocols.io.bft6jnre Colin Daniel  
Oregon Health and Science University

## MATERIALS

| NAME   | CATALOG #   | VENDOR            |
|--|-------------|-------------------|
| 100g Deoxycholate (Deoxycholic acid sodium salt) | DG090       | G-Biosciences     |
| NaCl   | 53014       | Sigma Aldrich     |
| Ethylenediaminetetraacetic acid                  | E9884       | Sigma – Aldrich   |
| Triton X-100                                     | BP151-100   | Fisher Scientific |
| Tris Base  | 604204      | Fisher Scientific |
| Sodium dodecyl sulfate solution                  | 71736-100ML | Sigma Aldrich     |

1 [M]20 Milimolar (mM) Tris pH7.5

[M]50 Milimolar (mM) NaCl

[M]0.5 % (v/v) Triton-X100

[M]0.5 % (v/v) DOC

[M]0.5 % (v/v) SDS

[M]1 Milimolar (mM) EDTA

2 On ice add protease and phosphatase inhibitors (Roche)