



Upload image

Version 4

Jun 26, 2020

Dye-terminator DNA sequencing V.4

Diep R Ganguly¹¹The Australian National University

1

Works for me

dx.doi.org/10.17504/protocols.io.bhxxkj7kw

Pogson Genomics Group



Diep Ganguly

The Australian National University



ABSTRACT

This protocol (based on the BigDye® Terminator v3.1 Cycle Sequencing Kit) is for performing terminator cycling sequencing reactions for Sanger sequencing of amplified PCR products or plasmid DNA on the 3130X genetic analyser (Applied Biosystems).

ATTACHMENTS

[BigDye Terminator v3.1.pdf](#)[wizard-sv-gel-and-pcr-clean-up-system-protocol.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bhxxkj7kw

PROTOCOL CITATION

Diep R Ganguly 2020. Dye-terminator DNA sequencing. **protocols.io**dx.doi.org/10.17504/protocols.io.bhxxkj7kw

KEYWORDS

Sanger sequencing, Dye-terminator sequencing

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 26, 2020

LAST MODIFIED

Jun 26, 2020

PROTOCOL INTEGER ID

38604

MATERIALS

NAME	CATALOG #	VENDOR
XRN-1 - 100 units	M0338L	New England Biolabs
Antarctic Phosphatase - 1,000 units	M0289S	New England Biolabs
96 well PCR Plate Non-skirted	MPS-499	Phenix Research
Wizard SV Gel and PCR Clean-Up System	A9281	Promega
Nuclease-free water (e.g. MilliQ or HPLC grade water)		
primers		
EDTA		

NAME	CATALOG #	VENDOR
10 mM dNTPs	10297-018	Life Technologies
Ethanol	100983	Merck Millipore
BigDye [®] Terminator v3.1 Cycle Sequencing Kit	4337454	Thermo Fisher
Exonuclease I (E. coli)	M0293S	NEB
Hi-Di [™] Formamide	4311320	Thermo Fisher Scientific

BEFORE STARTING

Optimize PCR cycling (if sequencing amplified PCR products) to ensure your reaction produces a single product. If needed, perform gel excision and clean-up to purify the target DNA fragment. Incubate with Antarctic phosphatase (SAP, AP, or CIP) and Exonuclease 1 to dephosphorylate and degrade unincorporated dNTPs prior to incorporating fluorescent nucleotides in the sequencing PCR (BigDye reaction).

Enzymatic PCR clean-up

- 1 If sequencing from a PCR amplified DNA fragment, gel purify target DNA band based on expected fragment size. Perform gel purification with Wizard SV Gel and PCR Clean-Up System (Promega, as per Manufacturer's instructions attached) followed by enzymatic clean-up (hydrolyze excess primers and nucleotides) with the following:

Component	Volume (μl)
Purified DNA fragment	6
Antarctic phosphatase (Genesearch)	0.5
XRN-1 (NEB)	0.25

Enzymatic clean-up of PCR products

Incubate the above in a thermal cycle for:

1. 37 °C for 15 minutes
2. 80 °C for 15 minutes.

Terminator cycling reaction

- 2 Perform sequencing PCR in PCR tubes (or 96-well plate) with BigDye Terminator cycling kit and forward or reverse^{30m} primers.

Component	Volume (μl)
v3.1 Ready reaction mix	1
5X Sequencing buffer	1.5
20 μM F/R Primer	0.5
Template (plasmid or cleaned PCR product)	50-150 ng DNA (plasmid or PCR product)
Nuclease-free water	to 10 μl

BigDye Terminator Cycling reaction

5x reaction buffer=400 mM TRIS, 10 mM MgCl₂

- 3 Run the following thermal cycling protocol: 4h
 1. 1 min at 96 °C
 2. 30-40 cycles: 96 °C for 10 seconds, 50 °C for 5 seconds, and 60 °C for 4 min.
 3. Hold at 4-12 °C.

Purification 1h 30m

- 4 Transfer PCR reaction to nuclease-free eppendorf tube. To the reaction, add 2.5 μL of 125 mM EDTA (make sure it touches bottom of tube).

- 5 Add 30 µl of 100% ethanol, *mix well* (inversion).
- 6 Incubate at room temperature for 15 minutes.
- 7 Centrifuge at 4 °C at max speed for 30 minutes.
- 8 Discard supernatant and add 50 µl of ice-cold 70% ethanol.
- 9 Centrifuge at 4 °C at max speed for 5 minutes.
- 10 Discard supernatant and allow to air-dry in the dark for >15 minutes.

Prepare for sequencing

- 11 Resuspend the pellet (likely transparent) in 7.5 µL HiDi Formamide (add to any empty wells). Incubate at RT for 5 minutes then transfer to plate. Spin down briefly.
- 12 Incubate plate at 95 °C for 3 minutes (denature) then place immediately on ice. Spin down briefly.
- 13 Submit for sequencing on 3130X genetic analyser (Applied Biosystems). Keep samples on ice.