



SEP 03, 2023

OPEN  ACCESS

Protocol Citation: José Manuel Villalobos-Escobedo 2023. Transformation Protocol ATMT Trichoderma atroviride. **protocols.io** <https://protocols.io/view/transformation-protocol-atmt-trichoderma-atroviride-b2edqba6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 28, 2021

Last Modified: Sep 03, 2023

PROTOCOL integer ID:
55461

Transformation Protocol ATMT Trichoderma atroviride

José Manuel Villalobos-Escobedo¹

¹UC Berkeley

ATMT




Maria Belen Mercado BS

ABSTRACT

This is the protocol for the transformation of conidia by Agrobacterium to perform the BarSeq technique in Trichoderma atroviride, which is reported by Villalobos-Escobedo et al., 2023.

DAY 1:

- 1 Inoculate 3 plates of *Trichoderma* @  27 °C Grow them for 4 days

Note

Sterilize 6 clean flasks
5 bottles
2 100 mL GC
2 250 ml GC

- 2 Prepare LB + Kan (500 mL)

Note

50 ug/mL is working concentration for Kan, amount added depends on stock concentration
Ex: 50 mg/mL stock concentration
1 L of LB would require 1 mL of 50 mg/ml Kan
-Usually, I put 1 ml in 500 ml of LB

DAY 3

- 3 Make liquid agro induction media (ABI)

For 100 transformations 5 mL of ABI are needed, water is autoclaved and all others are filtered and sterilized

200 mL to wash, 200 mL to resuspend (grow)

MIX AND FILTER INSIDE HOOD

Note

- CASamino acids (10% CASamino acids) (Heat up and stir)
- Glucose (36% Glucose) (Heat up and stir)
- Thiamine (10% Thiamine)

Note

Sterilized bottles

Dataset

ABI Media Spread Sheet

NAME

https://docs.google.com/spreadsheets/d/18ni9RM73H-b3UyHxGNHadtLOwJFGj_JT/edit?usp=sharing&ouid=113575406325457072750&rtpof=true&sd=true

LINK

Dataset

Chemical Mixes for ABI

NAME

https://docs.google.com/spreadsheets/d/1ezsYgXvGxfA8iqkbBXbE3URkroUd5RCiTNU_u6hufVA/edit?usp=sharing

LINK

For AB Buffer check pH, then filter

- 4 Prepare: 1 mL 3,5-dimethoxy-4-hydroxy-acetophenone (acetosyringone) *dissolved in dimethyl sulfoxide (DMSO)*. **(make this fresh every time)**


-- Acetosyringone is filtered through syringe and 0.22 um filter top in DMSO (0.098 g/ mL)

- 5
1. Autoclave 437ml ddH₂O with 20g Bacto Agar
 2. Add MIXED induction liquid to agar.
 3. Mix together inside hood for a few minutes
 4. Pour plates.

2d

Pour solid agro induction media plates and leave on bench to dry for 2 days.

- 6 Prepare medium for Agrobacterium LB + Kan.

7 Inoculate *A. tumefaciens* strain harboring the plasmid for Trichoderma transformation onto LB + Kan.  250 µL in 50 mL

8 Shake at 250rpm at 30C overnight

(Place around 7 pm- 8 pm)

Day 4

9 Measure the OD600 of the agro cell culture. (It will likely be around OD600 1.)

(Around 9 am)

9.1 In order to get an accurate measurement, dilute the sample 1:10 in H₂O in the cuvette.

(100 uL agro cell culture, 900 uL H₂O)

10 Back dilute the agro culture to OD600 0.5 in LB+Kan with twice the number of ml of media to the number of transformations you plan to do (i.e., if you are doing 100 transformations back dilute into 200ml LB).

100 mL from stock and 100 mL new LB+Kan


11 Shake at 250rpm at 30C for 2h

2h

12 Measure the OD600 of the agro cell culture after 2 hrs

Note

It should be around OD600 1. If it is not at least OD600 0.95, put it back on the 30C shaker for another 30min.

13 Centrifuge the agro culture at .  4000 rpm, 00:15:00

15m

14 Pour supernatant into waste container.


15 Wash cells in agro induction media

Note

The number of mL should correspond to twice the number of transformations that you are doing

ex: 20 ml for 10 transformations

New flasks are used

16 Centrifuge the agro cells at  4000 rpm, 00:15:00

15m

17 Pour supernatant into waste container.

18 Resuspend cells in agro induction media

Note

The number of mL should correspond to twice the number of transformations that you are doing

For 100 transformations 200 ml of induction media, and each transformation has 4×10^7 conidia,

it's 2×10^7 per ml. Concentration of agro.

We can use tubes with caps for roller.

18.1 Put cultures on shaker @ 25C and rpm 250, for 24 hrs

18.2

Note

Materials for next day:

- Filter top
- 0.45 um Filter
- Vacuum

DAY 5

6d 0h 5m

19

Note

(using plates of *T. atroviride* that are 4 days old)

NOTE: Make sure the big square plates are at room temperature.

Approximately 23h after putting agro cells into induction media, collect *Trichoderma* conidia from plates into 5-10ml of H₂O.

20 Shake and vortex vigorously to suspend conidia


21 Measure conidial concentration

Count four corner squares, and middle square

Avg spores: Sum of both sides / 10 squares

Math: (Avg. spores)(1x10⁴)(25)(100)

22 Aliquot 2e⁷ conidia into 2ml Eppendorf tubes

23 Spin cells at  10000 rpm, 00:05:00
or 4000 rpm for 20 min


5m

24 Pipet off supernatant

25 Resuspend cells (pellet) in 2ml of agro induction media cell culture (from yesterday)

26 Incubate at room temperature for at least 5 min

27 Put a 0.2um sterile PES bottle top filter onto a bottle and set up the vacuum



Equipment

Thermo Scientific™ Nalgene™ Rapid-Flow™ Sterile Disposable Bottle Top Filters with PES Membrane

NAME

0.2 µm Bottle Top Filter or equivalent

TYPE

Nalgene

BRAND

09-741-07

SKU


<https://www.fishersci.com/shop/products/nalgene-rapid-flow-sterile-disposable-bottle-top-filters-pes-membrane/0974107>

LINK

500 mL, 0.2 µm PES Bottle Top Filter

SPECIFICATIONS

- 28 Place a 0.45µm HAWP filter onto the bottle top filter using sterile forceps



Equipment

HAWP MF-Millipore Membrane Filter, 0.45 µm pore size

NAME

Membrane filter

TYPE

Millipore

BRAND

HAWP04700

SKU

[https://www.sigmaaldrich.com/catalog/product/mm/hawp04700?](https://www.sigmaaldrich.com/catalog/product/mm/hawp04700?lang=en®ion=US&gclid=CjwKCAjw8pH3BRAXEiwA1pvMsdoaQbbYstapLy8iGgQMUPbpUlubisFSK9v3zg7Ab-Uv1HEHZmOhSBoCPx8QAvD_BwE)

LINK

[lang=en®ion=US&gclid=CjwKCAjw8pH3BRAXEiwA1pvMsdoaQbbYstapLy8iGgQMUPbpUlubisFSK9v3zg7Ab-Uv1HEHZmOhSBoCPx8QAvD_BwE](https://www.sigmaaldrich.com/catalog/product/mm/hawp04700?lang=en®ion=US&gclid=CjwKCAjw8pH3BRAXEiwA1pvMsdoaQbbYstapLy8iGgQMUPbpUlubisFSK9v3zg7Ab-Uv1HEHZmOhSBoCPx8QAvD_BwE)

0.45 µm 47 mm

SPECIFICATIONS

29 Pipet the 2ml of Trichoderma/agro cell mixture onto 0.45µm HAWP filter

Note


Being careful to only get cell mixture on the HAWP filter.

30 Wait for the liquid induction media (ABI) to filter through the bottle top filter

31 Using sterile forceps, carefully transfer the filter with cells on it to the solid induction media petri dish, making sure to keep the side with the cells on it facing up!

32 Incubate induction media plates with filters at 20C in the dark for 6 days.

6d



Dried filters ready for incubation

32.1

Note

Materials for next part:

- 50 mL or 500 mL tubes with sterile beads and sterile water
- Layered PDA plates
- MM plates for quantifying
- Sterile water

Day 11

- 33 Label 2 falcon tubes for each HAWP filter

Note

25 - 50 filters fit into 500 mL tubes

- 34 Add 4ml sterile distilled H₂O to a 50ml conical tube

Note

500 mL conical tube

- 35 Using sterile forceps, carefully transfer the HAWP filter with cells on it to the  50 mL conical

36

Note

Use beads and VORTEX until it has all fallen off



- 37 Vortex vigorously to get the rest of the cells off of the HAWP filter

38 Transfer 2ml of cell suspension to sterile  2 mL Eppendorf tube.


38.1 Spin at  10000 rpm, 00:05:00 .

5m

39 Pipet off supernatant into a waste bottle

40 Transfer the rest of the cell suspension to the  2 mL Eppendorf tube and Spin at  10000 rpm, 00:05:00 .
Pipet off supernatant into a waste bottle

5m

41 Resuspend cell mixture into  100 μ L sterile distilled H₂O

42 To spread in large plates (150 mm) PDA.

- 1.- PDA with hyg and Carb (4 ml Hyg and 1 ml Carb per L): dispense 50 ml.
- 2.- PDA with Carb (1 ml per Liter): dispense 13.5 ml.
- 3.- Conidia transformants: 1 ml
- 4.- Low-melt PDA with Carb (1 ml per L). 13.5 ml. 37 C. we can mix with conidia.
- 5.- PDA layer carb and hyg. 13.5 ml, kept at 55 C.

42.1 1 plate with filter goes to 3 plates of top/bottom media

42.2
Plate cells on

5mL Top: MM + Carb

25 mL bottom: MM + Carb + Hyg

You don't need to dilute the control plate but for the experimental ones:

1. 100%:
2. 1:25: **4 ul of cells in 96 ul of water**

Note

Make up media in half the amount of water, and autoclave water separately, then combine


43 Place at 25C in the light for 7 days (checking efficiency with Triton)

43.1

Final Harvest (Day 18)

20m

44




Washing of plates with water

Note

Materials needed:


- 50% glycerol
- Sterile water
- Glass spreaders
- 1 ml Pipettes
- 25 ml pipettors and pipet tip
- Tips 100-1000 ul
- 15 mL Falcon tubes

Dispense  10 mL of sterile water onto the plate

45 Use the sterile glass spreader to detach as many as possible

46 Transfer into a  50 mL falcon tubes


46.1 Add another  10 mL of sterile water to plate and transfer into tube

47 Spin down  50 mL falcon tubes at 4,000 RPM for 20 minutes

20m

48 Resuspend pellets in 10 ml of 50% glycerol

48.1

Store in  1.5 mL tubes that are no more than half full, stored at an angle to avoid rupturing the tube.

Note

Where is it stored? What temperature?