

Sep 10, 2020

Mycoplasma

In 1 collection

Andrea Argouarch¹¹University of California, San Francisco
1 Works for me dx.doi.org/10.17504/protocols.io.8gphtvn

Neurodegeneration Method Development Community

Tech. support email: ndcn-help@chanzuckerberg.com

Andrea Argouarch

University of California, San Francisco

ABSTRACT

Detection of mycoplasma utilizing the bulldog-bio kit. Steps include, gDNA isolation from cells, nanodrop, PCR, and running an agarose gel to detect mycoplasma within a cell culture line.

Image from <https://cellculturedish.com/cell-culture-basics-mycoplasma-101-a-practical-guide-to-prevention-detection-and-elimination-of-mycoplasma-contamination/>

DOI

dx.doi.org/10.17504/protocols.io.8gphtvn

PROTOCOL CITATION

Andrea Argouarch 2020. Mycoplasma. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.8gphtvn>

COLLECTIONS ⓘ

Dural Cell Isolation and Culturing - Collection

KEYWORDS

dura mater, dural cells, mycoplasma, bulldog

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 19, 2019

LAST MODIFIED

Sep 10, 2020

PROTOCOL INTEGER ID

28911

PARENT PROTOCOLS

Part of collection

[Dural Cell Isolation and Culturing - Collection](#)

STEPS MATERIALS

NAME

CATALOG #

VENDOR

NAME	CATALOG #	VENDOR
Epicentre QuickExtract™ DNA Extraction Solution	QE09050	Epicentre
Trypsin-EDTA (0.05%), phenol red	25300062	Thermo Fisher
DPBS, no calcium, no magnesium	14190250	Thermo Fisher
e-Myco PLUS Mycoplasma PCR Detection Kit	25234	
Agarose	PBSA1705	
Ethidium Bromide (1% Solution/Molecular Biology)	BP130210	Fisher Scientific
Gel Loading Dye Blue (6X) - 4.0 ml	B7021S	New England Biolabs
GeneRuler 1 kb Plus DNA Ladder, ready-to-use	SM1333	Thermo Fisher

Observations

- 1 Culture should be in at least 3-6 days of subculture after split
- 2 During isolation and PCR, use clean gloves, spray down bench, and keep reagents and tubes cold
🧊 On ice

Isolation of gDNA

- 3 Isolate cells once confluent, at least 3-6 days after split
- 4 After trypsinization for freezing cells, rinse remaining cells in both flasks with 12 mls 🧴 12 mL PBS



Trypsin-EDTA (0.05%), phenol red

by Thermo Fisher

Catalog #: 25300062



DPBS, no calcium, no magnesium

by Thermo Fisher

Catalog #: 14190250

- 5 Spin for 5 mins 🕒 00:05:00 at 100 rpm 🌀 1000 rpm at 4C 🌡️ 4 °C
- 6 Aspirate supernate
- 7 Place in ice bucket 🧊 On ice on bench or store in -80C 🌡️ -80 °C

- 8 Add 40 µl  **40 µl** of Quick Extract to pellet













**Epicentre QuickExtract™ DNA
Extraction Solution**


by Epicentre

Catalog #: **QE09050**

- 9 Pipette up and down and transfer to PCR tubes
a. Label PCR tube with line, date, and P for pellet


- 10 Extract gDNA with thermocycler program
- Turn on Power button on the back of machine
 - Open lid with handle
 - Put PCR tubes in small slot and close lid
 - Select File, click Enter, select pcr chamber and start program
 - 65C  **65 °C** for 15 mins  **00:15:00**
 - 98C  **98 °C** for 10 mins  **00:15:00**
 - 4C  **4 °C** to hold
 - Run time around ~30mins  **00:30:00**
 - Click enter to exit after run
 - Turn off machine and place on ice  **On ice**





- 11 Centrifuge PCR tubes for 5 mins  **00:05:00** at max  **13.4 rpm** and transfer SN to new 1.5 mls tubes labeled gDNA, line, and date  **On ice**
- Tiny pellet ok left in PCR tube

- 12 Nanodrop to determine gDNA (pre-PCR) concentration on ice  **On ice**


- 13 Store at -20C  **-20 °C** or run PCR

Nanodrop

- 14 Items to bring to nanodrop machine
- p2 pipette, p2 tips, pen, kimwipes, ice bucket
 - Keep samples on ice  **On ice**
- 15 Open the NanoDrop 1000 software and select nucleic acid from the main menu





- 16 Clean the nanodrop pedestal with a kimwipe
 - 17 Add 2 μ l  2 μ l of nuclease free water to the center of the pedestal and gently close the lever arm
 - 18 Click 'OK' to initialize the spectrophotometer
 - 19 Select DNA from the dropdown menu
 - 20 Clean the nanodrop pedestal with a kimwipe, add 2 μ l  2 μ l of a blank sample, and click 'Blank'
 - a. Use quick extract solution as blank
 - 21 Clean the nanodrop pedestal with a kimwipe, add 2 μ l  2 μ l of sample and click 'Measure'
 - 22 Record the concentration (ng/ μ l) of the sample on the tube and a table
 - 23 Record the 260/280 ratio, and the 260/230 ratio
-  260/280 ratio should be above 1.8 (salt contamination)
260/230 ratio should be 1-2
- 24 Clean the nanodrop pedestal with a kimwipe in-between samples
 - 25 Click exit
 - 26 Leave kimwipe between the pedestal and the lever arm when finished

PCR Reaction


27 Make a 20 µl  20 µl of PCR reaction in provided eMyco Plasma PCR Tube  On ice



e-Myco PLUS Mycoplasma PCR
Detection Kit
Catalog #: 25234

- a. QS to 20 µl  20 µl with nuclease free water
- b. Add 40 ng  40 ng of sample
 - i. Can make working stock if concentration is too high
- c. Include a negative control with just water (20 µl)  20 µl
- d. Include a positive control with provided positive control (1 µl)  1 µl

28 Run PCR program

- a. Turn on machine at the back
- b. Tap to mix, place samples into thermocycler in small tube holes, and close lid
- c. Run myco program
 - Initial Denature** = 94C  94 °C for 1 min  00:01:00
 - ii. 35 Cycles
 - Denature** = 94C  94 °C for 30 secs  00:00:30
 - Anneal** = 58C  58 °C for 20 secs  00:00:20
 - Extend** = 72C  72 °C for 1 min  00:01:00
 - Final Extension** = 72C  72 °C for 5 mins  00:05:00
 - Hold** at 4C  4 °C
- d. Press enter to run program
- g. Run time is around ~1.5hrs  01:30:00

29 Can prep for agarose gel near the end of the PCR run

30 Store gDNA samples, at -20C  -20 °C

31 After run, select enter, turn off machine, and place PCR product on ice  On ice

32 Run gel or store at -20C  -20 °C

Agarose Gel

- 33 Wet rubber around gel chamber to help it slide into gel rig
- 34 Assemble with rubber on the black edges of the gel rig and place purple comb in the slot on the right
- 35 Make 2% Agarose in small beaker with screw cap for a mini rig size.



Agarose
Catalog #: PBSA1705

a. 1.2 g **1.2 g** of agarose with 60 mls **60 mL** of 1X TAE (Stock is 20X TAE)

- 36 Microwave for 2 mins **00:02:00** with a loose cap and swirl every 30 seconds **00:00:30** until dissolved
- 37 When temperature is cool to the touch, add 1:10,000 of Etbr



Ethidium Bromide (1%
Solution/Molecular Biology)
by Fisher Scientific
Catalog #: BP130210



Mutagen - may potentially cause carcinogenic or teratogenic effects

ie add 6 μ l **6 μ l** Etbr per 60 mls **60 mL** of agarose

- 38 Swirl to mix and add to chamber to set, allow for 10 mins **00:10:00** to cool in cold room **4 °C**
- 39 Prep samples for gel in PCR tubes on ice **On ice**
- a. 5 μ l **5 μ l** of PCR product
 - b. 3 μ l **3 μ l** of Nuclease free water
 - c. 2.6 μ l **2.6 μ l** of 6X loading buffer (2X final concentration from a 6X stock, purple)




Gel Loading Dye Blue (6X) - 4.0 ml

by New England Biolabs

Catalog #: B7021S

$$i.(6X)(a) = (2X)(8\mu l); a = 2.6\mu l$$

- 40 Rotate gel chamber when set, so the lanes are at the top of the chamber
- 41 Fill with 1X TAE until the top of the white knob in gel rig and remove comb slowly
- 42 Load 4 μ l  4 μ l of 1kb O'gene plus ladder into first lane



GeneRuler 1 kb Plus DNA Ladder, ready-to-use

by Thermo Fisher

Catalog #: SM1333

- 43 Load 11 μ l  11 μ l of each sample per well
- 44 Add lid, black at the top, red at the bottom, both at the right side
- 45 Run at 130V for 50 mins  00:50:00
- 46 Store PCR product at -20C  -20 °C
- 47 Dispose Etbr waste in proper container (liquid and solid)