



Jan 14, 2022

Nuclear extraction from endometrial tumors for single nuclei sequencing

gracefoley¹¹Northwestern Universitydx.doi.org/10.17504/protocols.io.b3vfqn3n

Northwestern University



gracefoley

Modified from: Slyper, M. *et al.* A single-cell and single-nucleus RNA-Seq toolbox for fresh and frozen human tumors. *Nat. Med.* 2020 26:792–802 (2020).

A optimized protocol for nuclear extraction from endometrial tumors. Performed with endometrial adenocarcinoma, endometrioid type, FIGO grade 1.

[Nuclear extraction optimization for single nuclei RNA-seq.pptx](#)

DOI

dx.doi.org/10.17504/protocols.io.b3vfqn3n<https://www.nature.com/articles/s41591-020-0844-1#Sec13>

gracefoley 2022. Nuclear extraction from endometrial tumors for single nuclei sequencing . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b3vfqn3n>



protocol ,

Jan 14, 2022

Jan 14, 2022

56967

Buffers:

2X ST Buffer

	Stock	Final	Volume for 10 mL
NaCl	5 M	292 mM	584 uL
Tris-HCl (pH 7.5)	1 M	20 mM	200 uL
CaCl ₂	1M	2 mM	20 uL
MgCl ₂	100 mM	42 mM	4.2 mL
Ultrapure Water	-	-	5 mL

NST Buffer

	Stock	Final	Volume for 10 mL
Nonidet P40 Substitute	10%		200 ul
BSA	5%		20 ul
ST Buffer	2X	1X	5 ml
Nuclease Free Water	-	-	4.75 mL

Nuclei Resuspension Buffer

	Stock concentration	Final	Volume Needed for 50 mL
BSA solution	5%	1%	10 ml
RNase Inhibitor	40 U/ul	0.2 U/ul	0.25 ml
1X PBS	-	-	39.75 ml

- 1 Place frozen tissue in a well of a 6 well plate containing 1 mL of cold NST disassociation buffer.

1.1 Ensure plate is on ice

- 2 Chop tissue using a scalpel until mostly homogenized

2m

- 3 Incubate solution on ice for 5 minutes 5m
- 4 Filter the homogenized solution through a 40 uM falcon cell strainer into a 50 mL falcon tube
- 5 Use an additional 1 mL of NST buffer to rinse the well and filter
- 6 Count nuclei
- 7 Centrifuge nuclei at 300 g for 10 minutes and aspirate supernatant completely 10m
- 8 Resuspend pellet and bring the volume up to 5 mL using ST buffer
- 9 Transfer solution to 15 mL conical tube and centrifuge for 10 minutes at 300 g at 4oC 10m
- 10 10. Resuspend pellets on ice in 1 ml of the ST buffer
- 11 11. Filter through 35 um falcon strainer
- 12 12. Count nuclei

13 Centrifuge nuclei at 300 g for 10 minutes and aspirate supernatant completely 10m

14 13. Resuspend nuclei pellet in resuspension buffer to get a concentration of 1,000 nuclei per uL

14.1 Start with .5 mL per 0.2 grams of endometrial tumor tissue

14.2 Look at nuclear membranes at 40x-60x to ensure high nuclear membrane quality with minimal blebbing