

Electron
Microscope
Sample
Preparation
Technique

Aug 02, 2022

Electron microscope sample preparation technique_V2

Forked from [Electron microscope sample preparation technique](#)Wai Kit Lam (Leo)¹¹Laboratory of Michael Lazarou (Team Hurley) at Monash University

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.14egn75zmv5d/v1 Wai Kit Lam (Leo)

ABSTRACT

There are eight main steps in preparing EM sample, including fixation, en bloc staining, agarose embedding, dehydration, infiltration, resin embedding, resin-embedded sample trimming and ultrathin sectioning.

This EM sample preparation protocol is established by Benjamin Padman and slightly modified by Wai Kit Lam and Runa Lindblom. Of note, this protocol is designed for processing cultured cell sample.

ATTACHMENTS

[_ASAP_Team_Hurley_Lazarou_lab_Leo_Lam_Electron_microscope_sample_preparation_technique.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.14egn75zmv5d/v1

PROTOCOL CITATION

Wai Kit Lam (Leo) 2022. Electron microscope sample preparation technique_V2.
protocols.io
<https://dx.doi.org/10.17504/protocols.io.14egn75zmv5d/v1>

FORK NOTE



FORK FROM

Forked from [Electron microscope sample preparation technique](#), Dorotea Fracchiolla

KEYWORDS

ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 02, 2022

LAST MODIFIED

Aug 02, 2022

PROTOCOL INTEGER ID

68051

MATERIALS TEXT



List of reagents

- 0.2 M Phosphate buffer (pH 7.4)
- 0.2 M sodium-cacodylate (NaCac) (pH 7.4)
- 16 % Paraformaldehyde solution, EM grade
- 25 % Glutaraldehyde solution, EM grade
- 2 % osmium tetroxide (OsO₄)
- Potassium ferricyanide (K₃Fe(III)(CN)₆)
- Uranyl acetate (UA)
- Lead citrate
- Agarose, low-gelling point
- 70 %, 80 %, 90 %, 95 %, 100 % Ethanol
- Propylene oxide
- Dodecenylsuccinic anhydride (DDSA)
- Araldite 502
- Procure 812
- Benzyldimethylamine (BDMA)

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

Fixation 1h

- 1 Aspirate cell culture media then fix cells with pre-warmed phosphate-buffered 4% Paraformaldehyde (PFA) at  **Room temperature** for  **01:00:00** .

1h

- 2  

Wash cells twice with **0.1 Molarity (M)** sodium-cacodylate (NaCac) **pH 7.4** buffer.

- 3 Post-fix the sample with **2.5 Mass Percent** glutaraldehyde in **0.1 Molarity (M)** NaCac buffer for **24:00:00** at **4 °C**. 1d



Wash cells three times with **0.1 Molarity (M)** NaCac buffer.

- 5 Post-fix the sample with **1 Mass Percent** osmium tetroxide (OsO₄), **1.5 Mass Percent** Potassium ferricyanide (K₃Fe(III)(CN)₆) in **0.1 Molarity (M)** NaCac buffer **On ice** for **01:00:00**. 1h

- 6 Bring the sample to BioWave Pro microwave (Pelco), and expose it to three microwave dutycycle (**00:02:00** on, **00:02:00** off) at 100W under vacuum. 4m

- 7 Rinse cells three times with MilliQ water.

en bloc staining

- 8 Add **2 Mass Percent** uranyl acetate (UA) in MilliQ water to sample, followed by exposing the sample to three microwave duty-cycle (120s on, 120 s off) at 100W under vacuum.

It is best to filter the 2% (w/v) UA using 0.22 µm pore size filter.

- 9 Rinse cells three times with MilliQ water.

Agarose embedding




10 

Scrape and pellet ( **10000 x g, 00:03:00**) cells in MilliQ water.

11 Resuspend the cell pellet in  **70 % volume** Ethanol.

12 

Centrifuge the sample again at  **10000 x g, 00:03:00** .

13 Gently remove the supernatant (i.e. 70 % Ethanol).

14 Add one drop of low-melting point agarose to the cell pellet.

15 

Vortex and centrifuge the sample immediately.


The low-melting point agarose must be super-hot (close to 100°C) to achieve step 15.


16 Cut the solidified agarose-cell pellet into 1 mm cubes using a razor blade.

17 Transfer the agarose-cell pellet cubes to a flat-bottom crew cap tube.


Dehydration

40s

18 Perform a microwave assisted dehydration by graduated ethanol series (80 %, 90 %, 95 %, ^{40s}100 %, 100 % (v/v); each at 150 W for  **00:00:40**).

- 19 Perform a microwave assisted dehydration by propylene oxide (100%, 100% (v/v); each at 150 W for  00:00:40).

Infiltration 3m

- 20 Infiltrate samples with Araldite 502/Procure 812 (Resin) by graduated concentration series in propylene oxide (25 %, 50 %, 75 % 1 100 %, 100 % (v/v)); each at 250 W for  00:03:00 under vacuum).

The composition of the Resin is 51.7 % (w/w) DDSA, 18.6 % (w/w) Araldite 502, 26.3 % (w/w) Procure 812 and 3.4 % (w/w) BDMA.

Resin embedding

- 21 Put the infiltrated agarose-cell pellet cubes and the appropriate paper label into the resin-block mould, then fill up the mould with resin.

It is important to remove all the bubbles. Usually, some bubbles are trapped under the paper label. It is best to use a pair of forceps to lift-up the paper label to release these bubbles, followed by using a 1 ml plastic transfer pipette to remove them.

- 22  

2d

Incubate the resin-embedded sample to a  60 °C oven for  48:00:00 or until the resin is fully polymerised.

Resin-embedded sample trimming

- 23 Use a razor blade to carefully trim the resin block to expose the underlying agarose-cell pellet cubes.
- 24 Create a mirror surface on the agarose-cell pellet cubes using Ultra UCT ultramicrotome (Leica Biosystems) equipped with a glass knife.

Before proceeding to the ultrathin sectioning, the block face must be cleaned, and the edges on the mirror surface must be parallel to each other.

Ultrathin sectioning

8m

- 25 Cut 70-90 nm sections on an ultramicrotome equipped with a diamond knife (Ultra 45° Diatome).

Make sure you check the angle of the knife holder. For the Ultra 45° Diatome, the angle needs to be at 6°.

- 26 Load the ultrathin sections onto the grid. (Copper side up)

- 27 The grid is ready for TEM imaging.