




May 28, 2020

Immunofluorescence staining of paraffin embedded cell block sections

Oncoimmunology

Erez Nissim Baruch¹, Rona Ortenberg²¹Department of Clinical Immunology and Microbiology, Sackler Faculty of Medicine, Tel Aviv University, Israel, The Ella Lemelbaum Institute for Immuno-Oncology, Sheba Medical Center, Tel-HaShomer, Israel,²The Ella Lemelbaum Institute for Immuno-Oncology, Sheba Medical Center, Tel-HaShomer, Israel1 Works for me dx.doi.org/10.17504/protocols.io.6pqhdmw
 Rona Ortenberg

ABSTRACT

This protocol is designated for Immunofluorescence staining of paraffin embedded cell block sections.

EXTERNAL LINK

<https://doi.org/10.1080/2162402X.2020.1741267>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Baruch EN, Ortenberg R, Avivi C, Anafi L, Dick-Necula D, Stossel C, Moshkovits Y, Itzhaki O, Besser MJ, Schachter J, Barshack I, Markel G, Immune co-culture cell microarray – a feasible tool for high-throughput functional investigation of lymphocyte–cancer interactions. Oncoimmunology 9(1). doi: [10.1080/2162402X.2020.1741267](https://doi.org/10.1080/2162402X.2020.1741267)

MATERIALS

NAME	CATALOG #	VENDOR
Ethanol 100%		
Tween 20	P1379-500ml	Sigma-aldrich
Xylene	XC9800.SIZE.1L	Bio Basic Inc.
Fluoromount	F4680	Sigma
DAPI	D3571	Invitrogen - Thermo Fisher
DPBS no calcium no magnesium	02-023-1A	
10x Citrate Buffer pH6.0	K035	

- 1
 1. After the slides were cut, let them dry in room temperature over-night, or 2-3 hours on a heating plate at 37°C.
 2. Before staining incubate the slides in oven at 55°C for 1 - 2 hours
- 2 In a chemical hood:
 1. Wash in Xylene for 10 minutes
 2. Wash in new Xylene for another 10 minutes
 3. Wash in Ethanol 100% for 3 minutes
- 3 On the bench, room temperature:
 1. Wash in new Ethanol 100% for 3 minutes.
 2. Wash in Ethanol 95% for 3 minutes.
 3. Wash in Ethanol 70% for 3 minutes.
 4. Wash in 0.05% Tween-20 in PBS for 5 minutes.
 5. Wash again in new 0.05% Tween-20 in PBS for 5 minutes.



- 4 Microwave:
1. Put the citrate buffer x 1 in a plastic container designated for microwave use.
 2. Put the slides in the container.
 3. Insert the container into the microwave and activate on maximal energy for 1 minute and on minimal energy for 10 minutes.
 4. Stop heating once the buffer boils.
 5. Take the container out of the microwave, and let it cool down for 15-30 minutes to room temperature. **Do not** transfer the slide directly from the boiled buffer to other fluids at room temperature. Wait until proper cool down.

- 5 Bench, room temperature:
1. Wash in a new 0.05% Tween-20 in PBS for 5 minutes.
 2. Wash again in a new 0.05% Tween-20 in PBS for 5 minutes.
 3. Take the slide out of the Tween PBS container, and put in the slide staining box.
 4. Prepare blocking solution - 3% normal serum of goat/sheep/donkey/other in PBS (make sure it's a different species from the secondary antibody).
 5. Gently drop 250-300 ul the blocking solution on the slide.
 6. Cover the slide with a paraffin film and incubate for 60 minutes.
 7. Remove the paraffin film from the slide. Gently tap the slide on a paper towel in order to dry it.



- 6 Bench, 4 °C:
1. Add 150-200 ul of a primary antibody diluted in a blocking solution
 2. Cover with a new paraffin film.
 3. Incubate for 60 minutes at room temperature or 4 °C.
 4. At 4°C, the antibody can stay on the slide over-night.

- 7 Bench, room temperature:
1. Remove the paraffin film.
 2. Wash in 0.05% Tween-20 in PBS for 5 minutes.
 3. Wash again in a new 0.05% Tween-20 in PBS for 5 minutes.
 4. Wash for a third time in new 0.05% Tween-20 in PBS for 5 minutes.
 5. Add 150-200 ul of Secondary Ab diluted in a blocking solution
 6. Cover with paraffin film and incubate 60 minutes in room temperature.
 7. Wash in 0.05% Tween-20 in PBS for 5 minutes.
 8. Wash again in a new 0.05% Tween-20 in PBS for 5 minutes.
 9. Wash for a third time in a new 0.05% Tween-20 in PBS for 5 minutes.
 10. Wash with DDW for 2 minutes

- 8 Stain with DAPI (1:2000 in DDW) for 2 minutes.
1. Wash with DDW for 5 minutes.
 2. Wash in a new 0.05% Tween-20 in PBS for 5 minutes.
 3. Gently tap the slide on a paper towel in order to dry it

- 9 Cover slides
1. Add a few (2-3) drops of a water-based glue
 2. Cover slide with cover-slip
 3. Cover slide with an aluminium foil and let them dry over night



- 10 Proceed to microscopy analysis