



SEP 28, 2023

OPEN  ACCESS**DOI:**[dx.doi.org/10.17504/protocols.io.rm7vzx6e2gx1/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzx6e2gx1/v1)

**Protocol Citation:** melissa.kalladeen 2023. Conventional methods for isolation and Identification of Streptococcus pyogenes.

**protocols.io**<https://dx.doi.org/10.17504/protocols.io.rm7vzx6e2gx1/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Sep 28, 2023

## Conventional methods for isolation and Identification of Streptococcus pyogenes

melissa.kalladeen<sup>1</sup><sup>1</sup>University of The West Indies, St. Augustine Campus

UG and UWI Collab



Paul Cheddie

**ABSTRACT**

The protocol outlines the steps for isolating and identifying Group A streptococci (Streptococcus pyogenes) from growth on primary culture media.

## Culture and identification of *Streptococcus pyogenes*

- 1 *S. pyogenes* isolates that were presumptively identified were collected from the different hospitals and outpatient centres that formed part of the study. Isolates were inoculated to 5% sheep blood agar plates that were incubated overnight at 35°C.
- 2 Upon receipt, a beta-haemolytic colony was streaked onto a new 5% blood agar (HiMedia Laboratories Limited, Mumbai, In.) plate to obtain isolated colonies and incubated at 35-37°C for 24 hours in an anaerobic jar.
  - 2.1 A bacitracin antibiotic disc was placed into the inoculum of each isolate. Any zone of inhibition to the antibiotic after incubation was considered to be an indication of drug susceptibility.
  - 2.2 For quality control purposes, *S. pyogenes* ATCC 19615 (positive control), *S. agalactiae* ATCC 12386 (negative control), *Enterococcus faecalis* ATCC 29212 (negative control), and an uninoculated blood agar plate were incubated alongside the test isolates.
- 3 Following incubation, the pyrrolidonylarylamidase (PYR) test was performed on all beta-haemolytic isolates based on the manufacturer's instructions.
  - 3.1 A PYR disc (HiMedia Laboratories Limited, Mumbai, In.) was placed on a glass slide on a petri dish and rehydrated with one drop of sterile deionized water.
  - 3.2 Using a sterile loop, 2-3 colonies were added to the suspension on the glass slide and the mixture incubated at room temperature for 2 minutes.

- 3.3** One drop of the PYR chromogenic solution was placed on the disc and observations recorded after 3 minutes. The formation of a red colour was considered a positive result
- 4** Beta-haemolytic colonies from the test samples and controls were also exposed to the Voges-Proskauer (VP) Test (BD, New Jersey, USA).
- 4.1** A colony was subcultured to tubes containing 4 mL of trypticase soy broth (HiMedia Laboratories Limited, Mumbai, In.) and incubated at 35-37°C for 24 hours.
- 4.2** Following incubation, 0.6 mL of alpha-naphthol solution and 0.2 mL of potassium hydroxide solution were added to the TSB tubes and mixed by vortexing for 10 sections.
- 4.3** Observations were recorded within 5 minutes. QC observations were made using *E. faecalis* ATCC 29212 (positive) and *S. pyogenes* ATCC 19615. Tubes that showed a yellow-brown colour were considered VP negative and indicative of *S. pyogenes*.
- 5** Confirmatory testing was performed using the Lancefield classification method with a Strep Grouping kit (Thermo Scientific Remel PathoDX, Remel Europe LTD, Kent, UK) following the manufacturer's instructions.