



Version 2

Nov 30, 2020

# Quick Protocol for Monarch® DNA Gel Extraction Kit (NEB #T1020) V.2

New England Biolabs<sup>1</sup><sup>1</sup>New England Biolabs

1

Works for me

[dx.doi.org/10.17504/protocols.io.bp9gmr3w](https://dx.doi.org/10.17504/protocols.io.bp9gmr3w)

New England Biolabs (NEB)

Tech. support phone: +1(800)632-7799 email: [info@neb.com](mailto:info@neb.com)Danielle Freedman  
NEB

## ABSTRACT

This is the quick version of the Monarch® DNA Gel Extraction Kit Protocol (NEB #T1020). For the full protocol, please click [here](#).

## EXTERNAL LINK

<https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-dna-gel-extraction-kit-t1020>

## DOI

[dx.doi.org/10.17504/protocols.io.bp9gmr3w](https://dx.doi.org/10.17504/protocols.io.bp9gmr3w)

## EXTERNAL LINK

<https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-dna-gel-extraction-kit-t1020>

## PROTOCOL CITATION

New England Biolabs 2020. Quick Protocol for Monarch® DNA Gel Extraction Kit (NEB #T1020).

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bp9gmr3w>

Version created by Lenny Teytelman



## EXTERNAL LINK

<https://www.neb.com/protocols/2015/12/08/quick-protocol-for-monarch-dna-gel-extraction-kit-t1020>

## KEYWORDS

Extracting, Monarch DNA Gel Extracting Kit, extracting dna

## LICENSE

— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 30, 2020

## LAST MODIFIED

Nov 30, 2020

## OWNERSHIP HISTORY

Nov 30, 2020



Danielle Freedman

NEB

PROTOCOL INTEGER ID

45064

GUIDELINES

For detailed protocol and more information, visit [www.neb.com/T1020](http://www.neb.com/T1020)

The full protocol is available [here](#).

A video protocol is available [here](#).

MATERIALS TEXT

MATERIALS

 **Monarch® DNA Gel Extraction Kit New England**

**Biolabs Catalog #T1020**

ABSTRACT

This is the quick version of the Monarch® DNA Gel Extraction Kit Protocol (NEB #T1020). For the full protocol, please click [here](#).

BEFORE STARTING

Please review the important information under the “Guidelines” tab before beginning.

Add 4 volumes of ethanol ( $\geq 95\%$ ) to one volume of DNA Wash Buffer.

- For 50-prep kit, add 20 ml of ethanol to 5 ml of Monarch DNA Wash Buffer
- For 250-prep kit, add 100 ml of ethanol to 25 ml of Monarch DNA Wash Buffer

All centrifugation steps should be carried out at 16,000 x g (~13,000 RPM).

Please note: column holds 800  $\mu$ l

- 1 Excise the DNA fragment from the agarose gel, taking care to trim excess agarose. Transfer to a 1.5 ml microfuge tube and weigh the gel slice.** Minimize exposure to UV light.

Trim/remove excess agarose and minimize exposure to UV light.

- 2 Add 4 volumes of Gel Dissolving Buffer to the gel slice** (e.g., 400  $\mu$ l buffer per 100  $\mu$ l or 100 mg agarose).

- 3 Incubate the sample between 37–55°C (typically 50°C), until the gel slice is completely dissolved** (generally 5–10 minutes). The time that takes a gel slice to melt depends on the size of the slice, the temperature used in the incubation as well as the percent agarose used in the gel. The time recommended above should be used just as a guideline.

For DNA fragments > 8 kb, an additional 1.5 volumes of water should be added after the slice is dissolved to mitigate the tighter binding of larger pieces of DNA (e.g., 100  $\mu$ l gel slice: 400  $\mu$ l Gel Dissolving Buffer: 150  $\mu$ l water).

- 4 Insert column into collection tube and load sample onto the column. Spin for 1 minute at 16,000 x g, then discard flow-through.**

 **00:01:00**

- 5 **Re-insert column into collection tube. Add 200 µl DNA Wash Buffer (with ethanol added) and spin for 1 minute at 16,000 x g.** Discarding flow-through is optional.

🕒 00:01:00

- 6 **Repeat Step 5** (Step 5: Re-insert column into collection tube. Add 200 µl DNA Wash Buffer and spin for 1 minute at 16,000 x g. Discarding flow-through is optional).

🕒 00:01:00

- 7 **Transfer column to a clean 1.5 ml microfuge tube.** Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for 1 minute.

Use care to ensure that the tip of the column does not come into contact with the flow-through. If in doubt, re-spin for 1 minute.

- 8 **Add ≥ 6 µl of DNA Elution Buffer to the center of the matrix. Wait for 1 minute, and spin for 1 minute at 16,000 x g to elute the DNA.**

🕒 00:02:00

Typical elution volumes are 6–20 µl. Nuclease-free water (pH 7–8.5) can also be used to elute the DNA. Yield may slightly increase if a larger volume of DNA Elution Buffer is used, but the DNA will be less concentrated. For larger size DNA (≥ 10 kb), heating the elution buffer to 50°C prior to use can improve yield.