



Sep 17, 2020

Batrachochytrium dendrobatidis electroporation

Scientific Reports

Andrew J. M. Swafford¹, Shane P. Hussey¹, Lillian K. Fritz-Laylin¹¹University of Massachusetts at Amherst**1** Works for me dx.doi.org/10.17504/protocols.io.88xhzxn

Andrew Swafford

ABSTRACT

Batrachochytrium dendrobatidis (BD) is a parasitic symbiote with many frogs. BD infection is devastating to most frog species, causing mass mortality in ~6.5% of amphibian species. To understand the cellular biology of BD and its interaction with amphibians, we have developed an method of electroporation, allowing for efficient deliver of extracellular constructs to infectious BD zoospores. The current method reliably loads 95% of treated cells with extracellular dextrans.

EXTERNAL LINK

<https://doi.org/10.1038/s41598-020-71618-2>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Swafford AJM, Hussey SP, Fritz-Laylin LK, High-efficiency electroporation of chytrid fungi. Scientific Reports doi: [10.1038/s41598-020-71618-2](https://doi.org/10.1038/s41598-020-71618-2)

DOI

dx.doi.org/10.17504/protocols.io.88xhzxn

PROTOCOL CITATION

Andrew J. M. Swafford, Shane P. Hussey, Lillian K. Fritz-Laylin 2020. Batrachochytrium dendrobatidis electroporation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.88xhzxn>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Swafford AJM, Hussey SP, Fritz-Laylin LK, High-efficiency electroporation of chytrid fungi. Scientific Reports doi: [10.1038/s41598-020-71618-2](https://doi.org/10.1038/s41598-020-71618-2)

EXTERNAL LINK

<https://doi.org/10.1038/s41598-020-71618-2>

KEYWORDS

chytrid, electroporation, bd, chytridiomycosis

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 12, 2019

LAST MODIFIED

Sep 17, 2020

GUIDELINES

Keep everything cold and work quickly when off ice.

Dispose of supernatant as quickly as possible as zoospores will swim up from out of the pellet after centrifugation.

MATERIALS

NAME	CATALOG #	VENDOR
Magnesium Chloride	AC223210010	Fisher Scientific
Sodium phosphate monobasic monohydrate	S9638	Sigma Aldrich
Gene Pulser Electroporation Cuvettes, 0.2 cm gap	1652086	Bio-rad Laboratories
D-Mannitol	MB0335.SIZE.100g	Bio Basic Inc.
Sodium phosphate, dibasic, anhydrous	S0404.SIZE.500g	Bio Basic Inc.
Sodium succinate, dibasic, hexahydrate	SB0889.SIZE.500g	Bio Basic Inc.
Paraformaldehyde	P6148	Sigma Aldrich
Potassium Chloride	P9541	Sigma Aldrich
sucrose	sucrose	Fisher Scientific
Dextran, Fluorescein, 3000 MW, Anionic	D3305	Thermo Fisher

SAFETY WARNINGS

This protocol requires the use of dangerous chemicals and biological hazards. Follow appropriate PPE, safety, and biological waste guidelines.

BEFORE STARTING

Synchronize cultures 2 hrs before counting zoospores

Media & Material Prep

10m

10m

- 1 Make 50 mL SM buffer (SM buffer should not be more than 7 days old).

1.1 SM Buffer:

- 5 mL [M]50 Milimolar (mM) KCL (Autoclaved)
- 750 µl pH7.2 [M]1 Molarity (M) Sodium Phosphate Buffer (Filter Sterilized)
- 7.5 mL [M]100 Milimolar (mM) MgCL2 (filter sterilized)
- .338 g Sodium Succinate dibasic acid hexahydrate
- .228 g D-Mannitol
- 36.75 mL Sterile water

10m

- 2 Make an appropriate volume of Paraformaldehyde fixation buffer.

2.1 Final Concentration of PFA fixation buffer:

- 4.8% PFA
- 9mM Sucrose
- 50mM Sodium Phosphate Buffer, pH 7.2

- 3 Make 2mg/mL dextrans in SM buffer (**Be sure it is Cat No. D3305 from Invitrogen by Thermo Fisher Scientific**) 5m
 - dextrans is sticky, measure in the tube you will be adding buffer to.
 - Spin down at high speed, then use supernatant.
- 4 Measure out the required amount of autoclaved growth media into an appropriate flask. 2m
- 5 Place fixation buffer, growth media, and bonner's salts on ice. 15m
- 6 Turn Centrifuge on to **4 °C**

Zoospore Prep

- 7 Transfer synced cells to 50mL conicals 5m
- 8 Spin down cells and discard supernatant **2500 rpm @ 4 °C** for **00:05:00** 6m
- 9 Resuspend in 10 mL of SM buffer, combine after resuspension if multiple tubes were needed to hold all the synced spores. 2m
- 10 **2500 rpm @ 4 °C** for **00:05:00** . Discard supernatant immediately. 6m
- 11 Resuspend in 5mL of SM 2m
- 12 **2500 rpm @ 4 °C** for **00:05:00** . Discard supernatant immediately. 5m
- 13 Resuspend cells in the appropriate volume of SM buffer needed for the # of cuvettes and tubes that will be used (**100 µl** SM for each cuvette) 2m
- 14 Count cell concentration. Aiming for 2×10^7 - 1×10^7 cells/mL 10m

14.1 If the spores are overly concentrated, dilute to the appropriate concentration. 2m

15 Prepare cuvettes & tubes.
E- = electroporation control (no voltage)
D- = Dextrans control (no dextrans)

15.1 Prepare E-D- tubes. 1.5mL eppendorfs: 1m

- 100 µl SM
- 100 µl Cell suspension

15.2 Prepare E-D+ tubes: 1.5ml eppendorfs 1m

- 100 µl Dextrans solution
- 100 µl Cell suspension

15.3 Prepare E+D- cuvettes: 0.2cm electroporation cuvettes 1m

- 100 µl SM
- 100 µl Cell suspension

15.4 Prepare E+D+ cuvettes: 0.2cm electroporation cuvettes 1m

- 100 µl Dextrans solution
- 100 µl Cell suspension

16 Set samples on ice for ~ 10 minutes 10m

Electroporation

17 Set electroporation protocol

17.1 2m

- Squarewave protocol
- 1000V
- 2 pulses
- 5 sec intervals
- 2mm cuvette

18 Gently pipette liquid in cuvette up and down to mix cells just before electroporation

19 Insert cuvette

- 20 Hit the PULSE button, record percent droop
- 21 Remove cuvette and place immediately on ice.
- 22 Repeat steps 18-21 for all cuvettes, be sure to pipette control eppendorf tubes even though they will not be exposed to a pulse.
- 23 ⌚ 00:10:00 allow samples to rest on ice 10m
- 24 Add 🧴 200 µl chilled growth media to each cuvette
- 25 ⌚ 00:10:00 allow samples to rest on ice

Wash cells

- 26 Transfer contents of each tube/cuvette to 15mL conicals.
- 27 🧴 5 mL chilled growth media to each sample
- 28 🌀 2500 rpm 🌡 4 °C ⌚ 00:05:00 discard supernatant
- 29 ⌚ go to step #27 Unless supernatant has been discarded three times.
- 30 add 🧴 5 mL Bonner's salts 🌡 4 °C Bonner's Salts
- 31 🌀 2500 rpm 🌡 4 °C ⌚ 00:05:00 discard supernatant
- 32 ⌚ go to step #30 until supernatant as been discarded 2 times

33 add  **600 µl Bonner's salts**  **4 °C**


34 Aliquot  **100 µl** cells for imaging

35 Aliquot  **100 µl** cells for counting


36 Fix remaining cells for flow cytometry


Cell Fixation

37 Add an equal volume of fixation buffer (see [go to step #2](#)) to each remaining volume of cells.

38 Immediately place tubes on ice.  **00:15:00** (centrifugation time in the following step is IN ADDITION to this time).
Leave alone for 15 minutes.

39
Centrifuge  **2200 rpm** @  **Room temperature**  **00:05:00**

40 Properly dispose of supernatant and resuspend in  **400 µl Bonner's salts**

41 Immediately quantify via flow cytometry or wrap in foil and store at  **4 °C** for up to 4 hours.