



Jul 10, 2021

## DNA extraction from frozen whole blood using the DNeasy Blood & Tissue kit (QIAGEN)

Luise A. Seeker<sup>1</sup>, Jennifer Fairlie<sup>1</sup>, Sarah L. Underwood<sup>1</sup>, Rachael V. Wilbourn<sup>1</sup>, Rebecca Holland<sup>1</sup>, Daniel H. Nussey<sup>1</sup>

<sup>1</sup>University of Edinburgh

1 Works for me



Share

[dx.doi.org/10.17504/protocols.io.i6cchaw](https://dx.doi.org/10.17504/protocols.io.i6cchaw)

Luise Seeker

### ABSTRACT

Even though DNA extraction from whole blood is a standard protocol of most commercial kits DNA quality and purity can be disappointing if those protocols are not optimised. Here, we present an optimised protocol for the extraction of high quality genomic DNA using the DNeasy Blood & Tissue kit (QIAGEN). Our downstream application of the DNA was the measurement of relative leukocyte telomere length by qPCR. The protocol might also be appropriate for other downstream applications. Changes reported here significantly improved DNA yield and 260/280 and 260/230 ratios measured on NanoDrop.

### ATTACHMENTS

[20170729\\_DNAExtraction  
Protocol\\_SpinColumns.docx](#)

### DOI

[dx.doi.org/10.17504/protocols.io.i6cchaw](https://dx.doi.org/10.17504/protocols.io.i6cchaw)

### PROTOCOL CITATION

Luise A. Seeker, Jennifer Fairlie, Sarah L. Underwood, Rachael V. Wilbourn, Rebecca Holland, Daniel H. Nussey 2021. DNA extraction from frozen whole blood using the DNeasy Blood & Tissue kit (QIAGEN). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.i6cchaw>

### KEYWORDS

high quality DNA

### LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

### CREATED

Jul 29, 2017

### LAST MODIFIED

Jul 10, 2021

### PROTOCOL INTEGER ID

7076

### MATERIALS TEXT

#### MATERIALS

[QIAGEN DNeasy Blood and Tissue Kit, 50](#)

[rxn Qiagen Catalog #69504](#)

[250ml RBC Lysis Buffer G-](#)

[Biosciences Catalog #786-650](#)

[1 Liter PBS \[10X\] \(Phosphate Buffered Saline\) \(80mM Na2HPO4, 1.5M NaCl, 20mM KH2PO4, 30mM KCl, pH 7.4\) G-](#)

[Biosciences Catalog #R027](#)

