



Jun 02, 2022

Assessment of PKC-dependent activation of LRRK1 in vitro

Athanasios Karapetsas¹, Asad Malik¹, Dario R Alessi¹¹Medical Research Council Protein Phosphorylation and Ubiquitylation Unit, School of Life Sciences, University of Dundee, Dow Street, Dundee DD1 5EH, UK

2

dx.doi.org/10.17504/protocols.io.5jyl89d5rv2w/v1

asap

Dario Alessi

We describe a non-radioactive assay that we deploy for analysing the kinase activity of recombinant LRRK1 following *in vitro* activation by Protein kinase C (PKC) isoforms. This assay can also be used to analyse the effect of PKC on LRRK1 immunoprecipitated from cells.

[393-856.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.5jyl89d5rv2w/v1

Athanasios Karapetsas, Asad Malik, Dario R Alessi 2022. Assessment of PKC-dependent activation of LRRK1 in vitro. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl89d5rv2w/v1>



PKC activation, lipid vesicles, Phosphorylation of LRRK1, immunoblotting analysis, ASAPCRN

protocol ,

Mar 28, 2022

Jun 02, 2022

Mar 28, 2022 renuka.s

Mar 31, 2022 Dario Alessi

Note: Once the *in vitro* kinase assay has been performed, we recommend analysing the reaction products by quantitative immunoblotting (as described in [XXXXXX](#)).

Note: This protocol can be adapted to analyse activation of LRRK1 that has been immunoprecipitated from cells (as described in [XXXXXX](#)).

Reagents:

1. Recombinant PKC isoform protein (available from MRC Reagents and Services: <https://mrcpureagents.dundee.ac.uk/>)
2. Recombinant Rab7A protein (available from MRC Reagents and Services: <https://mrcpureagents.dundee.ac.uk/>)
3. Recombinant LRRK1 wild type [27-2015] protein

Note: Recombinant LRRK1 protein is expressed and purified by following the protocol described in: [XXXXXX](#)

4. Kinase assay buffer:

A	B
HEPES pH 7.5	25 mM
2-mercaptoethanol	0.1% (v/v)
KCl	50 mM
CaCl ₂	1 mM
MgCl ₂	10 mM
ATP	1 mM

1. L- α -Phosphatidylserine (Avanti Polar Lipids, resuspended in methanol and chloroform at a 1:1 ratio for long-term storage)
2. L- α -Diacylglycerol (Avanti Polar Lipids, resuspended in methanol and chloroform at a 1:1 ratio for long-term storage)
3. **4X Loading buffer:**

[NUPAGE LDS sample buffer \(4x\)](#) Thermo Fisher

Scientific Catalog #NP0007

, or

4X SDS loading buffer

A	B
Tris-HCl, pH6.8	250mM
SDS	8% (w/v)
Glycerol	40% (v/v)
Bromophenol blue	0.02% (w/v)

 [Anti-Rab7 antibody Mouse monoclonal](#) **Sigma**

Aldrich Catalog #R8779

 [Recombinant Anti-PKC alpha](#)

antibody Abcam Catalog #ab11723

Equipment:

- Refrigerated bench-top centrifuge (Eppendorf microcentrifuge 5417R, or equivalent).

Refrigerated Centrifuge
Centrifuge

Eppendorf EP-5417R 

- Savant SpeedVac system (Thermo #SPD140DDA, or equivalent)
- Thermo mixer (Eppendorf ThermoMixer, or equivalent)
- Disposable Glass Culture Tubes (Fisherbrand Round Bottom Disposable Borosilicate Glass Tubes, or equivalent)

Preparation of lipid vesicles for PKC activation

1 

Clean a disposable glass culture tube by washing. Allow to air-dry.

1.1 Clean a disposable glass culture tube by washing with 100% methanol. (1/3)

1.2 Clean a disposable glass culture tube by washing with 100% methanol. (2/3)

1.3 Clean a disposable glass culture tube by washing with 100% methanol. (3/3)

2 

Pipette  **0.5 µL** of Diacylglycerol (stock concentration is  **10 mg/mL**) and  **5 µL** of

Phosphatidylserine (stock concentration is **10 mg/mL**) into the cleaned and dried glass tube.

Note: These quantities will provide sufficient lipid vesicles for 25 reactions at a volume of **20 µL** per reaction.

- 3 Vacuum dry lipids using a SpeedVac system for **00:10:00**. This should leave a visible, ^{10m} translucent lipid pellet.

Note: Ensure that lipids are completely dried as any residual chloroform or methanol will inhibit the kinase reaction.

- 4 Resuspend lipids in **50 µL** of **25 millimolar (mM)** HEPES **pH 7.4**, **50 millimolar (mM)** KCl. Vortex gently until pellet is no longer visible.

Kinase Reaction Step 1: Phosphorylation of LRRK1 by PKC

5

Note: If using immunoprecipitated LRRK1 from cells, perform immunoprecipitation and washes (as described in **XXXXXX**) before proceeding with Step 5.

Prepare a primary “2X master mix” containing

A	B
HEPES pH 7.5	50 mM
KCl	100 mM
2-Mercaptoethanol	0.2% (v/v)
MgCl ₂	20 mM
ATP	2 mM
CaCl ₂	2 mM
Phosphatidylserine	200 µg/ml
Diacylglycerol	20 µg/ml

6



For each reaction, add **10 µL** of the primary “2X master mix” to a clean Eppendorf tube.

7



5m

Add **5 µL** of **200 nanomolar (nM)** LRRK1 wild type protein (final concentration is **50 nanomolar (nM)**) to each reaction and allow equilibration **On ice** for **00:05:00**.

Note: If using LRRK1 immunoprecipitated from cells, add **10 µL** of the primary “2X master mix” and **5 µL** of **25 millimolar (mM)** HEPES **pH 7.5**, **50 millimolar (mM)** KCl, 0.1% (v/v) 2-Mercaptoethanol to each tube containing beads-bound immunoprecipitated LRRK1.

8



Start the kinase reaction by adding **5 µL** of **400 nanomolar (nM)** PKC Alpha protein (final concentration is **100 nanomolar (nM)**).

Note: The final reaction volume should be **20 µL**.

9

After **00:30:00**, transfer the Eppendorf tubes from Step 8 **On ice**.

30m

Kinase Reaction Step 2: Phosphorylation of Rab7A by PKC-activated LRRK1

10

Prepare a secondary “master mix” (=Master Mix B) containing

A	B
HEPES pH 7.5	25 mM
KCl	50 mM
MgCl ₂	10 mM
ATP	1 mM
Rab7A	1 μ M

11



Start the second step of the kinase reaction by adding **10 μ L** Master Mix B to the Eppendorf tubes from Step 5.

12



45m

Transferring the Eppendorf tubes to the thermo mixer set at **30 $^{\circ}$ C**, 1,000 rpm. Incubate for **00:45:00**.

13



Stop the kinase reaction by adding **10 μ L** of 4 \times LDS (supplemented with 5% (v/v) 2-Mercaptoethanol) loading buffer to the reaction mix to a final concentration of 1 \times .

14



10m

If using LRRK1 immunoprecipitated from cells, stop the kinase reaction by adding **30 μ L** of 4 \times LDS loading buffer to the reaction mix to a final concentration of 2 \times , incubate the mixture at **70 $^{\circ}$ C** on a heat block for **00:10:00** to elute LRRK1 from the resin, and collect the eluent by centrifugation through a 0.22- μ m-pore-size Spinex column.

15



5m

Incubate the samples for **00:05:00** at **70 $^{\circ}$ C** on a heat block before proceeding to quantitative immunoblotting analysis.

Note: If using LRRK1 immunoprecipitated from cells, supplement the samples from Step 14 with 2-Mercaptoethanol to 1% (v/v) before proceeding to Step 15.

Analysis of kinase reaction products by quantitative immunoblotting analysis:

1h 15m

16



The reaction products can be analysed by quantitative immunoblotting analysis (as described in XXXX). **Table 1** lists the primary antibodies that we recommend using, which include antibodies to detect Rab7A phosphorylation at Serine-72.

A	B	C	D	E
Antibody Target	Company	Cat. number	Host species	Dilution
pS72 Rab7A	Abcam Inc.	MJF-38, Clone 1	Rabbit	1 mg/ml
Rab7A (Total)	Sigma	R8779	Mouse	1.430556
LRRK1 (total) (C-terminus)	MRC-PPU Reagents and Services, University of Dundee	S405C	Sheep	1 mg/ml
PKC Alpha	Abcam Inc.	ab11723	Mouse	1.430556

17

1h 15m






Figure 1: PKC alpha dose-dependent activation of recombinant LRRK1 in vitro. Recombinant LRRK1 wild type [27-2015] was incubated with increasing concentrations of PKC Alpha (1 to 300 nM) at 30 °C for 00:30:00 with excess Mg-ATP. Reactions were subsequently incubated with 1 micromolar (μM) recombinant Rab7A and subjected to a

 **00:45:00** kinase reaction at  **30 °C** in the presence of excess Mg-ATP. Kinase reactions were subjected to immunoblot analysis with the indicated antibodies and the membranes were developed using the Odyssey CLx scan Western Blot imaging system.