



JUL 13, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.261genwoog47/v1

Protocol Citation: Mason Mackey, Mark Ellisman, Marianne Martone 2023. Karnovsky's Modified Fixative. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.261genwoog47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Apr 20, 2022

Last Modified: Jul 13, 2023

PROTOCOL integer ID:
 61106

Keywords: karnovsky, fixative, mice, rat, paraformaldehyde, Glutaraldehyde, cacodylate buffer

Karnovsky's Modified Fixative

Mason

Mackey¹, Mark Ellisman^{1,2}, Marianne Martone¹

¹National Center for Microscopy and Imaging Research, University of California, San Diego, La Jolla CA;

²Department of Neurosciences, University of California, San Diego, La Jolla CA

NCMIR@UCSD



NCMIR

University of California, San Diego

ABSTRACT

Fixative for mice and rat tissues to be used during perfusion and immersion fixation.

GUIDELINES

Preparing Fix solution (2% paraformaldehyde and 2.5% glutaraldehyde)

For 100 mL of 2% paraformaldehyde and 2.5% glutaraldehyde in 0.15M cacodylate buffer solution at 7.4 pH.

MATERIALS

- 100mL Erlenmeyer flask
- 100mL graduated cylinder
- #1 Whatman filter
- Paraformaldehyde PRILLS (Catalog no. 19202)
- Glutaraldehyde EMS (Catalog no. 16220)
- Sodium Cacodylate buffer TED PELLA, INC. (Catalog no. 18851)
- Stir bar

SAFETY WARNINGS



Always wear appropriate PPE.

BEFORE START INSTRUCTIONS

Make sure heating plate for dissolving paraformaldehyde is between 60-70° C.

- 1 Add 2.0 grams paraformaldehyde to a 100mL Erlenmeyer flask.
- 2 Add 35mL of DDH₂O to 100mL to a Erlenmeyer flask.
- 3 Heat paraformaldehyde in DDH₂O to 60-70°C. Use stir bar to mix.
- 4 Add 2-3 drops 1N NaOH.
- 5 Filter solution with #1 Whatman filter into a 100mL graduated cylinder.
- 6 Add 10.0 ml of 25% glutaraldehyde.
- 7 Dilute with DDH₂O to a total volume 49.0 ml.
- 8 Add 1.0 ml of 0.204 M CaCl₂·2H₂O.
- 9 Add 50 ml 0.3M cacodylate buffer stock solution at pH 7.4.

10 Volume should be 100 ml. Transfer fixative into a 100mL Erlenmeyer flask.

11 Keep in bath water 35°C.