



Oct 05, 2022

HILIC-MS analysis of central carbon metabolites in gram negative bacteria

Bashar Amer¹, Ramu Kakumanu¹, Edward EK Baidoo¹¹Lawrence Berkeley National Laboratory

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.4r3l2opzxv1y/v1

Edward EK Baidoo

Lawrence Berkeley National Laboratory, Joint BioEnergy Insti...

ABSTRACT

Central carbon metabolites are essential for energy metabolism, secondary metabolite and macromolecule biosynthesis in prokaryotes. Phosphorylated intermediates such as sugar phosphates are, however, difficult to resolve with conventional reversed phase (RP) chromatography. In addition, coenzyme A (CoA) compounds and nucleotide cofactors, which are also phosphorylated are difficult to resolve by RP chromatography. To this end, we have developed a hydrophilic interaction liquid chromatography quadrupole time-of-flight mass spectrometry (HILIC-QTOF-MS) method that to quantify these metabolic intermediates and cofactors in *Escherichia coli* DH1 and *Pseudomonas putida* K2440.

ATTACHMENTS

[fverbdeip.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.4r3l2opzxv1y/v1

PROTOCOL CITATION

Bashar Amer, Ramu Kakumanu, Edward EK Baidoo 2022. HILIC-MS analysis of central carbon metabolites in gram negative bacteria. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.4r3l2opzxv1y/v1>



FUNDERS ACKNOWLEDGEMENT

US Department of Energy
Grant ID: DE-AC02-05CH11231

KEYWORDS

Metabolomics, sugar phosphates, nucleotides, coAs, organic acids, HILIC, LC-MS

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 10, 2022

LAST MODIFIED

Oct 05, 2022

PROTOCOL INTEGER ID

64367

MATERIALS TEXT

Solvents used

LC-MS grade methanol (part number LC230-4, Honeywell Burdick & Jackson, Charlotte, NC, USA), LC-MS grade water (part number BJLC365-4, Honeywell Burdick & Jackson, Charlotte, NC, USA), and LC-MS grade acetonitrile (part number LC015-4, Honeywell Burdick & Jackson, Charlotte, NC, USA).

Reagents used

[M]5 millimolar (mM) InfinityLab Deactivator Additive - **25 mL** (part number 5191-3940, Agilent Technologies, Santa Clara, CA, USA), ammonium hydroxide (28–30% NH₃ basis, part number AX1303, EMD Chemicals, Billerica, MA, USA), and ammonium acetate (part number 17836, ACS reagent, from Sigma-Aldrich, St. Louis, MO, USA).

Analytical standards

Chemical standards were purchased from Sigma-Aldrich.

MS calibrants

Reference mass correction (100 mM trifluoroacetic acid ammonium salt (part number I8720243) and 2.5 mM HP-0921 (part number 18720241)) and ESI-L-Low concentration tuning mix (part number G1969-85000) were purchased from Agilent Technologies (Santa Clara, CA, USA).

UHPLC system, LC column, and guard column

The Agilent Technologies 1290 Infinity II UHPLC system was used throughout. An Agilent Technologies InfinityLab Poroshell 120 HILIC-Z guard column (**2.1 mm** inner diameter, **2.7 μm** stationary phase particle size, part number 821725-947) was connected to an Agilent Technologies InfinityLab Poroshell 120 HILIC-Z PEEK-lined column (**2.1 mm** inner diameter, **100 mm** length, **2.7 μm** stationary phase particle size, part number 675775-924).

QTOF-MS system

An Agilent Technologies 6545 (part number G6545B) Quadrupole Time of Flight (QTOF) LC/MS system was used throughout.

SAFETY WARNINGS

Wear the appropriate PPE protection (i.e., gloves, safety goggles, and lab coat) and prepare solvents and LC-MS mobile phases in a chemical fume hood. Store organic solvents in a flammable storage cabinet and peroxide-forming chemicals in the appropriate safety storage cabinets.

Sample preparation

- 1 Calibration curves: Pure analytical standards were dissolved in 50:50 methanol:water (v/v) to make a stock solution of **[M]100 micromolar (μM)** concentration. The stock solution was diluted 4-fold to make a standard solution of **[M]25 micromolar (μM)** (with 50:50

methanol:water, v/v). A seven-point calibration curve was produced via a series of 2-fold serial dilutions (with 50:50 methanol:water, v/v), which were conducted from the [M]25 micromolar (μM) standard solution to a standard solution ending in a concentration of [M]0.390625 micromolar (μM) ; 2-fold serial dilutions for CoA standards, however, were conducted from a [M]20 micromolar (μM) standard solution to standard solution ending in a concentration of [M]0.3125 micromolar (μM) . All solutions were prepared on ice and stored at 4 -20 °C until LC-MS data acquisition.

Extraction of intracellular metabolites: Wild type *Escherichia coli* DH1 (*E. coli*) and *Pseudomonas putida* K2440 (*P. putida*) were grown on 0.5% glucose in LB medium at 4 37 °C . Growth was monitored during mid exponential phase to early stationary phase. *E. coli* and *P. putida* reached average optical densities at 600 nm wavelength (OD 600 nm) of 1.67 and 2.16, respectively. 1.5 mL of the cell culture was harvested (in triplicate) and transferred to a 2 mL centrifuge tube. The cell culture was centrifuged at 14000 x g for 2 minutes at room temperature. After removing the supernatant, metabolic activity from the cell pellet was quenched by adding 250 μL of ice-cold methanol and mixing thoroughly by vortexing for 2 minutes at 4 4 °C . 250 μL of water was then added to the methanol lysate, mixed by vortexing, and centrifuged at 14000 x g for 5 minutes at 4 4 °C . 450 μL of the lysate was filtered via a 3k Da molecular weight cut-off centrifuge filter (Amicon-Ultra, part number UFC500396, MilliporeSigma, St. Louis, MO, USA) at 13000 x g at 4 4 °C , for 60 minutes prior to LC-MS analysis.

Baidoo EEK, Wang G, Joshua CJ, Benites VT, Keasling JD (2019). Liquid Chromatography and Mass Spectrometry Analysis of Isoprenoid Intermediates in *Escherichia coli*. Methods in molecular biology (Clifton, N.J.). https://doi.org/10.1007/978-1-4939-8757-3_11

Ultra high performance liquid chromatography (UHPLC) conditions

- 2 Chromatographic separation was performed via an Agilent Technologies 1290 Infinity II UHPLC system (Table 1). The sample tray and column compartment were set to 4 5 °C and 4 30 °C , respectively. A sample injection volume of 1 μL was used throughout. The UHPLC mobile phases were composed of solvents (A) [M]10 millimolar (mM) ammonium acetate, [M]0.2 % volume ammonium hydroxide, and [M]5 micromolar (μM) medronic acid in water

and (B) **10 millimolar (mM)** ammonium acetate, **0.2 % volume** ammonium hydroxide, and **5 micromolar (µM)** medronic acid in **90 % volume** acetonitrile, and water as the remaining solvent. The mobile phases were degassed by sonication for 5 minutes. The 1290 Infinity II UHPLC system pump was purged with the aforementioned mobile phases for up to 5 minutes (at **50 % volume** solvent A and **50 % volume** solvent B) prior to LC column installation. The analytes were separated on an Agilent Technologies InfinityLab Poroshell 120 HILIC-Z PEEK lined (**2.1 mm** internal diameter, **100 mm** length, and **2.7 µm** stationary phase particle size) column, which was connected to an Agilent Technologies InfinityLab Poroshell 120 HILIC-Z guard column (**2.1 mm** internal diameter and **2.7 µm** stationary phase particle size). The LC column was then equilibrated with the starting mobile phase composition at a flow rate of 0.1 mL/min until the system backpressure was stable. The UHPLC gradient is described in table 2.

A	B	C
Component	Description	Part number
1260 Infinity II Isocratic Pump	Isocratic pump for internal reference lock mass delivery	G7110B
1290 Infinity II MCT	Column compartment	G7116B
1290 Infinity II Multisampler	Autosampler for automatic sample injection	G7167B
1290 Infinity II High Speed Pump	UHPLC pump	G7120A
6545 Q-TOF LC/MS	Quarupole-time-of-flight mass spectrometer	G6545B

Table 1. The Agilent Technologies UHPLC-QTOF-MS system.

A	B	C	D
Time (min)	Percentage of mobile phase B (%B)	Flow rate (mL/min)	Maximum system pressure (Bar)
0	90	0.250	600
2	90	0.250	600
11.5	61.5	0.250	600
11.7	60	0.330	600
12.7	60	0.400	600
12.9	90	0.550	600
14.3	90	0.550	600

Table 2. UHPLC gradient. The maximum allowable system backpressure for the column was 600 bar.

QTOF-MS method parameters

- 3 The Agilent Technologies Infinity II UHPLC system was coupled to an Agilent Technologies 6545 QTOF-MS system. The LC column effluent was delivered to the Agilent Jet Stream (AJS) ion source. An Agilent Technologies Isocratic Pump was used to deliver the reference mass correction solution. AJS facilitated the production of gas-phase $[M - H]^-$ or $[M - 2H]^{2-}$ ions in the negative ion mode. High mass accuracy was achieved via reference mass correction. Reference mass correction was performed with the trifluoroacetic acid (anion) $[M - H]^-$ at 112.98559 m/z (from the LC mobile phase) and the trifluoroacetic acid anion adduct of HP-0921 (i.e. hexakis(1H, 1H,3H-tetrafluoropropoxy)phosphazine (C₁₈H₁₈F₂₄N₃O₆P₃, 921.23 Da, CAS NO. 58943-98-9) at 1033.98811 m/z . The reference ions were diluted to 2 μ M in 80:20 acetonitrile:water (v/v) and delivered to an electrospray ionization sprayer at a flow rate of 5 μ L/min via the 1260 Infinity II Isocratic Pump (Table 1). Prior to analysis, the Agilent 6545 QTOF-MS system was tuned via the ESI-L-Low calibration solution for an acquisition range of up to 1700 m/z . MS data acquisition parameters are described in table 3.

A	B
QTOF-MS parameters	Values
Acquisition range (m/z)	70-1100 m/z
Acquisition rate (spectra/s)	0.86
Nebulizer pressure (*Psi)	20
Drying gas temperature (°C)	300
Drying gas flow rate (L/min)	10
Sheath gas temperature (°C)	350
Sheath gas flow (L/min)	12
Capillary voltage (V)	3500
Fragmentor (V)	100
Skimmer (V)	50
OCT 1 RF Vpp (V)	300
Nozzle voltage [Expt] (V)	2000

Table 3. QTOF-MS parameters. *Psi is lb/in².

LC-MS data acquisition and analysis software

- 4 LC-MS data acquisition was performed via the Agilent MassHunter Workstation software (version 8). Data processing and analysis were performed via Agilent MassHunter Qualitative Analysis (version 6) and Profinder (version 8) software.

Method validation

5




Figure 1. Extracted ion chromatograms.

Compound	Mass-to-charge ratio (m/z)	Linearity (R^2 coefficient)	Limit of detection ($s/n=3$)	% RSD of standard retention time ($n=6$)	Standard retention time (min) ($n=6$)	Metabolite retention time (min) ($n=6$, * $n=3$)
Glucose 6-phosphate	259.022442	0.9964	0.22	0.18	8.68	8.76
Fructose 6-phosphate	259.022442	0.9994	0.23	0.08	8.27	8.30
Fructose 1,6-bisphosphate	338.988772	0.9996	0.24	0.12	9.64	9.66
Glyceraldehyde 3-phosphate	168.990748	0.9973	0.50	0.14	7.86	7.88
Phosphoglycerate	184.985663	0.9998	0.02	0.13	8.48	8.50
Phosphoenolpyruvate	166.975098	0.9985	0.07	0.03	8.75	8.75
Pyruvate	87.008768	0.9998	0.63	1.09	1.76	1.74
Succinate	117.019332	0.9955	0.16	0.19	7.54	7.51
6-Phosphogluconate	275.017357	0.9980	0.77	0.46	8.93	8.89
Sedoheptulose 7-phosphate	289.033007	0.9974	0.36	0.14	8.41	8.46
Erythrose 4-phosphate	199.001313	0.9933	0.70	0.14	8.19	8.34
Ribulose 5-phosphate	229.011877	0.9930	0.14	0.12	7.86	7.92
Ribose 5-phosphate	229.011877	0.9949	0.45	0.21	8.18	8.19
ATP	505.988469	0.9924	0.28	0.11	8.48	8.55
NAD	662.101845	0.9991	1.67	0.27	7.73	*7.76
NADH	664.117495	0.9999	0.35	0.17	6.82	*6.82
Acetyl-CoA	403.555610	0.9985	0.23	0.29	7.65	7.65
Free CoA	382.550328	0.9900	0.25	0.20	7.91	7.91

Figure 2. Method validation table. The mass-to-charge ratios shown are theoretical values. The CoA compounds were detected via $[M - 2H]^{2-}$ ions. All other compounds were detected via $[M - H]^-$ ions. Average standard and metabolite retention times are used in the last two columns of fig. 2. Phosphoglycerate data represents 2- and 3-phosphoglycerate, which could not be resolved by the method. ATP, NAD, and NADH are abbreviations for adenosine triphosphate, nicotinamide adenine dinucleotide, and nicotinamide adenine dinucleotide, reduced, respectively.

Biological data

6




Figure 3. Specific sugar phosphate concentrations in *E. coli* and *P. putida*.




Figure 4. Specific CoA and nucleotide cofactor concentrations in *E. coli* and *P. putida*.




Figure 5. Specific organic acid concentrations in *E. coli* and *P. putida*.

Conclusions

- 7 The HILIC-QTOF-MS method was able to separate the analytes tested, except for the isomers 2- and 3-phosphoglycerate (fig. 1). The method showed good linearity for all the analytes tested, with R^2 values of ≥ 0.99 . The limit of detection ranged from **20 nanomolar (nM)** to **1.67 micromolar (μ M)** for the analytes detected. Retention time repeatability measurements were ≤ 0.49 %RSD for the analytes tested, except for pyruvate, which had a %RSD value of 1.09. There was very little difference between standard and metabolite retention times, suggesting a minimal effect of the sample matrix on retention time stability (fig. 2). The method successfully quantified intracellular sugar phosphate (fig. 3), acetyl-CoA and free CoA (fig. 4), and organic acid metabolites (fig. 5) in *E. coli* and *P. putida* metabolite extracts. Additionally, the method was able to quantify the intracellular nucleotide cofactors ATP, NAD, and NADH (fig. 4) in these microbial metabolite extracts. The metabolite concentration data suggests that the method may be applicable for other phosphorylated metabolites and organic acids.