



Apr 05, 2021

# Neurolucida 360: Detecting muscle fiber orientation

Maci Heal<sup>1</sup><sup>1</sup>MBF Bioscience

1

Works for me

This protocol is published without a DOI.

SPARC

Tech. support email: [info@neuinfo.org](mailto:info@neuinfo.org)

sbaldwin

## ABSTRACT

Detecting muscle fiber orientation in 3D microscopy images using Neurolucida 360.

## PROTOCOL CITATION

Maci Heal 2021. Neurolucida 360: Detecting muscle fiber orientation . [protocols.io](https://protocols.io/view/neurolucida-360-detecting-muscle-fiber-orientation-br4qm8vw)  
<https://protocols.io/view/neurolucida-360-detecting-muscle-fiber-orientation-br4qm8vw>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Feb 04, 2021

## LAST MODIFIED

Apr 05, 2021

## PROTOCOL INTEGER ID

46960

### Setup

- 1 Launch SPARC-enabled Neurolucida 360.

#### Neurolucida 360

by MicroBrightField Bioscience

- 2 Open a microscopy image via the Open icon, File>Open, or dragging and dropping into the program window.

NOTE: MBF Bioscience software supports a variety of image file formats from common microscopy vendors such as JPX/JP2, TIFF, LIF, ND2, IMS, OIF/OIB, almost all CZI, and more.

- 3 The SPARC Vocabulary Services window will appear. The dialog is displayed so that you can specify subject-specific metadata and gain access to the [SciCrunch](#) database to retrieve anatomical terminology lists compliant with FAIR data principles to use during annotation.

- 3.1 Fill out the **Subject Information** and select your **Criteria for Anatomic Terms**. Then select **Begin**.

Welcome!  
To begin annotation, specify the following:

Subject information:

Species:

Subject ID:

Sex:

Age:

Select criteria for anatomic terms:

Organ:

Species:

Parcellation:

☒ Suppress this dialog this session

#### Using the Detect Orientation tool

- 4 When your image is loaded, close the 3D window. From the **TRACE** ribbon, select the **Detect Orientation** icon.




- 5 Specify the size of objects to detect in the **Detect Orientation** dialog box.

Detect Orientation

Detect orientation of objects of this size:   $\mu\text{m}$


- 5.1 NOTE: If you are unsure of the size of your object of interest, find the measure line tool from the **TRACE** ribbon.




- 6 Once your object size is specified, select **Detect Orientation**. The detection algorithm will then detect the orientation of the local region of image signal at that designated size.

## Detect Orientation

A progress indicator will display while Neurolucida 360 software determines muscle fiber orientation.



- 7 View results in the 3D window. The muscle fiber orientations are modeled with red arrows overlaid on the image data.



- 8 Save the data file in XML format.

The orientation data for all detected vectors will be written to the XML file for use in your individual analysis pipelines and/or MAP-Core segmentation registrations to organ scaffolds.

- 8.1 The algorithm detects the orientation of image regions based on a user-defined size (orientation *scale*). The orientation *score* for each detected vector, represented as an arrow is calculated; values range from 0 to 1, with 1 representing perfect alignment in a principal direction. Only the fiber direction vectors are displayed in Neurolucida 360 software, however, the orientation of the *sheet* vector (lateral direction normal to the fiber in the plane of the sheet) and the vector normal to these two are recorded in the XML data file.