



Mar 04, 2022

TSS, TXS, SDS Serial Extraction

Haley Geertsma¹¹University of Ottawa

1

dx.doi.org/10.17504/protocols.io.b5vvq666 **Haley Geertsma**
University of Ottawa

This protocol is used to serially extract protein from mouse brain tissue in increasingly stringent buffers.

DOI

dx.doi.org/10.17504/protocols.io.b5vvq666

Haley Geertsma 2022. TSS, TXS, SDS Serial Extraction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5vvq666>



protocol ,

Mar 03, 2022

Mar 04, 2022

59029

- 1 Weigh brain sample and add 3X volume of TSS buffer + 1X protease inhibitor. 2m
TSS Buffer: 140mM NaCl + 5mM Tris-HCl in H₂O
- 2 Homogenize in a dounce homogenizer then transfer to a clean tube. 2m
- 3 Centrifuge at 21130g for 30 minutes at 4°C then save the supernatant as 'TSS fraction'. 30m

Resuspend the pellet in the same volume of TSS buffer then centrifuge at 21130g for 30 35m

- 4 minutes at 4°C and discard the supernatant.
- 5 Resuspend the pellet in the same volume of TXS buffer then incubate on ice for 10 minutes.^{15m}
TXS Buffer: TSS buffer + 0.5% Triton X-100
- 6 Centrifuge at 21130g for 30 minutes at 4°C then save the supernatant as 'TXS fraction'.^{30m}
- 7 Resuspend the pellet in the same volume of TSS buffer then centrifuge at 21130g for 30 minutes at 4°C and discard the supernatant.^{35m}
- 8 Resuspend the pellet in the same volume of SDS buffer then incubate at room temperature for 10 minutes.^{15m}
SDS Buffer: TSS buffer + 1% sodium dodecyl sulfate (SDS)
- 9 Centrifuge at 21130g for 30 minutes at room temperature then save the supernatant as 'SDS fraction'.^{30m}