



Oct 02, 2022

🌐 Salting-out Protocol for Extracting HMW Genomic DNA from frozen octocorals

Santiago Herrera¹¹Lehigh University

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.bypypvpw

High molecular weight DNA extraction from all kingdoms

Long Read Club

 Santiago Herrera

DISCLAIMER

Tom Jenkins' protocol modified by Santiago Herrera, 20191016.

ABSTRACT

Salting-out protocol for extracting high-molecular weight (HMW) genomic DNA from frozen octocorals for genomic methods. Suitable for genome sequencing, RAD-seq and other methods that require HMW DNA. Starting material must be alive when flash-frozen in liquid nitrogen and always stored below -80C.

DOI

dx.doi.org/10.17504/protocols.io.bypypvpw

PROTOCOL CITATION

Santiago Herrera 2022. Salting-out Protocol for Extracting HMW Genomic DNA from frozen octocorals. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bypypvpw>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 30, 2021

LAST MODIFIED

Oct 02, 2022

PROTOCOL INTEGER ID

53720

GUIDELINES

Follow steps 1-12

MATERIALS TEXT

Reagents

Proteinase K (20 mg/ml; Invitrogen; optimal digestion at 65°C)

RNase A (100 mg/ml) optional

Cell lysis buffer (100mM Tris-Cl; 50 mM EDTA; 1 % SDS)

0.5M EDTA

β-Mercaptoethanol (very toxic, handle with extreme care in fume hood only)

7.5Mammonium acetate (cooled at 4°C)

100% isopropanol (cooled at 4°C)

70% ethanol (cooled at 4°C)

TE pH 8.0 or AE (Qiagen) buffer

Liquid N2

Dry Ice

Distilled water

Crushed ice

Supplies

Kim-wipes(small)

Aluminium foil

Whirl-pak bags (4oz)

1.5 ml microcentrifuge tubes (low-binding)

2.0 ml screw capped microcentrifuge tubes

Lighter

2x Microcentrifuge racks

Equipment

Liquid Nitrogen Cooled Mini Mortar & Pestle Set

3L Stainless Steel Dewar Flask

2x Small Styrofoam coolers

IncuShaker

Hot plate

Ultracentrifuge

Vortex

Fume hood

Small forceps

Squirt bottle for distilled water

Alcohol burner

Small beaker

SAFETY WARNINGS

Discard all liquids and solids downstream of addition of **β-Mercaptoethanol** (step 6) in appropriate hazardous waste containers. Wear full PPE when handling this chemical and be careful. Only manipulate in fume hood. Avoid any contact with skin. Re-seal stock bottle with parafilm.

DISCLAIMER:

Preparation

- 1 Make sure that there are enough Materials and Supplies available and at the right temperature conditions.
- 2 Make sure that the Ultracentrifuge, Heat Block and IncuShaker are set up in the fume hood.
- 3 Place tube rack in IncuShaker.
- 4 Turn on the Heat Block and the IncuShaker and set to 65°C.
- 5 Pre-label 3 sets of 1.5ml microcentrifuge tubes for each sample (label one of the three tubes [the last one to be used] with a green tough-spot and full SH ID in addition to labelling the side of the tube with the date of the extraction.
- 6 Prepare the following digestion mix for each tube (mix in a 2.0ml tube, scale accordingly for up to 4 samples) and pre-heat in hot plate: 350µl cell lysis buffer, 42µl 0.5 M EDTA.
- 7 Fill 3L Stainless Steel Dewar Flask with liquid N₂
- 8 Fill a small Styrofoam cooler with dry ice
- 9 Fill a small Styrofoam cooler with crushed ice
- 10 10. Take cryovials containing frozen tissues out of -80C freezer and place in styrofoam cooler with dry ice.

- 11 11. Take Proteinase K, 7.5M ammonium acetate, 100% isopropanol, and 70% ethanol stocks out of fridge and place in cooler with crushed ice.
- 12 12. Pour liquid N₂ from the Stainless Steel Dewar Flask into the Liquid Nitrogen Cooled Mini Mortar base up to the fill line and then place the metal microcentrifuge tube mortar holder in the base to chill. Digestion

Digestion

- 13 13. Sterilize the small forceps and the metal Pestle by exposing to the flame of the alcohol burner for ~10 seconds. Cool each tool by squirting distilled water on the tools (place collection beaker below to catch the water). Dry up the tools by shaking them vigorously.
- 14 14. Place forceps on top of a clean piece of aluminium foil.
- 15 15. Place the tip of the metal Pestle in a new Whirl-Pak bag
- 16 16. Place a 1.5 mL microcentrifuge tube (corresponding to the right sample) in the chilled MiniMortar.
- 17 17. Break a small piece of coral tissue (~3mm³) with the forceps (it should be very brittle and easy to break) and place in the microcentrifuge tube that is sitting in the MiniMortar.
- 18 18. Wipe forceps with clean Kim-wipe and place on the Aluminium foil.
- 19 19. Chill the Pestle by immersing the tip in liquid N₂ (~10s in the Dewar Flask) while keeping it in the Whirl-pak bag.
- 20 20. Grind the sample to a powder with the Pestle until you don't feel any more crunching.

- 21 21. Wipe Pestle with clean Kim-wipe and place on Aluminium foil.
- 22 22. Loosen pulverized tissue by quickly flicking the bottom end of the tube with your finger.
- 23 23. Add pre-heated 392µl of digestion mix and quickly mix by inverting and flicking. Place tube in heat block.
- 24 24. Add 10µl of Proteinase K to the mix. Place in tube rack in IncuShaker (minimize door opening).
- 25 25. Repeat steps 13-24 for each sample.
- 26 26. Incubate at 65oC for 2.5-3 hours. Mix digestion rack by inversion every 20-30min
- 27 27. Optional Add 2 µl RNase A. Mix by vortexing. Incubate at 37oC for 30 minutes.

Protein and Cellular Debris Precipitation

- 28 28. Add 140 µl 7.5 M ammonium acetate to each tube. Mix by inversion. Incubate at 4oC (in the fridge) for 10 minutes.
- 29 29. Centrifuge at 12,000xg for 10 minutes.
- 30 30. Transfer supernatant (500 µl) to a new 1.5 ml microcentrifuge tube. Be careful not to disturb waste pellet.

- 31 31. Repeat steps 28-29.
- 32 32. Transfer supernatant to a new 1.5 ml microcentrifuge tube (with the green tough spot label). Be careful not to disturb waste pellet.

DNA Precipitation

- 33 33. Add 680 µl cold isopropanol (volume ratio ~1:1). Mix by inverting gently ~50 times. Centrifuge at 8000 xg for 5 minutes.
- 34 34. Carefully discard the supernatant, by pouring out. Keep an eye on DNA pellet to avoid avoiding losing it.
- 35 35. Pulse centrifuge the tubes and pipette out remaining supernatant while avoiding contact with the DNA pellet. DNA Washing

DNA Washing

- 36 36. Add 800 µl 70 % ethanol. Invert the tube several times to wash the DNA pellet. Centrifuge at 8000 xg for 1 minute.
- 37 37. Carefully discard the supernatant, by pouring out. Keep an eye on DNA pellet to avoid avoiding losing it.
- 38 38. Repeat steps 36-37.
- 39 39. Pulse centrifuge the tubes and pipette out remaining ethanol while avoiding contact with the DNA pellet.
- 40 40. Allow pellet to air dry for 20 minutes. Avoid over-drying the DNA pellet, as the DNA will be difficult to dissolve. DNA Rehydration

DNA Rehydration

- 41 41. Re-suspend dried pellets with 150 μ l TE or AE buffer. Invert tube to mix and spin down using a centrifuge.
- 42 42. Incubate at room temperature for 60 minutes or incubate in the fridge overnight.
- 43 43. Pulse centrifuge the tubes and store at -20oC.