



Sep 06, 2022

Native-PAGE analysis of VCP hexamer

Itika Saha¹, F. Ulrich Hartl¹, Mark S. Hipp^{2,3}¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry, Am Klopferspitz 18, 82152 Martinsried, Germany;²Department of Biomedical Sciences of Cells and Systems, University Medical Center Groningen, University of Groningen, Antonius Deusinglaan, 1, 9713 AV Groningen, The Netherlands;³School of Medicine and Health Sciences, Carl von Ossietzky University Oldenburg, 26129 Oldenburg, Germany

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.36wgqjrzyvk5/v1

Felix Kraus

ABSTRACT

Valosin-containing protein (VCP) is a homo-hexameric AAA+ ATPase in eukaryotic cells. This protocol describes the analysis of myc-tagged versions of VCP transiently transfected in HEK293 cells (stably expressing and propagating aggregates of Tau repeat domain fused to YFP) for hexamer formation.

DOI

dx.doi.org/10.17504/protocols.io.36wgqjrzyvk5/v1

EXTERNAL LINK

<https://www.biorxiv.org/content/10.1101/2022.02.18.481043v1.full>

PROTOCOL CITATION

Itika Saha, F. Ulrich Hartl, Mark S. Hipp 2022. Native-PAGE analysis of VCP hexamer. **protocols.io**
<https://protocols.io/view/native-page-analysis-of-vcp-hexamer-cf7wtrpe>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

The AAA+ chaperone VCP disaggregates Tau fibrils and generates aggregate seeds Itika Saha, Patricia Yuste-Checa, Miguel Da Silva Padilha, Qiang Guo, Roman Körner, Hauke Holthausen, Victoria A. Trinkaus, Irina Dudanova, Rubén Fernández-Busnadiego, Wolfgang Baumeister, David W. Sanders, Saurabh Gautam, Marc I. Diamond, F. Ulrich Hartl, Mark S. Hipp bioRxiv 2022.02.18.481043; doi: <https://doi.org/10.1101/2022.02.18.481043>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 05, 2022

LAST MODIFIED

Sep 06, 2022

PROTOCOL INTEGER ID

69590

- 1 Plate 1.5×10^5 cells in 12-well plate.
- 2 Next day, transfect with plasmids expressing myc-tagged VCP variants (Saha et al. *BioRxiv*, 2022) using a standard transfection protocol.
- 3 Two days later, collect cells and lyse them in 50 μ L 0.5% Triton X-100/PBS supplemented with protease inhibitor cocktail (Roche) and DNase for 1 h on ice.
- 4 Centrifuge lysates at 10,000 x g for 2 min and collect supernatant. 2m
- 5 Determine protein concentration in the supernatant and normalize across all samples.
- 6 Add 2x native sample buffer (40 % glycerol, 240 mM Tris pH 6.8, 0.04 % bromophenol blue) to 40 μ g lysate.
- 7 Run samples on a Native PAGE gel (e.g. Novex Value 4 to 12% Tris-glycine gels (Thermo)) in ^{1h} 20 mM Tris 200 mM Glycine buffer at pH 8.4.

- 8 Transfer proteins to nitrocellulose membrane in standard Tris-glycine buffer, block in 5% low-fat dry milk for 1 h at room temperature (RT).^{2h}

NOTE: Nitrocellulose membranes produce less background than PVDF membranes with fluorescent secondary antibodies.

- 9 Dilute anti-myc (9E10) and anti-VCP (1:2000, Novus Biologicals) primary antibodies together in blocking solution and incubate membrane overnight.
- 10 Next day, wash membrane 3 times with TBST and incubate with anti-mouse (LI-COR Biosciences Cat# 926-68070, RRID:AB_10956588; 1:10,000) and anti-rabbit (LI-COR Biosciences Cat# 926-32211, RRID:AB_621843; 1:10,000) fluorescent secondary antibodies for 2 h at RT.^{2h}
- 11 Wash membrane 3 times with TBST.
- 12 Detect fluorescent myc and VCP signal on a fluorescent imager.