



Mar 25, 2021

Human Lung Digestion: Deriving a single cell suspension

Morrisey Lab¹¹University of Pennsylvania

1

Works for me

This protocol is published without a DOI.

LungMap2 Consortium

Tech. support email: lungmap2dcc@gmail.com[Click here to message tech. support](#)

Chelseahortman

SUBMIT TO PLOS ONE

ABSTRACT

Human Lung Digestion: Deriving a single cell suspension

Digestion solution: (add PBS (without Ca,Mg) up to 50mL)

- 84mg Collagenase Type I (Life Tech, cat# 17100017, powder 285U/mg)
- 5mL Dispase (Fisher, cat# CB40235, 50U/mL)
- 3mL DNase1 (Roche, cat# 10104159001/165U/ml)

FACS Buffer:

- 1L dPBS (Cell Center, cat# 14190136, without Ca++ & Mg++)
- 2mL EDTA (Invitrogen, cat# 15575020, UltraPure, 0.5M, pH 8.0, 100mL)
- 10mL of 10%FBS (Invitrogen, cat# 10437-028, 500mL)

0.04% PBS+BSA (resuspending sample for sequencing)

- 10mg BSA (Jackson, cat# 001-000-162, IgG free, protease free)
- 25mL PBS (Cell Center, cat# 14190136, without Ca++ & Mg++)

MACS Accessories: (Miltenyi Biotec)

- GentleMACS Octo Dissociator w/ heaters (cat# 130-096-427)
- GentleMACS C Tubes (cat# 130-096-334)
- MACS- Smart Strainers 70uM (4x25pk) (cat# 130-110-916)

ATTACHMENTS

Morrissey Human Lung
Digestion Protocol March
.pdf

PROTOCOL CITATION

Morrissey Lab 2021. Human Lung Digestion: Deriving a single cell suspension. **protocols.io**
<https://protocols.io/view/human-lung-digestion-deriving-a-single-cell-suspen-btntnmen>



LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Mar 25, 2021

LAST MODIFIED

Mar 25, 2021

PROTOCOL INTEGER ID

48563