

Jul 07, 2022

Deep Dye Drop Protocol

Chiara Victor¹, Ben Gaudio¹, Mirra Chung¹, Mario Niepel¹, Marc Hafner¹, Luca Gerosa¹, Clarence Yapp¹, Kartik Subramanian¹, Peter Sorger², Caitlin Mills², Ajit Johnson Nirmal², Nicholas Clark²

¹Harvard Medical School; ²Harvard University

1 Works for me

Share

This protocol is published without a DOI.

Laboratory of Systems Pharmacology

Caitlin Mills

ABSTRACT

High-throughput measurement of cells perturbed using libraries of small molecules, gene knockouts, or different microenvironmental factors is a key step in functional genomics and pre-clinical drug discovery. However, it remains difficult to perform accurate single-cell assays in 384-well plates, limiting many studies to well-average measurements (e.g. CellTiter-Glo®). Here, a public domain “Dye Drop” method that uses sequential density displacement and microscopy to perform cell count and viability assays is described. Cell viability and DNA replication assays are followed by immunofluorescence imaging to collect single-cell dose-response data in the “Deep Dye Drop” version of the protocol. The resultant data can be used to calculate growth rate inhibition (GR) values and metrics. Dye Drop is rapid, reproducible, customizable, and compatible with manual or automated laboratory equipment. Dye Drop improves the tradeoff between data content and cost, enabling the collection of information-rich perturbation-response datasets.

PROTOCOL CITATION

Chiara Victor, Ben Gaudio, Mirra Chung, Mario Niepel, Marc Hafner, Luca Gerosa, Clarence Yapp, Kartik Subramanian, Peter Sorger, Caitlin Mills, Ajit Johnson Nirmal, Nicholas Clark 2022. Deep Dye Drop Protocol. **protocols.io**
<https://protocols.io/view/deep-dye-drop-protocol-96zh9f6>

FUNDERS ACKNOWLEDGEMENT

NIH
Grant ID: U54-CA225088

NIH
Grant ID: U54-HL127365

NIH
Grant ID: U24-DK116204

KEYWORDS

High-throughput microscopy, live-cell assays, dose response, small molecule perturbation, viability, cell cycle, GR metrics

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 05, 2019

LAST MODIFIED

Jul 07, 2022

OWNERSHIP HISTORY

Dec 05, 2019  Chiara Victor Harvard Medical School


Jul 06, 2022  Caitlin Mills

PROTOCOL INTEGER ID

30649

MATERIALS TEXT

MATERIALS

 [Optiprep \(Iodixanol\)](#) Sigma

Aldrich Catalog #D1556-250ML

 [Odyssey® Blocking Buffer \(PBS\)](#) LI-

COR Catalog #927-40000 927-40100

 [Microseal® 'F' Foil](#) BioRad

Sciences Catalog #MSF-1001

 [Hoechst 33342](#) Contributed by

users Catalog #H3570

 [Copper \(II\) sulfate pentahydrate](#) Sigma –

Aldrich Catalog #209198

 [Triton X-100](#) Sigma

Aldrich Catalog #X100

 [LIVE/DEAD™ Fixable Far Red Dead Cell Stain Kit, for 633 or 635 nm excitation](#) Thermo

Fisher Catalog #L34974

 [1X PBS](#) VWR

Scientific Catalog #75800-986

 [Formaldehyde solution 37%](#) Sigma

Aldrich Catalog #F1635-500ML

 [EdU Lumiprobe](#) Catalog #10540

 [Sulfo-Cy3](#)

azide Lumiprobe Catalog #B1330

 [Phospho-Histone H3 Alexa 488](#) Cell Signaling

Technology Catalog #3465S

 [Ascorbic acid](#)

Sigma Catalog #A4544

Stock solutions

- LIVE/DEAD Red- prepare according to directions (add  **50 µL** DMSO per tube), store at  **-20 °C** , limit

freeze/thaws

- EdU- prepare [M]**10 millimolar (mM)** working solution in DMSO, store at **-20 °C**
25 mg EdU
10 mL DMSO
- Sulfo-cy3 azide- prepare [M]**4 Molarity (m)** working solution in DMSO, store at **-20 °C**
3 mg sulfo-cy3-azide
1 mL DMSO
- Ascorbic acid- prepare **200 mg/ml** working solution in water ****fresh each time!****
- CuSO₄·5H₂O- prepare [M]**200 millimolar (mM)** in water, store at room temperature, protected from light
250 mg CuSO₄S
5 mL water
- Triton X-100- prepare 10% working solution in 1X PBS
10 mL Triton X-100
90 mL PBS

1 Pulse cells with EdU+ stain dead cells with LDR

For 10 ml:

- 1 mL** optiprep (10% final)
- 9 mL** PBS
- 5 µL** LDR (1:2000 final)
- 10 µL** EdU (10 µM final)

384-well plate: Add **15 µL** per well along the edge of the wells using a multi-channel pipette and incubate for **01:00:00** (or desired pulse duration) @ **37 °C**

96-well Plate: add **60 µL** per well.

1.1

This file can be used for all solution calculations:

 **DDD.xlsx**

2 Fix cells

For 10 ml:

- 2 mL** Optiprep (20% final)
- 6.9 mL** PBS
- 1.1 mL** formaldehyde (4% final)

384-well plate: Add **20 µL** per well along the edge of the wells using a multi-channel pipette and incubate for **00:30:00** @ **Room temperature** in the dark (cover with foil)

96-well Plate: add **80 µL** per well.

- 2.1 Run DeepDyeDrop1 protocol on the plate washer (aspirate all but **10 µL** , dilute with **80 µL** PBS, stop, store @ **4 °C** . Aspirate all but **10 µL** if continuing with permeabilization right away).

3 Permeabilize:

For 10 ml:

- 1 mL** Optiprep (10% final)
- 8.5 mL** PBS
- 500 µL** 10% Triton X-100 (0.5 % final)

384-well plate: Add **15 µL** per well along the edge of the wells using a multi-channel pipette and incubate for **00:20:00** @ **Room temperature** in the dark.

96-well plate: add **60 µL** per well.

4 Click Reaction:

For 10 ml:

Combine **in order**

- 7 mL** PBS
- 2 mL** optiprep (20% final)
- 100 µL** 200mM CuSO₄ (2mMfinal)
- 10 µL** 4mM sulfo-cy3-azide (4µM final)
- 1 mL** 200 mg/ml ascorbic acid (20 mg/ml final)

384-well plate: Add **20 µL** per well along the edge of the wells using a multi-channel pipette and incubate for **00:30:00** @ **Room temperature** in the dark.

96-well Plate: add **80 µL** per well.

- 4.1 Run DeepDyeDrop1 protocol on the plate washer (aspirate all but **10 µL** , dilute with **80 µL** PBS).

Important note: Stopping point if not performing immunofluorescence staining. Include **2 µL** Hoechst 33342 (1:5000 final) in the click reaction mix

5 Immunofluorescence:

384-well plate: Add **40 µL** Odyssey blocking buffer per well and incubate @ **Room temperature** for **01:00:00** in the dark on a plate rocker at the slowest setting.

96-well plate: Add **160 µL** per well.

- 5.1 Run DeepDyeDrop2 protocol on the plate washer (aspirate all but **10 µL**)

Antibodies:

5.2

For 10 ml:

▢ **10 mL** odyssey blocking buffer

▢ **2 µL** Hoechst 33342 (1:5000 final)

▢ **5 µL** pH3-A488 antibody (1:2000 final)

384-well plate: Add ▢ **15 µL** antibody solution per well, incubate 🕒 **Overnight** @ 🌡 **4 °C** in the dark, on a plate rocker at the slowest setting.

96-well plate: add ▢ **60 µL** per well.

5.3

Run DeepDyeDrop3 protocol on the plate washer (wash once with PBST, twice with PBS, leave ▢ **80 µL** PBS per well). Seal plates with Microseal 'F' foil seal and store @ 🌡 **4 °C** until imaging.

6 Image Acquisition:

On the Operetta:

Hoechst (ex 360-400 em 410-480)
EdU-cy3 (ex 520-550 em 560-630)
pH3-A488 (ex 460-490 em 500-550)
LDR-A647 (ex 620-640 em 650-700)

On the IXM-C:

DAPI
FITC
TRITC
Cy5

Important note: Focus height will vary with plate type and cell line, but expect to image the pH3 higher than the others since cells ball up during mitosis

7 Image Analysis

1. Apply flat field correction if available
2. Segment nuclei based on Hoechst signal
3. Define a ring around the nuclei, be sure this does not include any of the nuclear area (it will be used to subtract the local background from each channel)
4. calculate the **nuclear area**
5. for each channel, calculate the average intensities within the nuclear mask, and within the surrounding ring
6. **for each channel, subtract the ring intensity from the nuclear intensity**
7. multiply the nuclear area by the Hoechst intensity for the **DNA content**
8. calculate **SER spot pixel size 8 texture feature** for the LDR signal within the nuclear area
9. output features (nuclear area, DNA content, Hoechst, EdU, pH3, LDR background corrected intensities, LDR texture feature) in a .tsv file
10. identify cells in mitosis through pH3 intensity
11. identify dead cells through LDR intensity and/or texture- threshold independent of cell line
12. identify additional dead cells through size (small) and intensity (high)- optimize per cell line
13. identify population of cells in each phase of the cell cycle based on DNA content and EdU intensity
14. identify additional dead cells based on size (small) and intensity (high)

