



Sep 05, 2022

Integra Total Nucleic Acid Extraction

Sophana Chea¹, Sreyngim Lay¹, Mengheng Oum¹, Gechlang Tang¹,
Cheata Hou¹, Manu Vanaerschot², Christina Yek³, Jessica Manning³,
Vida Ahyong²

¹International Center of Excellence in Research, National Institute of Health, Cambodia;

²Chan Zuckerberg Biohub; ³National Institute of Health

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.81wgbydq1vpk/v1



Sophana Chea

International Center of Excellence in Research, National Ins...

ABSTRACT

This SOP is the process of extracting Total Nucleic Acid (TNA) from Sera and/or Nasopharyngeal or Nasal swabs. The isolated high-quality nucleic acids are suitable for Next-Generation Sequencing (NGS).

DOI

dx.doi.org/10.17504/protocols.io.81wgbydq1vpk/v1

PROTOCOL CITATION

Sophana Chea, Sreyngim Lay, Mengheng Oum, Gechlang Tang, Cheata Hou, Manu Vanaerschot, Christina Yek, Jessica Manning, Vida Ahyong 2022. Integra Total Nucleic Acid Extraction. **protocols.io**
<https://protocols.io/view/integra-total-nucleic-acid-extraction-cf2ztqf6>



KEYWORDS

Integra, Total Nucleic Acid, DNA, RNA

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 01, 2022

LAST MODIFIED

Sep 05, 2022

PROTOCOL INTEGER ID

69433

GUIDELINES

Adapted from the Quick-DNA/RNA Pathogen MagBead Kit Manual (Zymo Research, Cat# R2145)

MATERIALS TEXT

1. RNase away sprays for RNase decontaminants (Thermo Scientific, Cat#7002).

[🔗RNase Away spray](#) **Molecular BioProducts**

2. 1ml deep well sterile plate.

3. 2ml deep well sterile plate.

4. Hard-shell PCR Plates 96-well (Bio-Rad, Cat# HSP9601).

5. PCR Plate Seal, foil (Bio-Rad, Cat# MSF1001).

6. Quick-DNA/RNA Pathogen MagBead kit (Zymo Research, Cat# R2145)

[🔗Quick-DNA/RNA Pathogen MagBead kit](#) **Zymo**

Research Catalog #R2145

7. Molecular-Grade Isopropanol (100% isopropanol).

8. Molecular-Grade Absolute ethanol (100% ethanol).

9. Proteinase K w/ Storage buffer 20mg (Zymo Research, Cat# D3001-2-20).

[🔗Proteinase K w/ Storage buffer 20mg set](#) **Zymo**

Research Catalog #D3001-2-20

10. DNase/RNase-free water

[🔗Nuclease-free water](#)

Ambion Catalog #AM9932

11.

VIAFLO

96 channel pipette

Integra

VIAFLO 96



12. DNase I (Zymo Research, Cat# E1010)

[🔗DNase I set](#) **Zymo**

Research Catalog #E1010

13. 96-well magnetic stand

SAFETY WARNINGS

All steps should be performed at **🔗 Room temperature**.

***NOTE:** When reusing tips, make sure to include a bit of extra air aspiration to avoid drops at the bottom of tips when aspirating volumes, and also a bit extra air blows out

at the end of dispensing steps in plates.

Perform the TNA extraction in the extraction room separate from the PCR room.
Respect the Laboratory safety guideline for all steps of the protocol.
Wearing PPE is recommended.

Reagent preparation (required with new kit) 15m

1. Add **20 mL** of isopropanol to the MagBead DNA/RNA Wash 1 concentrate.
 2. Add **30 mL** of isopropanol to the MagBead DNA/RNA Wash 2 concentrate.
 3. Add **1.2 mL** of Proteinase K Storage Buffer per vial to reconstitute the lyophilized Proteinase K at **20 mg/mL**.
- Vortex to dissolve. STORE AT **-20 °C** Freezer.

Preparation the buffer plate (before starting protocol) 1h

1. Pre-make pathogen buffer plate with **880 µL** Pathogen DNA/RNA buffer in 1ml deep well plate.
2. Pre-make bead plate with **25 µL** MagBinding beads into 96V-well PCR plate.
*Make immediately before starting, <1h prior to starting the protocol, to ensure the beads are mixed.
3. Pre-make DNA/RNA Wash 1 plate with **550 µL** Wash 1 buffer into a 1ml deep well plate.
4. Pre-make DNA/RNA Wash 2 plate with **550 µL** Wash 2 buffer into a 1ml deep well plate.
5. Pre-make 100% ethanol plate with **1 mL** 100% ethanol into a 2ml deep well plate.
6. Pre-make 80% ethanol plate with **600 µL** 80% ethanol into a 1ml deep well plate.
7. Pre-make water plate with **50 µL** DNAase/RNase-free water in a 96 V-well PCR plate.
8. Pre-make water plate with **30 µL** DNAase/RNase-free water in a 96 V-well PCR plate.

Spin all plates down for **00:01:00** except the bead plate. Perform a quick pulse spin down of the bead plate, just enough to get all liquid down. Centrifuge the rest of the plates at 12 000 rpm for **00:01:00**.

Sample preparation and Proteinase K 16m

3. 1. Create a plate map so you know which sample you are adding to each well. Add **400 µL** of your samples to 2ml deep well plate (Plate 1).
2. Manually add **65 µL** of Proteinase K to each well of 8 well PCR strip tubes. Using a manual multichannel pipet, aliquot **4 µL** of Proteinase K into each sample (Plate 1).
3. Load a set of Integra tips (tip set 1) onto the Integra.
4. **Program: Pipet/Mix 250ul, 15 cycles, speed 10.** Program the Integra to pipet

▢ **250 µL** of your samples up and down for ⌚ **00:01:00** (15 cycles), then incubate at
🌡 **Room temperature** for ⌚ **00:15:00** .

*Keep tips.

Sample binding and washing

1h 20m

1h 20m

- 4 5. **Program: Pipet 300ul.** add ▢ **800 µL** total of Pathogen DNA/RNA Buffer to the sample plate (Plate 1).
6. **Program: Pipet/Mix 250ul, 30 cycles, speed 10.** Program the Integra to mix samples and buffer for ⌚ **00:02:00** . Keep tips.
7. **Program: Pipet/Mix 20ul, 20 cycles.** Program the Integra to mix the MagBinding beads plate so the beads are fully resuspended. If beads are not fully to the bottom of the plate, perform a short pulse spin.
8. **Program: Pipet 20ul.** Program the Integra to aspirate ▢ **20 µL** from the MagBinding beads plate. Check that all wells have beads!
9. **Program: Pipet/Mix 250ul, 30 cycles x 5 (10 min total), speed 7.** Program the Integra to mix the beads with the sample for ⌚ **00:10:00** total. Keep tips.
10. Place the 96-well magnetic stand underneath the sample plate for ⌚ **00:05:00** until a bead ring forms.
11. **Program: Manual Pipet 300ul.** Aspirate and discard the cleared supernatant into a 2ml deep well waste plate.
Discard tips, load new tips.
12. Remove magnetic stand.
13. **Program: Pipet 250ul.** Dispense a total of ▢ **500 µL** Wash 1 into the sample plate.
Keep tips.
14. **Program: Pipet/Mix 250ul, 30 cycles, speed 7.** Program the Integra to mix the beads with the Wash buffer for ⌚ **00:02:00** total. Keep tips.
15. Place the 96-well magnetic stand underneath the sample plate for ⌚ **00:02:00** until bead ring forms.
16. **Program: Manual Pipet 300ul.** Aspirate and discard the cleared supernatant into the 2ml deep well waste plate. Keep tips.
17. **Program: Pipet 250ul.** Dispense a total of ▢ **500 µL** Wash 2 into the sample plate.
Keep tips.
18. **Program: Pipet/Mix 250ul, 30 cycles, speed 7.** Program the Integra to mix the beads with the wash buffer for ⌚ **00:02:00** total. Keep tips.
19. Place the 96-well magnetic stand underneath the sample plate for ⌚ **00:02:00** until a bead ring forms.
20. **Program: Manual Pipet 300ul.** Aspirate and discard the cleared supernatant into the 2ml deep well waste plate. Discard tips, load new tips.
21. Remove magnetic stand.
22. **Program: Pipet 250ul.** Dispense a total of ▢ **500 µL** 100% Ethanol into the sample

plate.

23. **Program: Pipet/Mix 250ul, 30 cycles, speed 7.** Program the Integra to mix the beads with Ethanol for 🕒00:02:00 total. Keep tips.

24. Place the 96-well magnetic stand underneath the sample plate for 🕒00:02:00 until a bead ring forms.

25. **Program: Manual Pipet 300ul.** Aspirate and discard the cleared supernatant into the 2ml deep well waste plate.

26. Remove magnetic stand.

27. Repeat steps 22 to 26 for another round of Ethanol wash. Discard tips after 2nd Ethanol wash. Load new tips.

28. Keep the sample plate on the magnetic stand after the final Ethanol wash until the beads are dry (~ 🕒00:10:00).

29. **Program: Pipet/Mix 33ul, 30 cycles, speed 7.** Pipet 📏33 µL of Nuclease-free water (Heat the water at 🔥55 °C for 🕒00:10:00 before elution) into the dried beads and mix the beads for 🕒00:02:00 total. Keep tips.

30. Place the 96-well magnetic stand underneath the sample plate for 🕒00:02:00 until a bead ring forms.

31. **Program: Manual Pipet 30ul.** Pipet 📏30 µL from the sample plate and dispense into a new 96 V-bottom PCR plate.

32. Store TNA sample immediately at 🔥-80 °C .

DNase treatment post-TNA purification

15m

- 5
1. Aliquot SPRI beads (well mixed and resuspended) 📏50 µL into a new 96 V-bottom PCR plate.^{15m}
 2. Prepare DNase I solution
 - a. Add 📏275 µL DNase/RNase-free water to reconstitute lyophilized DNase I and mix.
 3. Aliquot 📏20 µL of TNA into a new 96 V-bottom PCR plate.
 4. Prepare DNase master mix:
 - a. 1x rxn = 20ul sample + 2.5ul DNase I + 2.5ul Digestion Buffer.
 - b. Using a manual multichannel pipet, aliquot 📏5 µL of DNase I master mix to each sample and mix.
 - c. **Program: Pipet/Mix 20ul, 20 cycles, speed 7.**
 5. Incubate 🕒00:15:00 🔥Room temperature

SPRI bead clean-up

30m

- 6
1. **Program: Pipet/Mix 45ul, 20 cycles, speed 7.** Add 📏45 µL SPRI beads (1.8x ratio) to the DNase I treated sample and mix by pipetting.^{30m}
 2. Incubate 🕒00:05:00 🔥Room temperature .

3. Place the 96-well magnetic stand underneath the sample plate for ⌚00:02:00 until a bead ring forms.
4. **Program: Manual Pipet 100ul.** Aspirate and discard the cleared supernatant into the 2ml deep well waste plate. Keep tips.
5. **Program: Pipet 200ul. (1st Ethanol Wash).** Aspirate and dispense 📏200 µL of freshly made 80% Ethanol into sample plate.
6. Incubate for ⌚00:01:00 at 🌡 Room temperature .
7. **Program: Manual Pipet 200ul.** Aspirate and discard the cleared supernatant into the 2ml deep well waste plate. Keep tips.
8. **Program: Pipet 200ul. (2nd Ethanol Wash).** Aspirate and dispense 📏200 µL of freshly made 80% Ethanol into sample plate.
9. Incubate for ⌚00:01:00 at 🌡 Room temperature .
10. **Program: Manual Pipet 200ul.** Aspirate and discard the cleared supernatant into the 2ml deep well waste plate. Discard tips, load new tips.
11. Incubate for ⌚00:05:00 at 🌡 Room temperature until beads are dry.
12. Remove plate from magnetic stand.
13. **Program: Pipet/Mix 22ul, 10 cycles, speed 7.** Add 📏22 µL DNase/RNase-free water to each well and mix beads.
14. Place the 96-well magnetic stand underneath the sample plate for ⌚00:02:00 until a bead ring forms.
15. **Program: Manual Pipet 20ul.** Aspirate 📏20 µL of purified TNA sample to a new 96 V-bottom PCR plate.
16. Store sample at 🌡 -80 °C immediately.