



Jan 01, 2021

# Cell Harvesting Protocol

Elizabeth Fozo<sup>1</sup><sup>1</sup>In-house protocol

1

Works for me

This protocol is published without a DOI.

Eadewunm

## ABSTRACT

A general procedure for harvesting bacterial cells for flow cytometry or fluorescence microscopy  
Write up by S. Shore 11/17/2020 based on B. Bogati communications, proofread by Bikash Bogati

## PROTOCOL CITATION

Elizabeth Fozo 2021. Cell Harvesting Protocol. **protocols.io**  
<https://protocols.io/view/cell-harvesting-protocol-bq4wmyxe>

## KEYWORDS

cell harvesting, cell harvesting protocol, protocol, harvesting cells, cells, harvesting

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jan 01, 2021

## LAST MODIFIED

Jan 01, 2021

## PROTOCOL INTEGER ID

45942

## DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## ABSTRACT

A general procedure for harvesting bacterial cells for flow cytometry or fluorescence microscopy  
Write up by S. Shore 11/17/2020 based on B. Bogati communications, proofread by Bikash Bogati

## Procedure

- 1 Dilute the overnight culture to OD<sub>600</sub> 0.01 and grow cells in an appropriate medium to the desired time point/OD<sub>600</sub> in at least 10mL of liquid culture.

- 2 Harvest cells directly into pre-labeled falcon tubes (if running flow same-day) or Eppendorf tubes. Once cells have been harvested, minimize exposure to direct light if looking for fluorescence.
  1. For *E. coli* from overnight LB, at least 25uL
  2. For *E. coli* from exponential phase (OD~0.3) LB, at least 50uL
  3. For *E. coli* from the overnight or exponential phase in minimal media, at least 500uL
- 3 For membrane depolarization assay:
  1. Wash the cells twice with 1x PBS and resuspend in 1 ml PBS
  2. Add DiBAC<sub>4</sub>-3 reagent- final concentration 10mg/mL (4ml of 2.5 mg/mL working stock in 1 mL)
  3. Mix and incubate in dark for 20 minutes. Wash twice in PBS.
  4. Fix using paraformaldehyde and glutaraldehyde (see below)
- 4 Optional: Fix with desired fixation method in the chemical fume hood
  1. Required if using the LSR-II in Mossman. Special permission needed to run live-cells through a flow cytometer.
  2. Potential fixation methods:
    - 0.5% paraformaldehyde for 10 min in a cold room in dark (T. Hancock)
    - 3% glutaraldehyde for 1hr at room temp in dark (R. Johnston)
    - paraformaldehyde to final concentration 2.8% (vol/vol) and glutaraldehyde to 0.04% (vol/vol) for 15 min at room temp in dark (Dr. J.Mannik)
- 5 Wash cells 3 times
  1. Centrifuge at 4 degrees until pellets form (tabletop for 3 min, if in falcon tubes, 10 min in a large centrifuge)
  2. Decant supernatant in an appropriate chemical waste container in the hood
  3. Flood sample with PBS and vortex
  4. Repeat steps 1-3 2 additional times
- 6 Decant liquid and add PBS to the desired volume.
  1. For flow and *E. coli*, 1 mL for minimal media samples, 0.5 mL for LB samples if the followed recommendation in step 2.
  2. If mounting cells on slides for microscopy, do not add PBS during the last round. Aspirate remaining liquid out instead.
- 7 If applicable, store cells in the fridge (4 degrees).
  1. If using the J.Mannik fixation method, cells can be stored up to a month in the dark.