



AUG 01, 2023

OPEN  ACCESS

Protocol Citation: Juliana Alcoforado Diniz, Paul Robson, Dylan Baker, Jessica Garofalo 2023. KAPP-Sen TMC: Dissociation of Pancreatic Islets (non-recovered) Protocol.

protocols.io

<https://protocols.io/view/kapp-sen-tmc-dissociation-of-pancreatic-islets-non-cxz5xp86>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Aug 01, 2023

Last Modified: Aug 01, 2023

KAPP-Sen TMC: Dissociation of Pancreatic Islets (non-recovered) Protocol

Juliana Alcoforado Diniz¹, Paul Robson², Dylan Baker²,
Jessica
Garofalo²

¹The Jackson Laboratory; ²The Jackson Laboratory for Genomic Medicine

Cellular Senescence Network (SenNet) Method Development
Community



Ashley M Raynock

UConn Health, UConn Center on Aging

ABSTRACT

The dispersed samples were shipped cold from **PRODOLABS**.

Prior to scRNA-seq, dispersed samples from brain dead donor's pancreatic acinar and ducts were dissociated as follows.

Abstract

- 1 The dispersed samples were shipped cold from **PRODOLABS**. Prior to scRNA-seq, dispersed samples from brain dead donor's pancreatic acinar and ducts were dissociated as follows.

Cell Dissociation with TrypLE

- 2

NOTE: Before beginning cell dissociation coat all the materials (pipettes, tubes, etc.) with stopping media to prevent sticking. Stopping Media: DMEM + 10% FBS + 1:100 Glutamate.

 1. Distribute specimens into 50 ml conical tubes, try to make similar cell concentrations across all tubes. *** Receive acinar in ~300 ml, should distribute into 12 tubes with 25 ml each***
 2. Centrifuge acinar samples for 2 min at 230 g and 4C. Leave duct samples to settle without disturbance.
 3. Wash samples 3X with 25 ml PBS.
 4. Add 10 ml TrypLE E and incubate at 37C, mixing with a pipette approximately every 5 min.
 - i. Acinar: in TrypLE E at 37C for approximately 30 mins.
 - ii. Ducts: in TrypLE E at 37C for approximately 15 mins.
 5. To stop reaction, add 30 ml per tube of stopping media.
 6. Centrifuge acinar and ducts at 1300 rpm, 4C, 2 mins.
 7. To stop: add 30 ml per tube of stopping media.
 8. Centrifuge acinar and ducts at 1300 rpm, 4C, 2 mins.
 9. Aspirate TrypLE E and wash 1X with stop media. Aspirate then resuspend in stop media to (Acinar) 5 ml or (Ducts) 2 ml.
 10. Transfer cell suspension through 70 µm filter then combine to 2 ml (Acinar) and 1ml (Ducts). Filter again through a 40 µm Flowmi.
 11. Count cells using AO/PI (acridine orange/propidium iodide) Cell Viability Kit for Luna-FL automated cell counter.
 12. Proceed to cell fixation.

Fixation of Cells & Nuclei for Chromium Fixed RNA Profiling

- 3 Cells were fixated prior to scRNAseq according to [https://dx.doi.org/10.17504/protocols.io.\[...\]/v1](https://dx.doi.org/10.17504/protocols.io.[...]/v1)