

May 05, 2020

# Separation and purification of human PBMC from FRESH BLOOD

Marco Cosentino<sup>1</sup>, Elisa Storelli<sup>2</sup>, Alessandra Luini<sup>1</sup>, Massimiliano Legnaro<sup>3</sup>, Emanuela Rasini<sup>3</sup>, Marco Ferrari<sup>3</sup>, Franca Marino<sup>3</sup>

<sup>1</sup>Center for Research in Medical Pharmacology, University of Insubria (Varese, Italy),

<sup>2</sup>Center for Research in Medical Pharmacology, University of Insubria, <sup>3</sup>University of Insubria

**1** Works for me [dx.doi.org/10.17504/protocols.io.bfxmjpk6](https://dx.doi.org/10.17504/protocols.io.bfxmjpk6)



Elisa Storelli

Center for Research in Medical Pharmacology, University of I...

## ABSTRACT

Separation and purification of PBMC from FRESH BLOOD: list of published work using this protocol

Kustrimovic, N., Comi, C., Magistrelli, L., Rasini, E., Legnaro, M., Bombelli, R., Aleksic, I., Blandini, F., Minafra, B., Riboldazzi, G., Sturchio, A., Mauri, M., Bono, G., Marino, F., & Cosentino, M. (2018). Parkinson's disease patients have a complex phenotypic and functional Th1 bias: cross-sectional studies of CD4+ Th1/Th2/T17 and Treg in drug-naïve and drug-treated patients. *Journal of neuroinflammation*, 15(1), 205. <https://doi.org/10.1186/s12974-018-1248-8>

Kustrimovic, N., Rasini, E., Legnaro, M., Bombelli, R., Aleksic, I., Blandini, F., Comi, C., Mauri, M., Minafra, B., Riboldazzi, G., Sanchez-Guajardo, V., Marino, F., & Cosentino, M. (2016). Dopaminergic Receptors on CD4+ T Naïve and Memory Lymphocytes Correlate with Motor Impairment in Patients with Parkinson's Disease. *Scientific reports*, 6, 33738. <https://doi.org/10.1038/srep33738>

Cosentino M., Ferrari M., Kustrimovic N., Rasini E., Marino F. (2015). Influence of dopamine receptor gene polymorphisms on circulating T lymphocytes: A pilot study in healthy subjects. *Human immunology*, 76, 10, 747-752. <https://doi.org/10.1016/j.humimm.2015.09.032>

## MATERIALS

NAME	CATALOG #	VENDOR
Ficoll Paque PLUS	17144003-500 ml	Ge Healthcare
Fetal Bovine Serum (FBS)	ECS0180L-500 ml	EuroClone
RPMI 1640	ECM 0495L- 500 ml	EuroClone
Trypan Blue solution 0.4%	T8154- 100 ml	Sigma Aldrich

## MATERIALS TEXT

### Instrumentation required:

- Laminar flow hood
- Centrifuge
- Cellometer (automated cell counter) or Optical Microscope (manual cell count)
- Flow Cytometer
- Autoclave

## EQUIPMENT

NAME	CATALOG #	VENDOR
Allegra AVANTI 30	Beckman Italy	Beckman Coulter
BD FACS Celesta	Milan Italy BD	
Cellometer Auto T4	Euroclone	

#### BEFORE STARTING

If you need to obtain **PBMC for cell culture**, make sure you are using **sterile PBS, culture medium, filtered Lysis Buffer and sterile plastic disposables as well**. Moreover, work under laminar flow hood when you are processing samples. Otherwise, use non-sterile solutions and plastic disposables, and process samples in cell isolation laboratory.

#### ALL REAGENTS USED IN THIS PROTOCOL MUST BE AT ROOM TEMPERATURE!

1 Put the needed amount of blood sample into a **50 ml** conical tube.

2 Add an equal volume of **PBS 1X** and mix well.




PBS 1X  
by Marco Ferrari,  
University of Insubria

3 Place **3 ml** of FICOLL in a **15 ml** conical tube.



4 Carefully layer **12 ml** of diluted blood on FICOLL with a glass Pasteur Pipette to a final volume of 15 ml as shown in the figure below.




5 Centrifuge samples **400 x g 00:40:00** at room temperature (RT) without break.



Allegra AVANTI 30  
Centrifuge  
Beckman Coulter Beckman Italy

- 6 After centrifugation, take out the tubes carefully to not disturb the mononuclear cell layer that appears as a white, cloudy band between the plasma and FICOLL as shown in the figure below.





- 7 Carefully with a glass Pasteur pipette transfer the mononuclear lymphocyte cell layer to another 15 ml conical tube.

- 8 Wash the isolated PBMC with **PBS/FBS 2%** to a final volume of **10 ml** and centrifuge at **600 x g 00:10:00** at RT.







RPMI - FBS and PBS - FBS  
by Elisa Storelli,  
Center for Research in Medical Pharmacology, University of Insubria

- 9 Remove supernatants, resuspend pellet in **1 ml** of **Lysis Buffer** and add another **9 ml** of **Lysis Buffer**. Immediately centrifuge tubes at **300 x g 00:10:00** at RT.





Lysis Buffer  
by Elisa Storelli,  
Center for Research in Medical Pharmacology, University of Insubria

- 10 Remove supernatant and resuspend pellet in  **10 ml** of **PBS/FBS 2%** and centrifuge at  **600 x g 00:10:00** at RT.




RPMI - FBS and PBS - FBS  
by Elisa Storelli,  
Center for Research in Medical Pharmacology, University of Insubria


- 11 Remove supernatant and resuspend the obtained pellet in  **10 ml** of **RPMI/FBS 10%** for cell counting.



RPMI - FBS and PBS - FBS  
by Elisa Storelli,  
Center for Research in Medical Pharmacology, University of Insubria

- 12 **For manual cell count use Türk solution for checking purity.**

Mix  **10 µl** of cell suspension with an equal amount of Türk solution (dilution factor = 2), allow mixture 3 min at room temperature.

Take  **10 µl** of the mixture and place it inside a Bürker chamber and view under an optical microscope using 40X magnification.

Count the cells in each square found in the four corners and in the central square (see figure 1 below), including those that lie on the bottom and left-hand perimeters, but not those that lie on the top and right hand perimeters (see figure 2 below).

**Total number of cells per ml = mean number of cells x dilution factor x 10<sup>4</sup> (hemacytometer volume).**




Figure 1




Figure 2




Trypan blue and Turk solution  
by Marco Ferrari,  
University of Insubria




## 13 OPTIONAL STEP

**For automatic cell count with Cellometer machine use Trypan Blue.** The machine will calculate the n° of cells/ml and the % of viability.

Take **10  $\mu$ l** of cell suspension and add an equal amount of Trypan Blue. Use all the volume to place it in a counting chamber. Place the chamber inside Cellometer and count.




Cellometer Auto T4  
Automated Cell Counter  
Nexcelom Bioscience Euroclone



Trypan blue and Turk solution  
by **Marco Ferrari**,  
University of Insubria


- 14** If needed, check the purity of PBMC suspension by using morphological parameter of the flow cytometer.

For this test  $0.5 \times 10^6$  PBMC in 500  $\mu$ l of PBS are enough.



BD FACS Celesta  
Flow Cytometer  
Becton Dickinson Milan Italy BD

## 15 Expected results



**VIABILITY** - The expected viability by Trypan Blue should be  $\geq 90\%$ .

**PURITY** - The PBMC suspension obtained should contain at least 80% of lymphocytes, 10-15% of monocytes and few contaminant PMN cells ( $\leq 5\%$ ) as confirmed by flow cytometry.

**YIELD** - The expected amount of PBMCs should be  $\pm 28,5 \times 10^6$  starting from 25 ml of fresh blood.