




Oct 27, 2020

Protein expression in *Bacillus subtilis*

Kristoffer Bach Falkenberg¹, Cristina Hernandez Rollan¹, Maja Rennig¹, Andreas Birk Bertelsen¹, Morten Norholm¹

¹Technical University of Denmark

1 Works for me dx.doi.org/10.17504/protocols.io.bdmui46w

 Kristoffer Bach Falkenberg
Technical University of Denmark

ABSTRACT

B. subtilis is a gram-positive bacteria used by both academia and industry as a protein production workhorse. This is due to its' their excellent fermentation properties, high production titers, and capacity to secrete proteins into the extracellular medium.

This protocol describes how to express proteins in *B. subtilis*. The protocol is developed using KO7-S, although it might also work for other strains as well. The method is adapted from Rasmussen, M. D.; Bjoernvad, M. E.; Diers, I. Pectate Lyase Fusion for Expression and Secretion of Polypeptides. WO 00/75344, 2000 and Jensen, K.; Østergaard, P. R.; Wilting, R.; Lassen, S. F. Identification and Characterization of a Bacterial Glutamic Peptidase. *BMC Biochem* **2010**, 11 (1), 47. <https://doi.org/10.1186/1471-2091-11-47>.

DOI

dx.doi.org/10.17504/protocols.io.bdmui46w

PROTOCOL CITATION

Kristoffer Bach Falkenberg, Cristina Hernandez Rollan, Maja Rennig, Andreas Birk Bertelsen, Morten Norholm 2020. Protein expression in *Bacillus subtilis*. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bdmui46w>

KEYWORDS

B. subtilis, Bacillus, *Bacillus subtilis*, Protein expression

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Mar 13, 2020

LAST MODIFIED

Oct 27, 2020

PROTOCOL INTEGER ID

34196

MATERIALS TEXT

MATERIALS

 [Sodium molybdate dihydrate](#) Contributed by users

 [manganese sulfate](#) Contributed by users

 [Iron\(III\) chloride hexahydrate](#) Contributed by

users Catalog #236489

 [Zinc sulfate heptahydrate](#) Sigma

Aldrich Catalog #204986

 [Magnesium sulfate heptahydrate](#) Contributed by users

 [Sodium Phosphate dibasic](#) Fisher

Scientific Catalog #S373-500

 [Copper \(II\) sulfate pentahydrate](#) Sigma –

Aldrich Catalog #209198

 [Yeast extract](#) Contributed by users

 [Nalgene™ Rapid-Flow™ Sterile Disposable Bottle Top Filters with PES Membrane, 150mL, 0.45μm pore, 45mm neck](#) Thermo

Fisher Catalog #296-4545

 [Maltodextrin \(DE 13.0-17.0\)](#) Sigma

Aldrich Catalog #419680

 [Pluronic L-61](#) Sigma

Aldrich Catalog #435422

SAFETY WARNINGS

Be sure to wear protective equipment when adjusting the pH of the media. Follow local safety regulations

ABSTRACT

B. subtilis is a gram-positive bacteria used by both academia and industry as a protein production workhorse. This is due to its' their excellent fermentation properties, high production titers, and capacity to secrete proteins into the extracellular medium.

This protocol describes how to express proteins in *B. subtilis*. The protocol is developed using KO7-S, although it might also work for other strains as well. The method is adapted from Rasmussen, M. D.; Bjoernvad, M. E.; Diers, I. Pectate Lyase Fusion for Expression and Secretion of Polypeptides. WO 00/75344, 2000 and Jensen, K.; Østergaard, P. R.; Wiltng, R.; Lassen, S. F. Identification and Characterization of a Bacterial Glutamic Peptidase. *BMC Biochem* **2010**, 11 (1), 47. <https://doi.org/10.1186/1471-2091-11-47>.

BEFORE STARTING

Make sure you have your expression strain freshly streaked on an agar plate.

Cal18-2 media preparations

- 1 Prepare a stock solution of 2.0g/L Na₂MoO₄. Sterilize by filtration
- 2 Prepare a trace metal solution consisting of
 - 4.48g/L MnSO₄·H₂O
 - 3.33g/L FeCl₃·6H₂O
 - 0.625g/L CuSO₄·5H₂O
 - 7.12g/L ZnSO₄·7H₂OSterilize by filtration
- 3 Fill a blue cap bottle to ~80% of the desired final volume with MQ water.
- 4 Add a magnetic stirrer to the blue cap bottle and place the bottle on a stirring plate. Turn on the stirring, and make sure it's mixing well.
- 5 Add the following to the bluecap bottle:
 - 40g/L yeast extract
 - 1.3 g/L MgSO₄·7H₂O
 - 50 g/L maltodextrin (DE ~ 12)
 - 20 g/L NaH₂PO₄·2H₂O
 - 6.7mL/L 2.0g/L Na₂MoO₄ stock solution
 - 6.7mL/L Trace metal solution
 - 100μL/L Pluronic L-61
- 6 Make sure that all of the ingredients are dissolved
- 7 Adjust to pH 6 with 5M NaOH
- 8 Add MQ water to the desired final volume
- 9 Sterilize by filtration



The media easily clogs filters, so choose a 0.45μM vacuum bottle top filter for this step and be prepared to

use a few filters per liter

- 10 Store the media at **4 °C** until needed

Overnight culture - Day 1

- 11 Inoculate between **3 mL** to **50 mL** LB media with a single colony of the expression strain. Depending on the expression volume and overnight OD. The culture can be grown in a 24-deepwell plate, a falcon tube or a shake flask
- 12 Grow the strain at **37 °C** **Overnight**



Make sure to not incubate the overnight culture for longer than **16:00:00**. Using an overnight culture that has been incubating for longer than this, often results in non-reproducible results

Expression - Day 2 - 4/5

- 13 Prepare the desired volume of expression media in the desired vessel



Vessel and media volume has a massive influence on the final protein titers. In our experience, 24-deepwell plates in an incubator with a large shaking amplitude and CSTRs work well, while shake flasks give low protein yields. This is likely highly dependant on the equipment (e.g. incubators and the shape and baffle pattern of the shake flasks), and thus should be optimized for the individual labs

- 14 Inoculate the expression media to an OD₆₀₀ of 0.1

- 15 Incubate the expression culture at **20 °C** with 250 RPM shaking between **48:00:00** and **72:00:00**



The expression temperature and duration is dependant on the target protein, although the specified values seem to be a good starting point in our experience

Harvesting - Day 4/5

- 16 Harvest the culture by centrifuging at **6000 x g, 4°C, 00:05:00**



Depending on what the samples are for, the supernatant from the first centrifugation can additionally be centrifuged at **16000 x g, 4°C, 00:30:00** to clear out cell debris and other smaller contaminants

- 17 Keep the sample **On ice** when working with it and at **-20 °C** for storage