

JAN 18, 2024

OPEN BACCESS



DOI:

dx.doi.org/10.17504/protocols.io.x 54v9p4nmg3e/v1

Protocol Citation: Sarthak M. Singhal, thnasko 2024. Preparation of Recording ACSF (artificial cerebrospinal fluid) for slice electrophysiology. **protocols.io**

https://dx.doi.org/10.17504/protocols.io.x54v9p4nmg3e/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Preparation of Recording ACSF (artificial cerebrospinal fluid) for slice electrophysiology

Sarthak M. Singhal¹, thnasko¹

¹University of California San Diego



Sarthak M. Singhal

ABSTRACT

Hnasko lab- preparation of Recording ACSF Recording ACSF is used for incubating brain slices during cell-attached or whole-cell electrophysiological recordings of neurons.

MATERIALS

NaCl (Fisher Scientific, cat# BP358) KCl (Fisher Scientific, cat# BP366) NaH₂PO₄ (Fisher Scientific, cat# BP329)

D-glucose (Sigma Aldrich, cat# G8270) CaCl₂ (Fisher Scientific, cat# BP510)

MgSO₄ (Fisher Scientific, cat# M80)

NaHCO₃ (Fisher Scientific, cat# BP328)

HCI (Fisher Scientific, cat# A481-212)

Created: Jan 18, 2024

Last Modified: Jan 18, 2024

PROTOCOL integer ID: 93770

- **1** Prepare 10X stock solution.
 - 1.1 Composition of L 1 L of 10X Recording ACSF (in mM): 1250 NaCl, 25 KCl, 12 NaH₂PO₄, 125 d-glucose, 20 CaCl₂ and 20 MgSO₄.
 - 1.2 Composition of \square 1 L of 10X Recording ACSF (in g): 73.05 NaCl, 1.86 KCl, 1.44 NaH₂PO₄, 22.52 d-glucose, 2.94 CaCl₂ and 4.93 MgSO₄.
- 2 Dissolve all salts in \square 1 L of ddH₂O to make 10X stock solution of Recording ACSF.
- 3 Dilute 10X solution to 1X for use in experiments.
 - 3.1 Dilute \bot 100 mL of 10X solution to \sim \bot 800 mL with ddH₂O and add (in mM): 26 NaHCO₃; (in g): 2.18 NaHCO₃.
 - 3.2 Add ddH₂O to final volume of $\boxed{4}$ 1 L

Oct 18 2024

- 4 Bubble the solution with carbogen (95% O_2 5% CO_2) for ~10 min.
- After bubbling, measure and adjust pH of solution to \sim \bigcirc 7.3 with 12 N HCl (<1 ml).
- 6 Refrigerate until use.