



Nov 30, 2021

NOMe-seq of fixed cells

Florian Noack¹, Boyan Bonev¹

¹Helmholtz Pioneer Campus

1



protocol .



Boyan Bonev
Helmholtz Pioneer Campus

Protocol is based on the NOMe-seq protocol (Kelly et al.; 2012) with major modifications in the incubation times to work with fixed cells (Nordström et al.; 2019) and at the library preparation step.

Florian Noack, Boyan Bonev 2021. NOMe-seq of fixed cells. **protocols.io**
<https://protocols.io/view/nome-seq-of-fixed-cells-brdwm27e>



protocol ,

Jan 13, 2021

Nov 30, 2021

46230

Required Kits

Accel-NGS Methyl-Seq DNA Library Kit (Swift Bioscience, Cat. N.: 30024)
 Methyl-Seq Set A Indexing Kit (12 indices, 24 rxns) (Swift Bioscience, Cat. N.: 36024)
 EZ DNA Methylation-Gold Kit (Zymo Research, Cat. N.: D5005)
 Qubit dsDNA High Sensitivity Assay (ThermoFisher, Cat. N.: Q32851)

Buffers and Mastermixes

GpC methylation mix I (for control DNA)

- 5µl 10x GpC buffer (New England Biolabs, Cat. N.: M0227S)
- 1µl 32mM SAM (New England Biolabs, Cat. N.: M0227S)
- 5µl M.CviPI (New England Biolabs, Cat. N.: M0227S)
- 19µl nuclease free water

2% Formaldehyde solution

Dilute 1ml of 16% Formaldehyde solution (ThermoFisher, Cat. N.: 28908) with 7ml PBS.

2M Glycine solution

Mix 30.024g of Ultrapure Glycine (Invitrogen, Cat. N.: 15527013) with 200ml of PBS.

Lysis buffer (1ml)

10mM Tris-HCl, pH 8.0 (ThermoFisher, Cat. N.: 15567027)
10mM NaCl (ThermoFisher, Cat. N.: AM9760G)
0.2% Igepal-CA630 (Sigma-Aldrich, Cat. N.: I3021)
1x cOmplete™, EDTA-free Protease Inhibitor Cocktail (Roche, Cat. N.: 11873580001)
0.05% RNasin plus RNase inhibitor (Promega, Cat. N.: N261A)

For 1ml

- 10µl Tris-HCl, pH 8.0 (1M)
- 20µl Igepal-CA630 (10% stock solution, freshly prepared)
- 2µl NaCl (5M)
- 20µl cOmplete™, EDTA-free Protease Inhibitor Cocktail (50x stock)
- 0.5µl RNasin plus RNase inhibitor
- 947.5 µl nuclease free water

GpC wash buffer (250µl)

- 25µl GpC buffer (New England Biolabs, Cat. N.: M0227S)
- 50µl 5% Bovine Serum Albumin (ThermoFisher, Cat. N.: AM2618)
- 175µl nuclease free water

GpC reaction buffer (50µl)

- 7.5µl M.CviPI (New England Biolabs, Cat. N.: M0227S)
- 5µl 10x GpC buffer (New England Biolabs, Cat. N.: M0227S)
- 1µl 32mM SAM (New England Biolabs, Cat. N.: M0227S)
- 36.5µl nuclease free water

Library amplification mix (75µl)

- 2.5µl indexing primers of the Methyl-Seq Set A Indexing Kit (Swift Bioscience, Cat. N.: 36024)
- 25µl 5x EpiMark Hot Start Taq Reaction Buffer (New England Biolabs, Cat. N.: M0490)
- 2.5µl 10mM dNTPs (New England Biolabs, Cat. N.: N0447L)
- 0.65µl EpiMark Hot Start Taq (New England Biolabs, Cat. N.: M0490)
- 44.35µl nuclease free water

Prepare control DNA

1 NOTE: Control DNA has to be prepared only once and can be reused.

To prepare GpC methylated control DNA, mix 10µl of CpG methylated pUC19 DNA (Zymo Research, Cat. N.: D5017) with 10µl of unmethylated lambda DNA (Promega, Cat. N.: D1521).

2 Perform GpC methylation by mixing 20µl of the pUC19/lambda DNA mix with 30µl of the GpC methylation mix I and incubate for 1h @ 37 °C followed by 20 minutes at 65 °C.

3 Perform a 1x AmpureXP (Agencourt, Cat. N.: A63881) purification by adding 50 µl of room temperature AmpureXP beads to the samples followed by a incubation of 10 minutes at room temperature under slow rotation.

- 4 Pellet the beads on a magnetic-rack, remove the supernatant and wash the beads twice with freshly prepared 80% Ethanol.
- 5 After the last wash remove as much as possible of the 80% Ethanol and air-dry the beads for approximately 5 minutes at room temperature.
- 6 Resuspend the beads with 16µl nuclease-free water, incubate for 10 minutes at room temperature pellet the beads on a magnetic-rack and transfer the DNA containing supernatant into a Covaris microTUBE-15 AFA Beads Screw-Cap tubes (Covaris, Cat. N.: 520145).
- 7 Shear the DNA to a targeted size of 550bp using following parameter on a Covaris S220 sonicator:
 - Peak Incident Power (W) 18
 - Duty Factor 20%
 - Cycles per Burst 50
 - Treatment Time (s) 22
- 8 Perform another 1x AmpureXP purification as described above und elute the purified DNA in 30µl nuclease-free water.
- 9 Quantify the yield using the Qubit dsDNA High Sensitivity Assay (ThermoFisher, Cat. N.: Q32851). Low amounts of DNA expected. Control DNA can be stored at -20°C.

Fixation of cells

- 10 Resuspend dissociated cells in PBS to reach a maximal cell concentration of 2×10^6 cells/ml and add freshly prepared 2% Formaldehyde solution to reach a final concentration of 1%. Incubate for 10 minutes at room temperature with slow rotation.
- 11 Add 2.0M glycine solution (Invitrogen, Cat. N.: 15527013) to a final concentration of 0.2M to quench the reaction. Incubate at room temperature for 5 minutes with slow rotation.
- 12 Centrifuge for 5 minutes at 500xg at 4°C. Discard supernatant into an appropriate collection container.
- 13 Resuspend cells in 1ml of cold 1x PBS with 0.5% BSA (ThermoFisher, Cat. N.: AM2618) and count the cell number.

Fixed cells can be either FAC-sorted, directly used for NOME-seq or pelleted (2500xg for 5 minutes at 4°C) and then snap frozen in liquid nitrogen for storage at -80°C.

NOME-seq: Cell Lysis and GpC methylation

14 NOTE: this protocol was optimized for 1×10^4 cells.

Pellet fixed cells for 5 minutes at 2500xg (4°C) and carefully resuspend them with 100µl of ice-cold Lysis Buffer.

If cells were frozen, thaw the cell pellet first on ice.

15 Incubate on ice for exactly 10 minutes and centrifuge at 2500xg for 5 minutes at 4°C and remove the supernatant.

16 Wash the nuclei pellet by adding 250µl GpC wash buffer followed by a centrifugation at 2500xg for 5 minutes at 4°C and the removal of the supernatant.

17 Carefully resuspend nuclei with 50µl GpC reaction buffer.

18 Incubate for 3 hours at 37 °C with slight shaking (500 RPM) and add to the reaction every hour 0.5µl 32mM SAM (New England Biolabs, Cat. N.: M0227S) and 1µl of M.CviPI (New England Biolabs, Cat. N.: M0227S).

NOME-seq: Decrosslinking and shearing

19 Add 478µl freshly prepared Proteinase-K digestion mix and incubate at 55°C for 30 minutes.

20 Add 53µl of 5M NaCl (ThermoFisher, Cat. N.: AM9760G) and incubate at 68°C overnight (900rpm, 10sec, every 5 min.).

- 21 Cool tubes at room temperature, add 1:1 volume Phenol/Chloroform/Isoamylalcohol (Invitrogen, Cat. N. 15593031), vortex for at least 30 seconds and spin with maximal speed for 10 minutes at 4°C.
- 22 Take the upper aqueous phase into a 2ml tube and resuspend with 2x volumes of ice cold 100% ethanol, 1/10th of the volume 3M NaAc (Ambion, Cat. N.: AM9740) and 1µl Glycogen (ThermoFisher, Cat. N.: 10814010). Store the mix at -20°C at least for 1 hour or -80°C for 10-15 minutes.
- 23 Centrifuge at maximal speed for 20 minutes at 4°C. Carefully remove the supernatant by pipetting.
- 24 Wash with 800µl freshly prepared 70% cold ethanol followed by centrifuge at maximal speed for 5 minutes at 4°C.
- 25 Wash with 400µl freshly prepared 70% cold ethanol and transfer everything to a 0.5ml tube (make sure glycogen pellet is transferred). Centrifuge at maximal speed for 5 minutes at 4°C.
- 26 Remove as much as possible of the ethanol and air-dry the pellet. Dissolve pellet in 16µl of 10mM Tris buffer pH 7.5 (ThermoFisher, Cat. N.: 15567027) and incubate at 37°C for 15 minutes to fully dissolve DNA.
- 27 Transfer 16µl of the sample into a Covaris microTUBE-15 AFA Beads Screw-Cap tube (Covaris, Cat. N.: 520145) and sheare DNA to a targeted size of 550bp using following parameter on a Covaris S220 sonicator:
 - Peak Incident Power (W) 18
 - Duty Factor 20%
 - Cycles per Burst 50
 - Treatment Time (s) 22
- 28 Quantify the yield using the Qubit dsDNA High Sensitivity Assay.

NOMe-seq: Bisulfite Conversion and Lib. preparation

- 29 **NOTE: for the highest yield, the libraries should be prepared directly after the bisulfite conversion and not more than 100ng of sheared DNA should be used.**

Prior the bisulfite conversion add roughly 0.05% M.CviP-treated control DNA to the sample.
 Proceed with the bisulfite conversion using the EZ DNA Methylation-Gold Kit (Zymo Research)

according to the manual instructions until the elution step.

- 30 Elute the bisulfite converted DNA in 16µl low EDTA TE water provided with the Accel-NGS® Methyl-Seq DNA Library Kit (Swift Bioscience, Cat. N.: 30024).
- 31 Use 15µl of the bisulfite converted DNA for the library preparation using the Accel-NGS® Methyl-Seq DNA Library Kit according to the manual instruction until the final amplification step. Elute the DNA after the last purification step in 50µl low EDTA TE water.
- 32 Final amplification is performed with the with the EpiMark Hot Start Taq (New England Biolabs, Cat. N.: M0490S) in 5 separate reactions in order to achieve a higher complexity. To do so, pipette 10µl of the purified and adapter ligated DNA in 5 separate PCR tubes and add 15µl of the Library amplification mix containing sample specific Methyl-Seq Set A Indexing Primers (Swift Bioscience, Cat. N.: 36024).
- 33 Amplify library using the following PCR program: 95°C 30s; {95°C 15s, 61°C 30s, 68 °C 60s} x10-14; 68°C 5min; Hold at 10°C.
- 34 After the PCR is finished pool the different reactions of a sample together (~125µl total volume) and perform a 0.8x AmpureXP (Agencourt, Cat. N.: A63881) purification by adding 100µl of room temperature AmpureXP beads to the samples followed by a incubation of 10 minutes at room temperature under slow rotation.
- 35 Pellet the beads on a magnetic-rack, remove the supernatant and wash the beads twice with freshly prepared 80% Ethanol.
- 36 After the last wash remove as much as possible of the 80% Ethanol and air-dry the beads for approximately 5 minutes at room temperature.
- 37 Resuspend the beads with 15µl nuclease-free water, incubate for 10 minutes at room temperature, pellet the beads on a magnetic-rack and transfer the DNA containing supernatant into a new tube. Quantify the yield using the Qubit dsDNA High Sensitivity Assay (ThermoFisher, Cat. N.: Q32851)