



MAY 23, 2023

AAV viral DNA and whole RNA recovery for AAV pool experiments in rhesus macaque

Miguel Chuapoco¹¹California Institute of Technology

Miguel Chuapoco

California Institute of Technology

ABSTRACT

This protocol outlines procedures to extract viral DNA and whole RNA from rhesus macaque tissue that had been treated with AAV *in vivo*

OPEN  ACCESS**DOI:**

dx.doi.org/10.17504/protocols.io.3byl4jo68lo5/v1

Protocol Citation: Miguel Chuapoco 2023. AAV viral DNA and whole RNA recovery for AAV pool experiments in rhesus macaque.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.3byl4jo68lo5/v1>




License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 08, 2023**Last Modified:** May 23, 2023**PROTOCOL integer ID:**
81617**DNA/RNA extraction**

3h 8m

- 1 Add  100 mg tissue sample (brain or liver) and  1 mL TRIzol reagent  TRIzol Reagent Thermo Fisher Scientific Catalog #15596026 to bead homogenizer tubes. Use prefilled tubes with 1.5 mm Zirconium beads or 2.8 mm stainless steel beads.

Equipment

Prefilled 2.0ml tubes, Zirconium Beads, 1.5mm Triple-Pure - High Impact, 50pk	NAME
Homogenizer tubes (1.5 mm Zirconium beads)	TYPE
Benchmark Scientific	BRAND
D1032-15	SKU
https://www.benchmarkscientific.com/product/d1032-group/	LINK

Equipment

Prefilled 2.0ml tubes, Stainless Steel, 2.8mm Acid Washed, 50pk	NAME
Homogenizer tubes (2.8 stainless steel)	TYPE
Benchmark Scientific	BRAND
D1033-28	SKU
https://www.benchmarkscientific.com/product/d1032-group/	LINK

- 2 Homogenize tissue in using the following settings:

- Speed: 5.0 m/s
- Time: 30 seconds
- Pause: 1 minute
- Cycles: 2

Incubate for  00:05:00 .

5m

Equipment

BEADBUG 6, SIX POSITION HOMOGENIZER, 115V

NAME

Tissue homogenizer (6 position)

TYPE

Benchmark Scientific

BRAND

D1036

SKU

<https://www.benchmarkscientific.com/product/d1036/>

LINK

Equipment

BEADBLASTER 24 MICROTUBE HOMOGENIZER, 115V

NAME

Tissue Homogenizer (24 position)

TYPE

Benchmark Scientific

BRAND

D2400


SKU

<https://www.benchmarkscientific.com/product/d2400/>

LINK

Note

Samples can be stored at -20 °C for up to year in TRIzol.

- 3 Centrifuge the homogenizer tubes containing the TRIzol solution and homogenized tissue using the following parameters:  12000 x g, 4°C, 00:05:00. Transfer the supernatant to a new tube (microcentrifuge tube or similar).

5m

Equipment

Centrifuge 5425/5425 R - Microcentrifuge

NAME

Refrigerated centrifuge

TYPE

Eppendorf

BRAND

2231000909

SKU

<https://www.eppendorf.com/us-en/eShop-Products/Centrifugation/Microcentrifuges/Centrifuge-5425-5425R-p-PF-934144>

LINK

Equipment

DNA LoBind® Tubes

NAME

Microcentrifuge tubes

TYPE

Eppendorf








BRAND

022431021

SKU


<https://www.eppendorf.com/us-en/eShop-Products/Laboratory-Consumables/Tubes/DNA-LoBind-Tubes-p-PF-56252>




LINK

- 4 Add  200 µL chloroform to each tube for every  1 mL TRIzol used for lysis, vortex briefly, and incubate for  00:03:00 . 3m
- 5 Add 1 equivalent volume of isopropanol, 1/10 volume of sodium acetate, and co-precipitant (e.g.  500 µL isopropanol ,  50 µL sodium acetate ,  2-3 µL co-precipitant) and vortex briefly. Incubate for  00:10:00 . 10m





Sodium Acetate 3M, pH 5.2 Thermo Scientific Catalog
#R1181






 GlycoBlue™ Coprecipitant Thermo Scientific Catalog
#AM9516


- 6 Centrifuge  12000 x g, 4°C, 00:10:00 to pellet nucleic acids. Discard supernatant and wash pellet with  1 mL 75% ethanol. Centrifuge again  7500 x g, 4°C, 00:05:00 and discard supernatant.


15m

- 7 Air dry pellet and resuspend in  84 µL PCR clean water


 UltraPure Distilled Water Invitrogen - Thermo Fisher Catalog #10977-015



- 8 To isolate DNA, treat half of the sample with RNase. Remove RNA by digestion with  1.5 µL RNase cocktail and digest with  1.5 µL SmaI. Supplement reaction with  5 µL CutSmart. Incubate at  Room temperature for 2-3 hours and  37 °C overnight.

Purify with  Zymo DNA Clean & Concentrator - 5 Zymo Research Catalog
#D4014




 RNase Cocktail™ Enzyme Mix Thermo Fisher Catalog
#AM2286

 SmaI - 2,000 units New England Biolabs Catalog
#R0141S

 CutSmart Buffer - 5.0 ml New England Biolabs Catalog
#B7204S

- 8.1 To obtain cDNA, take  1 µg RNA from sample (measured by high sensitivity RNA Qubit) and treat with  DNase I, Amplification Grade Thermo Fisher Catalog
#18068015

15m

Combine  1 µg RNA with  1 µL 10X DNase I reaction buffer,  1 µL DNase I, and

 UltraPure DNase/RNase-Free Distilled Water Thermo Fisher Scientific Catalog
#10977023

to 10 µL.

Incubate reaction Room temperature 00:15:00

8.2

Inactivate DNase I by adding 1 μ L 25 mM EDTA. Heat 65 °C 00:15:00

15m

8.3

To convert DNase I treated RNA to cDNA, take 1-5 μ L sample and combine with 1 μ L oligo(dT), 1 μ L dNTP and fill to 10 μ L using UltraPure water.

SuperScript™ III Reverse Transcriptase Thermo Fisher Catalog #18080093

9

Incubate 65 °C 00:05:00 and place on ice.

5m

10

Prepare cDNA synthesis mix according to manufacturer's specifications. Add 10 μ L cDNA synthesis mix to each RNA primer mixture and mix by gently flicking the tubes.

11

Incubate as follows:

- 50 °C 00:50:00
- 25 °C 00:10:00
- 50 °C 00:50:00
- 85 °C 00:05:00

1h 55m

12

Store samples at -80 °C until ready to use.

PCR amplification

3h 8m

13

Use Zymo DNA purification of PCR product according to manufacturer's suggested protocol.




Zymo DNA Clean & Concentrator - 5 Zymo Research Catalog
#D4014

- 14 Dilute PCR product 1:100 and use as template for an additional round of PCR amplification around the variable region with primers containing Read1 and Read2 sequences by 10 cycles of

50s

98 °C for 00:00:10 , 60 °C for 00:00:30 , and 72 °C for 00:00:10

using  Q5 High-Fidelity DNA Polymerase - 500 units New England Biolabs Catalog
#M0491L

and the following primers:

- Forward: 5'- TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGcgatgtccagattacgcttgag -3'
- Reverse: 5'- GTCTCGTGGGCTCGGAGATGTGTATAAGAGACAGattttgtaatccagaggttgattatcg - 3'


- 15 Use Zymo DNA purification of PCR product according to manufacturer's suggested protocol.



Zymo DNA Clean & Concentrator - 5 Zymo Research Catalog
#D4014

- 16 Append Illumina flow cell adapters and unique indices by PCR amplification with

50s

 NEBNext Multiplex Oligos for Illumina (Dual Index Primers Set 1) - 96 rxns New England Biolabs Catalog #E7600S

by 10 cycles of 98 °C for 00:00:10 , 59 °C for 00:00:30 , and 72 °C for 00:00:10 using



Q5 High-Fidelity DNA Polymerase - 500 units New England Biolabs Catalog
#M0491L

- 17 Run PCR products on a freshly-prepared 2%



UltraPure™ Low Melting Point Agarose Thermo Fisher Scientific Catalog
#16520050

gel

and gel purify amplified 200 bp PCR product.