



Jul 23, 2020

# Bacterial Transformation Protocol

Olivier George<sup>1</sup>, Sierra Simpson<sup>1</sup><sup>1</sup>University of California, San Diego**1** Works for me This protocol is published without a DOI.

George Lab

Tech. support email: [olgeorge@ucsd.edu](mailto:olgeorge@ucsd.edu)

Sierra Simpson

University of California, San Diego

## ABSTRACT

Transformation of bacteria to amplify DNA for cloning, virus production, or other molecular biology techniques.

## PROTOCOL CITATION

Olivier George, Sierra Simpson 2020. Bacterial Transformation Protocol . [protocols.io](https://protocols.io)  
<https://protocols.io/view/bacterial-transformation-protocol-82yhyfw>

## KEYWORDS

Bacteria, molecular biology

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 05, 2019

## LAST MODIFIED

Jul 23, 2020

## PROTOCOL INTEGER ID

29496

## MATERIALS

NAME	CATALOG #	VENDOR
LB agar plates with the proper antibiotic (e.g. Kanamycin)		
MACH1 or DH5a or TOP10 or NEB stable or etc for DNA purification/miniprep		
Luria Broth Base (Miller's LB Broth Base)&trade;, powder	12795027	Thermo Fisher

## MATERIALS TEXT

### 14-956-9C Culture Tubes

## BEFORE STARTING

- Prepare Luria broth (LB) agar plates and allow to set.
- If pre-poured plates are being used, ensure the plates are warmed to 37° C.
- Depending on the antibiotic marker present in the plasmid DNA, incorporate appropriate antibiotic in the LB agar.
- Heat the water bath or heat block to 42° C.
- Warm the sterile LB medium to room temperature (or 20-25° C in water bath).

- 1 Transfer competent cells (DH5a) from -80 to wet ice for 5-10 min or until thawed.

- 2 Add 1ng to 50ng of DNA directly to cells. Incubate for 10 minutes on ice.
- 3 Heat shock cells for 45 seconds at 42C in a heat block. (Do not go over 45 seconds!) You can kill bacteria by keeping them at high temps for too long.
- 4 Transfer cells to ice and incubate for 5 minutes.
- 5 Add cells to sterile culture tubes with 5ml of LB + selection antibiotic. Ensure lid is not tight to ensure proper aeration of the cultures.
- 6 Incubate starter culture for 1h at 230rpm in a shaker at 37C.
- 7 Add 100ul of starter culture after incubation to warmed LB plate. Spread using roller beads.
- 8 Incubate plates o/n at 37C with the plate facing down to avoid desiccation of the agar.
- 9 In the morning, select colonies and culture in 3-5ml of Luria Broth (LB) overnight for Mini-prep in following morning.
- 10 Isolate the DNA from the culture.
- 11 Digest the plasmid with appropriate restriction enzyme and visualize with gel electrophoresis to determine if the plasmid insert is correct.
- 12 Send out the plasmid for sequencing before using to make virus, or functional studies.