

Aug 10, 2022

Ligation (Instructor Protocol)


Brian Teague¹¹University of Wisconsin - Stout

1 Works for me

Share

This protocol is published without a DOI.

Yeast ORFans CURE

 Brian Teague
University of Wisconsin - Stout

ABSTRACT

This is the instructor protocol for [the student Ligation protocol](#).

The abstract for the student protocol explains the basics. I pre-digest the L2-01 backbone for my students, but I generally do not gel-purify it. I find that the ligation is generally pretty robust, with three caveats:

- Because I don't gel-purify, we need to bias the ligation towards ligating the annealed oligos (not the GFP that was digested out). To do so, I dilute the backbone to 10 fmol/ul and instruct the students to dilute their annealed oligos to 200 fmol/ul (which is **200 nanomolar (nM)**).
- Ligase is EXPENSIVE, and this kind of ligation does not need much -- so I dilute 1:5 into small tubes. Thus, if a beginning student screws something up (throws the tube away, or contaminates it, etc), the expensive master reagent tube isn't lost.
- The ligase buffer contains ATP, which is heat-labile. (When a ligation fails, it's often because the buffer is bad, not the enzyme!) Thawing and re-freezing is not good for it. So, I aliquot the ligase buffer into single-use 5 ul aliquots in PCR tubes.

PROTOCOL CITATION

Brian Teague 2022. Ligation (Instructor Protocol). **protocols.io**
<https://protocols.io/view/ligation-instructor-protocol-cev9te96>



KEYWORDS

ligase, ligation, oligonucleotides

LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

By Madprime - Own work, CC BY-SA 3.0,
<https://commons.wikimedia.org/w/index.php?curid=2161789>

CREATED

Aug 05, 2022

LAST MODIFIED

Aug 10, 2022

PROTOCOL INTEGER ID

68257

MATERIALS TEXT

- E. coli freezer stock transformed with the L2-01 plasmid
- [LB agar Contributed by users](#) + 50 µg/ml
[Kanamycin Research Products International \(rpi\) Catalog #K22000-25.0](#)
- [LB Broth Contributed by users](#) + 50 µg/ml
[Kanamycin Research Products International \(rpi\) Catalog #K22000-25.0](#)
[Monarch Plasmid Miniprep](#)
- [Kit NEB Catalog #T1010](#) Step 3
- [Nuclease-free Water Contributed by users](#) In 2 steps
[Esp3I New England](#)
- [Biolabs Catalog # R0734S](#) Step 5
[CutSmart® Buffer New England](#)
- [Biolabs Catalog #B7204S](#) Step 5
- [Monarch DNA Elution Buffer - 25 ml New England](#)
[Biolabs Catalog #T1016L](#) Step 8
(optional)
- [10X NEB T4 DNA ligase buffer New England Biolabs](#) Step 10
[T4 DNA Ligase - 20,000 units New England](#)
- [Biolabs Catalog #M0202S](#) Step 11
[Diluent A - 5.0 ml New England](#)
- [Biolabs Catalog #B8001S](#) Step 11

SAFETY WARNINGS

Some components of the miniprep kit are hazardous; wear appropriate PPE.

TE, oligos, ligase & buffer, etc. are not hazardous. HOWEVER, we are shedding nucleases -- enzymes that degrade DNA -- all the time. Wear lab coats and gloves to keep your samples nuclease-free.

Grow & miniprep L2-01

- 1 At least 48 hours before the lab, strike out the L2-01 E. coli strain from a frozen stock on an LB+Kan plate.



Sanity check -- because the L2-01 plasmid has a GFP cassette, the colonies should be bright green!

- 2 At least 24 hours before the lab: pick a colony of L2-01 into **5 mL** LB+Kan liquid media. Grow in a round-bottomed test-tube overnight on a shaker, **200 rpm, 37°C, 16:00:00**



This is a pretty big plasmid, so the culture may not be as turbid as you're used to.



Sanity check -- the culture should be bright green!

- 3 Miniprep the culture using the [Monarch Plasmid Miniprep Kit](#) **NEB Catalog #T1010** or another suitable miniprep kit.
- 4 Analyze the eluate on a Nanodrop for DNA concentration and purity.



Typical elution: 50 ul @ 130 ng/ul (20 fmol/ul)

Linearize L2-01 3h 20m

- 5 We're going to linearize the ENTIRE miniprep. To a 200 ul PCR tube, add:
 - **50 µL** miniprep
 - **38 µL** [Nuclease-free Water Contributed by users](#)
 - **10 µL** [CutSmart® Buffer New England](#)
 - **2 µL** [Esp3I New England](#)
 - **2 µL** [Biolabs Catalog # R0734S](#)

Flick several times to mix well, then pulse down in a microcentrifuge.

- 6 Incubate δ **37 °C** for \odot **03:00:00** , then inactivate the reaction at δ **65 °C** for \odot **00:20:00** ^{3h 20m}

- 7 Optional: gel-purify the backbone. Run the entire miniprep on a preparative agarose gel, then use a kit such as the

[☒ Monarch DNA Gel Extraction Kit New England](#)

Biolabs Catalog #T1020S

to purify

the *larger* band.

- 8 Optional: dilute with

[☒ Monarch DNA Elution Buffer - 25 ml New England](#)

Biolabs Catalog #T1016L

to a

final concentration of 10 fmol/ μ l.

Aliquot / dilute reagents

- 9 Aliquot \square **1 mL** of [☒ Nuclease-free Water Contributed by users](#) per 4 students. Store at δ **Room temperature** .

- 10 Aliquot [☒ 10X NEB T4 DNA ligase buffer New England Biolabs](#) into \square **5 μ L** single-use aliquots in 200 μ l PCR strip tubes. Store at δ **-20 °C** .

I like to store these in a 200 μ l tip box.

- 11 [☒ T4 DNA Ligase - 20,000 units New England](#)

Dilute **Biolabs Catalog #M0202S**

1:5 in

[☒ Diluent A - 5.0 ml New England](#)

Biolabs Catalog #B8001S

. I make tubes containing

\square **2 μ L** of ligase and \square **8 μ L** of buffer, one per 4 students. Store at δ **-20 °C** .

12 Instructor Tips

- The first step of the student protocol requires another dilution. Again, I'm hands-off for this one. A common student question is "what is the concentration of my annealed oligos?" and I suggest that they might be able to figure it out by reviewing the annealing protocol.
- I usually instruct students to run this protocol and the URA3 PCR protocol at the same time: set up the ligation, then while the ligation is incubating, set up their PCR.

13 Common Student Errors

- Omitting the dilution step (particularly if they need to retry the protocol.)
- Small-volume pipetting.