



VERSION 2  
JAN 06, 2023

## OPEN ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v2](https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v2)

**Protocol Citation:** Andreas Sagen 2023. Lysogeny Broth (LB) medium. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1o262lr2/v2> Version created by [Andreas Sagen](#)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Jan 06, 2023

**Last Modified:** Jan 06, 2023


**PROTOCOL integer ID:**  
74863

**Keywords:** LB, bacteria, Lennox broth, E. coli, Escherichia coli

# Lysogeny Broth (LB) medium V.2

Andreas Sagen<sup>1</sup>

<sup>1</sup>University of Oslo



Andreas Sagen  
[University of Oslo, The National Institute of Occupational H...](#)

## ABSTRACT

**Lysogeny broth (LB)** is a nutritionally rich medium which is primarily used for the growth of bacteria<sup>[1]</sup>. The initialism is also commonly, albeit incorrectly, taken to mean **Luria broth**, **Lennox broth**, or **Luria-Bertani medium**. According to its creator Giuseppe Bertani, the abbreviation **LB** was actually intended to stand for **lysogeny broth**. The formula of the **LB medium** was published in 1951 in the first paper of Bertani on lysogeny<sup>[2]</sup>. There are several common formulations of **LB**. Although they are different, they generally share a somewhat similar composition of ingredients used to promote growth, including the following: Peptides and casein peptones, Vitamins (including B vitamins), Trace elements (e.g. nitrogen, sulfur, magnesium) and Minerals.

Sodium ions for transport and osmotic balance are provided by sodium chloride (NaCl). Tryptone is used to provide *essential amino acids* such as peptides and peptones to the growing bacteria, while the yeast extract is used to provide a plethora of organic compounds helpful for bacterial growth. These compounds include vitamins and certain trace elements.

## GUIDELINES

Follow step by step, unless stated otherwise. Equipment needed should be standard to a microbiology lab.

## MATERIALS

Analytical scale, autoclave, bottle, weight vessel, LAF bench

## SAFETY WARNINGS





⚠ When removing autoclaved components, be sure to take care as this can be very hot. If using antibiotics, use sufficient PPE to protect yourself, as some can be toxic to humans.

## BEFORE START INSTRUCTIONS

Prepare glassware by cleaning it, and ensure that scale is sufficiently calibrated

- 1 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

## 500 mL LB-Lennox (broth) medium


- 1.1 Fill the bottle with  400 mL double-distilled water
- 1.2 Measure  5000 mg Tryptone,  2500 mg Yeast extract and  2500 mg Sodium chloride

Powdered compounds:


 Tryptone **Millipore Catalog #T9410**


 Yeast Extract **Sigma-aldrich Catalog #Y0875**


 Sodium chloride **Sigma-aldrich Catalog #S9625**

- 1.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for  00:05:00

5m

- 1.4 Adjust pH while mixing to  6.7 using concentrated sodium hydroxide

- 1.5 Add distilled water to a total of  500 mL

- 1.6 Autoclave liquid at  121 °C for  00:15:00


15m





### Note

Cool to 50°C and supplement with antibiotics as appropriate

- 2 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

## 500 mL LB-Lennox (agar) medium

2.1 Fill the bottle with  400 mL double-distilled water

2.2 Measure  5000 mg Tryptone,  2500 mg Yeast extract,  2500 mg Sodium chloride and  7500 mg agar


Powdered compounds:

 Tryptone **Millipore Catalog #T9410**

 Yeast Extract **Sigma-aldrich Catalog #Y0875**

 Sodium chloride **Sigma-aldrich Catalog #S9625**

 Agar **Sigma-aldrich Catalog #A1296**

2.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for  00:05:00

5m

2.4 Adjust pH while mixing to  6.7 using concentrated sodium hydroxide

2.5 Add distilled water to a total of 500 mL

2.6 Autoclave liquid at  121 °C for  00:15:00

15m

#### Note

Cool to 50°C and supplement with antibiotics as appropriate

Agar can be stored, then reheated to 50°C to be poured