

Sep 26, 2024

pGEM-T Cloning

DOI

dx.doi.org/10.17504/protocols.io.yxmvme389g3p/v1

Carolina Lopez¹

¹Washington University



Cecilia Escudero

wustl

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.yxmvme389g3p/v1**

Protocol Citation: Carolina Lopez 2024. pGEM-T Cloning. protocols.io **<https://dx.doi.org/10.17504/protocols.io.yxmvme389g3p/v1>**

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 28, 2024

Last Modified: September 26, 2024

Protocol Integer ID: 100760

Abstract

Protocol for cloning in pGEM-T vector



Materials

Reagents:

- pGEM-T Easy Kit (Promega, A1360)
- Competent E. coli cells (NEB, C3040H)
- LB plates with necessary antibiotic

Before start

PRINCIPLES BEHIND THE PROCEDURE MUST BE UNDERSTOOD. PLEASE CONSULT WITH EXPERIENCED LAB MEMBER THE FIRST TIME YOU USE THIS PROCEDURE. UPDATE AS A GENERAL PROCEDURE AS NECESSARY BUT DO NOT MODIFY WITH SPECIFICS TO YOUR PROJECT, INSTEAD DOWNLOAD AND PASTE A MODIFIED COPY IN YOUR NOTEBOOK. Updated by MH, April 2023



Insert Generation

1h 2m

- 1 1. Perform PCR to generate bands of interest in total volume of 50 μL .

1h 2m

	Volume (μL)
NF H ₂ O	37
Fwd Primer (10 μM)	1
Rev Primer (10 μM)	1
dNTPs	1
Plasmid (20ng)	1
Buffer	5
MgCl ₂	3
Taq Polymerase	1

2. Make 1% low melt-agarose gel.

- Mix 1 g of Agar with 100 mL of TAE Buffer.
- Microwave to boil agarose and let cool until you can touch bottle, but gel is not solid.
- Add 1.5 μL of EtBr to agarose and pour into DNA gel mold with 10 well comb.
- Let gel solidify.

3. Add 10 μL of loading buffer to PCR reaction then Load 60 μL of reaction into well of gel

4. Run gel for 00:45:00 at 120V.

5. Visualize band with UV light and cut out section of gel with band with new razor blade and place in 1.5 mL tube.

6. Purify Band from gel with Promega Wizard SV Gel and PCR Purification Kit (A9282)

- <https://www.promega.com/products/nucleic-acid-extraction/clean-up-and-concentration/wizard-sv-gel-and-pcr-clean-up-system/?catNum=A9281>

- Value
- Weight DNA gel fragment and add 10 μL of Membrane Binding Solution per 10 mg of gel slice.

- Value
- Incubate mixture at 65 $^{\circ}\text{C}$ for 00:10:00 or until gel is completely melted.


- Value
- Value
- Add melted gel mixture to SV minicolumn in Collection Tube and incubate at room temperature for 00:01:00.

- Value
- Value
- Centrifuge at max speed for 00:01:00. Discard flowthrough and reinsert column into tube.

- Value
- Value
- Value
- Add 700 μL of Membrane Wash Solution. Centrifuge at max speed for




- Value
- Value
- Value
- Value
- 00:01:00. Discard flowthrough and reinsert column into tube.




g) Add  500 μL of membrane wash solution. Centrifuge at max speed for

 00:01:00 . Discard flowthrough and reinsert column into tube.

h) Spin empty column for  00:01:00 at max speed to remove excess ethanol.


i) Transfer column to labelled  1.5 mL tube and add  35 μL of NF H₂O. Incubate for  00:01:00 at room temperature.

j) Centrifuge at max speed for  00:01:00 .

k) Keep eluate and store at  -20 °C .

pGEM Ligation Protocol

2h 24m

2 1. Set up Ligation Reaction in  0.2 mL PCR tube


2h 24m

Component	Volume (uL)
2x Rapid Ligation Buffer	5
pGEM-T Vector (50ng)	1
PCR Product	X*
H ₂ O	3-X
T4 DNA Ligase	1



a) *Volume of PCR product should be 2 insert:1 pGEM vector molar ratio

b) For example:


A	B	C	D	E	F	G
Component	Length of DNA (bp)	Molar ratio	ng of DNA	Volume of 50ng/ul solution		
pGEM-T Vector	3015	1	50	1 ul		
PCR Fragment	500	2	16.58	0.33 ul		
H ₂ O				2.67 ul		


2. Incubate for  01:00:00 at room temperature

3. Transform Product into E. coli






a) Add  2 μL of product to  50 μL of TOP10 cells and mix up and down slowly with a pipette.

b) Incubate for  00:20:00 on ice.

c) Heat shock bacteria in  42 °C waterbath for  00:01:00 .

d) Incubate on Ice for  00:03:00 .



- e) Add  100 μ L of SOC media and shake in warm room for  01:00:00 .
- f) Plate bacteria onto LB-Antibiotic/Xgal Plate and incubate overnight at  37 °C
- i. Add  40 μ L of Xgal to plate and spread with plate spreader. Let dry
- ii. Then plate  100 μ L of bacteria onto plate.
- iii. Spread cells with plate spreader to get individual colonies
- g) Pick white colonies for miniprep growth and sequencing.

Protocol references

pGEM-T Cloning:

<https://www.promega.com/products/pcr/pcr-cloning/pgem-t-easy-vector-systems/?catNum=A1360>