

9 ▼

Mar 22, 2022

# Total RNA and DNA from Microalgae (24 samples per day) V.9

Yingyu Hu<sup>1</sup>, Zoe V Finkel<sup>1</sup>

<sup>1</sup>Dalhousie University

1



protocol .

**Marine Microbial Macroecology Lab**  
Tech. support email: [ruby.hu@dal.ca](mailto:ruby.hu@dal.ca)



Yingyu Hu  
Dalhousie University

Here we describe a protocol for extracting and quantifying bulk RNA and DNA from microalgae, which is adapted from Berdalet E. et al. (2005).

RNA and DNA are extracted from microalgae samples and then quantified by fluorochrome SYBR Green II.

The level of sensitivity of this method was set at ca. 40 ~300 ng RNA or 10 ~ 100 ng DNA (ml assay)<sup>-1</sup>.

Berdalet E, Roldán C, Olivar MP, Lysnes K. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part A. Optimisation of the assay. *Scientia Marina*.  
<https://doi.org/10.3989/scimar.2005.69n11>

Berdalet E, Roldán C, Olivar MP. Quantifying RNA and DNA in planktonic organisms with SYBR Green II and nucleases. Part B. Quantification in natural samples. *Scientia Marina*.  
<https://doi.org/10.3989/scimar.2005.69n117>

Yingyu Hu, Zoe V Finkel 2022. Total RNA and DNA from Microalgae (24 samples per day) . **protocols.io**  
<https://protocols.io/view/total-rna-and-dna-from-microalgae-24-samples-per-d-b6j2rcqe>  
Yingyu Hu



Emphasize the incubation time

RNA, DNA, SYBR Green II, DNase, RNase, microalgae, fluorescence

protocol ,

Mar 22, 2022

Mar 22, 2022

59738

Estimation of RNA/DNA in the collected microalgae samples:

Under replete condition, RNA and DNA is about 5.7% and 1% in total dry mass, while Chl-a is about 1.1% in total dry mass. Therefore,  $\text{RNA}_{\text{ug/L}} = \text{Chl-a}_{\text{ug/L}} \times (5.7/1.1)$ ,  $\text{DNA}_{\text{ug/L}} = \text{Chl-a}_{\text{ug/L}} \times (1/1.1)$ .

Common dilution from sample collected on the filter to assay is 1/40.

#### STEP MATERIALS

[☒ Nuclease decontamination](#)

[solution IDT Catalog #11-05-03-01](#) Step 2

[☒ Ribonuclease A from bovine pancreas Sigma](#)

[Aldrich Catalog #R6513-50MG](#) Step 6.1

[☒ DEOXYRIBONUCLEASE1 RNase and Protease](#)

[Free Bioshop Catalog # DRB002.10](#) Step 7.1

[☒ Magnesium chloride solution Sigma](#)

[Aldrich Catalog #63069-100ML](#) Step 40.1

[☒ Calcium chloride solution Sigma](#)

[Aldrich Catalog #21115-100ML](#) Step 42.1

[☒ SYBR™ Green II RNA Gel Stain, 10,000X concentrate in DMSO Thermo](#)

[Fisher Catalog #S7564](#) Step 49

[☒ Tris\(hydroxymethyl\)aminomethane hydrochloride 1M pH 8.0 RNase free Fisher](#)

[Scientific Catalog #AAJ60080AK](#) Step 3.1

[☒ Deoxyribonucleic acid from calf thymus Sigma](#)

[Aldrich Catalog #D4522-1MG](#) Step 5.1

[☒ N-Lauroylsarcosine sodium salt solution \(20% RNase/DNase free\) Sigma](#)

[Aldrich Catalog #L744-50mL](#) Step 10.1

[☒ EDTA buffer \(0.5M DNase/RNase](#)

[free\) Bioshop Catalog #EDT333.100](#) Step 10.2

[☒ UltraPure™ DNase/RNase-Free Distilled](#)

[Water ThermoFisher Catalog #10977023](#) Step 3.2

[☒ E. coli Total RNA Thermo Fisher](#)

[Scientific Catalog #AM7940](#) Step 4.1

No data is available addressing the mutagenicity or toxicity of SYBR® Green II Nucleic Acid Gel Stain. Because this reagent binds to nucleic acids, it should be treated as a potential mutagen and used with appropriate care. The DMSO stock solution should be handled with particular caution as DMSO is known to facilitate the entry of organic molecules into tissues. We strongly recommend using double gloves when handling the DMSO stock solution. As with all nucleic acid stains, solutions of SYBR Green II Nucleic Acid Gel Stain should be poured through activated charcoal before disposal or collected in waste container to be treated later. The charcoal must then be incinerated to destroy the dye.

#### Day 1: Freeze-dry samples

- 1 Freeze dry samples and blank filters. Freeze at **-80 °C** until processed.

1. Freeze-drying should be as short as possible to reduce sample degradation.
2. The exact duration of freeze-drying depends on size of filter, quantity of sample and the size of container.

FreeZone® 2.5 L Benchtop Freeze Dryers  
Labconco® 700202000

#### Day 1: Prepare primary solutions

- 2 Turn on UV light in biosafety cabinet for ⌚00:15:00 and clean working surface with decontamination solution.

☒Nuclease decontamination

solution IDT Catalog #11-05-03-01

- 3 Prepare Tris buffer [M]5 mM pH8.0

- 3.1 Pour [M]1 M pH8.0 Tris into an RNase free 15 mL Falcon tube

☒Tris(hydroxymethyl)aminomethane hydrochloride 1M pH 8.0 RNase free Fisher Scientific Catalog #AAJ60080AK

Falcon® Centrifuge Tubes  
Polypropylene, Sterile, 15 mL  
Corning® 352096

- 3.2 Directly add [V]2.5 mL [M]1 M pH8.0 Tris into 500 mL RNase free water in its original package.

☒UltraPure™ DNase/RNase-Free Distilled

Water Thermofisher Catalog #10977023

BT Barrier Pipet Tips  
Pre-Sterile  
Neptune® BT1250, BT100, BT10

- 4 RNA primary standard solution ([M]200 ug/ml )

- 4.1 In the original package, the E. Coli Total RNA is of 1 mg/mL, in which total RNA is 200 ug.

☒E. coli Total RNA Thermo Fisher

Scientific Catalog #AM7940

[https://assets.thermofisher.com/TFS-Assets/LSG/manuals/sn\\_7940.pdf](https://assets.thermofisher.com/TFS-Assets/LSG/manuals/sn_7940.pdf)

<https://doi.org/10.1016/j.jm.2017.05.001>

- Cap the package and vortex for a thorough mix.

Finnpipette Stepper Pipette  
Thermo Scientific™ 4540000

Finntip stepper pipette tips  
500 ul (sterile)

Thermo Scientific      Thermo Scientific™  
9404173

Microcentrifuge Tubes  
1.7 mL/0.6 mL  
Axygen Scientific MCT-175-C/MCT-060-L-C

5 DNA primary standard solution ( $\approx$  **[M]500 ug/ml**)

- (**[M]5 mM** , **pH8.0** ).

⊗ Deoxyribonucleic acid from calf thymus **Sigma**

Aldrich Catalog #D4522-1MG

[https://www.sigmaaldrich.com/content/dam/sigmaaldrich/docs/Sigma/Product\\_Information\\_Sheet/d4522pis.pdf](https://www.sigmaaldrich.com/content/dam/sigmaaldrich/docs/Sigma/Product_Information_Sheet/d4522pis.pdf)

- Cap the package. Do not vortex or sonicate.

5.4 Aliquot 10 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at  $-80^{\circ}\text{C}$ .

Finntip stepper pipette tips  
500 uL (sterile)  
Thermo Scientific Thermo Scientific™  
9404173

5.5 Dilute 5 uL primary DNA standard solution with 95 uL Tris buffer ( $10\text{ mM}$ ,  $\text{pH } 8.0$ ) in a microtube (600 uL).  
Measure DNA concentration by using  $\mu$ drop plate (sample volume: 4 uL)  
Use Tris buffer ( $10\text{ mM}$ ,  $\text{pH } 8.0$ ) as blank.

$\mu$ Drop™ Plates  
Thermo Scientific N12391

Varioskan LUX Multimode Microplate  
Reader  
Thermo Fisher VL0L00D0

The dilution is to avoid saturated observation at 260 nm.

5.6  $\text{DNA concentration } (\mu\text{g/ml}) = (\text{Abs}_{260} - \text{Abs}_{260}(\text{blank})) \times 50 \mu\text{g/ml} \times (10\text{mm}/0.5\text{ mm}) \times \text{DF}$   
Where,  $\text{DF}=20$ .

If the measured DNA concentration is not close to  $500 \mu\text{g/ml}$ , check reverse pipetting technique.

6 RNase primary stock solution ( $10\text{ mg/ml}$ )

6.1 Uncap the original package of Ribonuclease A from bovin pancreas and add  $5\text{ mL}$  Tris buffer ( $10\text{ mM}$ ,  $\text{pH } 8.0$ ).  
Cap the package and vortex for a thorough mix.

☒ Ribonuclease A from bovine pancreas Sigma

Aldrich Catalog #R6513-50MG

- 6.2 Aliquot 30 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at  $-20^{\circ}\text{C}$ .

Finntip stepper pipette tips  
500 ul (sterile)

Thermo Scientific Thermo Scientific™  
9404173

Finntip™ Stepper Pipette Tips  
500 ul (Sterile)

Thermo Scientific 21-377-149

- 7 DNase primary stock solution ( $1\text{ mg/ml} = 10,000\text{ U/mL}$ )

- 7.1 Uncap the original package of Deoxyribonuclease1 and add  $1\text{ mL}$  Tris buffer ( $5\text{ mM}$ ,  $\text{pH}8.0$ )

Cap the package and vortex for a thorough mix.

☒ DEOXYRIBONUCLEASE1 RNase and Protease

Free Bioshop Catalog # DRB002.10

- 7.2 Aliquot 100 uL by stepper with sterile tip to 600 uL RNase free microtubes. Keep frozen at  $-20^{\circ}\text{C}$

Finntip stepper pipette tips  
1250 ul (sterile)

Thermo Scientific Thermo Scientific™  
9404183

## Day 2: Exact RNA and DNA

- 8 Turn on UV light in biosafety cabinet for  $00:15:00$  and clean working surface with decontamination solution.
- 9 Prepare falcon tubes and tube rack in biosafety cabinet

A	B
Volume of tube (mL)	Contents in the tube
5	0.5 M EDTA
5	20% sarcosine
15 or 50	5 mM Tris
15 or 50	1% STEB

Falcon® Centrifuge Tubes  
Polypropylene, Sterile, 15 mL  
Corning® 352096

Falcon® Centrifuge Tubes  
Polypropylene, Sterile, 50 mL  
Corning® 352070

## 10 Prepare STEB (1% )

Use the following formula to determine the total volume of 1% STEB required:  
(# samples + # blank filters) X (500 ul) + (500 ul) = total volume of 1% STEB required

10.1 Pour sarcosine (20% ) into an RNase free 5 mL falcon tube.

[N-Lauroylsarcosine sodium salt solution \(20% RNase/DNase free\)](#) **Sigma**  
**Aldrich Catalog #L744-50mL**

10.2 Pour EDTA (0.5 M ) into an RNase free 5 mL falcon tube.

[EDTA buffer \(0.5M DNase/RNase free\)](#) **Bioshop Catalog #EDT333.100**

10.3 Pour Tris buffer (5 mM , pH 8.0 ) into an RNase free 15 or 50 mL falcon tube.

For large volume of Tris buffer usage, use sterile serological pipet

10.4 Mix 500 µL sarcosine (20% ), 10 µL EDTA (0.5 M ) and 9 mL + 490 µL Tris buffer (5 mM , pH 8.0 ) to obtain STEB (1% ).

11 Prepare ice bath

12 Remove freeze-dried samples from -80°C freezer and place them [On ice](#).

13 Add [500 µL](#) Tris buffer ( [5 mM](#) , [pH 8.0](#) ) and [500 µL](#) STEB ( [1 %](#) ) to the bead tube. Place tubes [On ice](#).

LYSING TUBES  
MATRIX D 2 mL/15 mL  
MP BIOMEDICALS 116913500/116933050

14 Rinse forceps by [70 % volume](#) ethanol and air dry.

Filter forceps  
blunt end, stainless steel  
Millipore XX6200006P

15 Transfer sample/blank filter into the bead tube by using clean forceps.

16 Vortex immediately then put back [On ice](#).

20s

VWR ANALOG VORTEX MIXER  
VWR 10153-838  
With tube insert

17 Disrupt samples on the bead mill at 6.5 m/s.






30s

Fastprep-24 5G™ Sample Preparation  
Instrument  
MP Biomedicals 116005500


18 Keep tubes [On ice](#). Check the label on each tube, restore the label if it fades.

1m 30s



- 19 Disrupt samples on the bead mill at 6.5 m/s. 30s
- 20 Keep tubes  **On ice** . Check the label on each tube, restore the label if it fades. 1m 30s
- 21 Disrupt samples on the bead mill at 6.5 m/s 30s
- 22 Keep tubes  **On ice** . Check the label on each tube, restore the label if it fades. 1m 30s
- 23 Disrupt samples on the bead mill at 6.5 m/s. 30s
- 24 Continuously shake homogenate in a multi-head vortex at the highest speed for  **01:00:00**  **Room temperature** 1h
- Vortex mixer should be able to remain stable on the bench under this vortex speed.
- 25 In the biosafety cabinet, transfer 150 uL of homogenate into RNase free 600 uL tube.
- 26 Freeze both aliquote (150 uL) and the rest of homogenate (in bead mill tube) at  **-80 °C** until analyzed.

#### Day 3: Run the assay

- 27 Prepare ice bath.
- 28 Turn on UV light in biosafety cabinet for  **00:15:00** and clean working surface with decontamination solution.
- 29 Prepare falcon tubes, microtubes and tube racks in biosafety cabinet  
\* Maximum number of samples (including blanks) per assay is 20.

A	B	C
Number of tubes	Type of tubes	Contents
3	5 mL falcon tubes	1 M MgCl <sub>2</sub> , 1 M CaCl <sub>2</sub> , Sybr Green II working solution (SG-II WS)
1	50 mL falcon tube	5 mM Tris buffer
4	15 mL falcon tubes	0.05% STEB, Working solution A (WS-A), Working solution B (WS-B), Working solution (WS-C)
5	1.7 mL RNase free tubes	RNase working solution, Secondary RNA standard solution, Secondary DNA standard solution, 900mM MgCl <sub>2</sub> , 900 mM CaCl <sub>2</sub>
33	1.7 mL RNase free tubes	RNA standard solutions for RNA standard curves, DNA standard solutions for DNA standard curves
N= total number of samples and blanks	1.7 mL RNase free tubes	Samples and blanks
3XN	1.7 mL RNase free tubes	Diluted samples and blanks
5	Microtube racks	Tubes of 1.7 mL
1	Tube racks	Falcon tubes

Screw-Cap Centrifuge Tube  
5 mL  
VWR 10002-738

### Day 3: Run the assay

## 30 Organize and label the tubes as shown below

### Set 1:

In microtube rack, label 1.7 mL tubes for samples and blanks to be further diluted.

1	2	3	4	5	6	7	8	9	...	Blank
---	---	---	---	---	---	---	---	---	-----	-------

### Set 2:

In microtube rack, label 1.7 mL tubes for RNA (marked in pink) and DNA (marked in blue) standard solutions to be used as standard curves.

Tubes A is for standard solutions treated with working solution A (WS-A)

Tubes B is for standard solutions treated with working solution B (WS-B)

Tubes C is for standard solutions treated with working solution C (WS-C)

Tubes A	R1A	R2A	R3A	R4A	R5A	R6A	D1A	D2A	D3A	D4A	D5A
Tubes B	R1B	R2B	R3B	R4B	R5B	R6B	D1B	D2B	D3B	D4B	D5B
Tubes C	R1C	R2C	R3C	R4C	R5C	R6C	D1C	D2C	D3C	D4C	D5C

### Set 3:

In microtube rack, label 1.7 mL tubes for diluted samples and blanks.

Tubes A is for diluted samples and blanks treated with working solution A (WS-A)

Tubes B is for diluted samples and blanks treated with working solution B (WS-B)

Tubes C is for diluted samples and blanks treated with working solution C (WS-C)

Tubes A	1A	2A	3A	4A	5A	6A	7A	8A	9A	...	BlankA
Tubes B	1B	2B	3B	4B	5B	6B	7B	8B	9B	...	BlankB
Tubes C	1C	2C	3C	4C	5C	6C	7C	8C	9C	...	BlankC

- 31 Label tubes for reagents as following.  
Follow the sheet, add Tris buffer (1M 5 mM , pH 8.0 ) to the reagent tubes:

A	B
Reagent type	5 mM Tris (uL)
0.05% STEB	9X1000 + 500
RNA	990+495
DNA	998
RNase	380
900 mM MgCl <sub>2</sub>	40
900 mM CaCl <sub>2</sub>	40
WS-A	5X1000 + 640
WS-B	5X1000 +640
WS-C	5X1000 + 880
SG-II WS	2X1000 + 500

- 32 Add 900 µL Tris buffer (1M 5 mM , pH 8.0 ) to each tube in **Set 1**

Depending on the dilution of extracted sample

- 33 Follow the sheet, add Tris buffer (1M 5 mM , pH 8.0 ) to each tube in **Set 2**. The unit of volume is uL.

Tubes A	650	640	625	600	550	500	640	630	610	580	550
Tubes B	650	640	625	600	550	500	640	630	610	580	550
Tubes C	600	590	575	550	500	450	590	580	560	530	500

- 34 Follow the sheet, add Tris buffer (1M 5 mM , pH 8.0 ) to each tube in **Set 3**. The unit of volume is uL.

Tubes A	650	650	650	650	650	650	650	650	650	...	650
Tubes B	650	650	650	650	650	650	650	650	650	...	650
Tubes C	600	600	600	600	600	600	600	600	600	...	600

- 35 Prepare STEB (1M 0.05 % )

Add 500 µL STEB (1M 1 % ) to 0.05% STEB tube, and vortex.

- 36 Add 250 µL STEB (1M 0.05 % ) to each tube in **Set 2** by reverse pipetting.

- 37 Place RNase and DNase primary stock solutions, RNA and DNA primary standard solutions and samples **On ice** .

- 38 Turn on refrigerated centrifuge and set the temperature to **4 °C** .

CENTRIFUGE 5430 R  
Eppendorf MP2231000510

39 Turn on shaker/incubator and set temperature to  $37^{\circ}\text{C}$ .

SHAKING INCUBATOR  
71L  
Corning® LSE™ 6753

40 Prepare  $900\text{ mM}$   $\text{MgCl}_2$

40.1 Pour  $1\text{ M}$   $\text{MgCl}_2$  solution into 5 mL RNase free Falcon tube

[☒ Magnesium chloride solution Sigma](#)

[Aldrich Catalog #63069-100ML](#)

40.2 Transfer  $360\text{ }\mu\text{L}$   $1\text{ M}$   $\text{MgCl}_2$  solution into  $900\text{ mM}$   $\text{MgCl}_2$  tube

41 Add  $120\text{ }\mu\text{L}$   $900\text{ mM}$   $\text{MgCl}_2$  to WS-A and WS-B

42 Prepare  $900\text{ mM}$   $\text{CaCl}_2$

42.1 Pour  $1\text{ M}$   $\text{CaCl}_2$  solution into 5 mL RNase free Falcon tube

[☒ Calcium chloride solution Sigma](#)

[Aldrich Catalog #21115-100ML](#)

42.2 Transfer  $360\text{ }\mu\text{L}$   $1\text{ M}$   $\text{CaCl}_2$  solution into  $900\text{ mM}$   $\text{CaCl}_2$  tube

43 Add  $120\text{ }\mu\text{L}$   $900\text{ mM}$   $\text{CaCl}_2$  to WS-A and WS-B

44

Lunch break!

45 Prepare RNase working solution [M]0.5 mg/ml  
Add ▢20 µL RNase primary stock solution ([M]10 mg/ml) to RNase tube

46 Add ▢120 µL [M]0.5 mg/ml RNase to WS-B and WS-C.  
Keep WS-B and WS-C ⚡ On ice .

47 Add ▢120 µL DNase primary stock solution ([M]5 mg/ml) to WS-A.  
Keep WS-A ⚡ On ice .

Two microtubes of aliquot (70 ul/tube)

48 Centrifuge extracted samples 🌀10000 x g, 4°C, 00:04:00

49 Take one tube of  
▢SYBR™ Green II RNA Gel Stain, 10,000X concentrate in DMSO Thermo  
Fisher Catalog #S7564  
out of the freezer, keep it at ⚡ Room temperature .  
If open a new package, wrap 1.7 mL microtube with foil and aliquot 1000 ul to each tube, store at ⚡ -20 °C .

50 Prepare RNA secondary standard solution [M]2 ug/ml  
Add ▢15 µL RNA primary standard solution to RNA standard tube and vortex.  
Keep ⚡ On ice .

51 Prepare DNA secondary standard solution [M]1 ug/ml  
Add ▢2 µL DNA primary standard solution to DNA standard tube and vortex.  
Keep ⚡ On ice .


52 Load ▢50 µL WS-A to Tubes A in Set 2 and Set 3.


From Step 51 to 54: Reverse pipetting  
Decontaminate pipet between different WS.

53 Load ▢50 µL WS-A to Tubes C in Set 2 and Set 3.

54 Load ▢50 µL WS-B to Tubes B in Set 2 and Set 3.

55 Load ▢50 µL WS-C to Tubes C in Set 2 and Set 3.




56 Add  **100 µL** centrifuged samples to its corresponding tubes in **Set 1**.  
Vortex each tube.

57 From Set 1, transfer  **250 µL** of diluted samples to each corresponding tubes in **Set 3**.

In order to avoid cross contamination from RNase or DNase, use one tip for each dispensing.  
Pipette solution in the tube up and down for mixing.

58 Follow the sheet:  
Add RNA secondary standard to tubes (marked in pink) in Set 2.  
Add DNA secondary standard to tubes (marked in blue) in Set 2.  
The unit of volume is µL.

Tubes A	0	10	25	50	100	150	10	20	40	70	100
Tubes B	0	10	25	50	100	150	10	20	40	70	100
Tubes C	0	10	25	50	100	150	10	20	40	70	100

59 Vortex each tube for  **00:00:02** and place all tubes into the shaker/incubator at  **37 °C**, continuously shaking at <sup>20m</sup> 200 RPM for  **00:20:00**.



Incubation time is critical. Temperature might be disturbed by door open/close. Don't start the timer until temperature returns to 37°C.


#### Read fluorescence

60 Prepare SYBR Green II working solution (SG-II WS)

60.1 Each 96-well microplate requires 1 mL of SG-II WS.

For 24 samples:  
2.5 mL Tris  
17.5 µL SG-II WS

60.2 Centrifuge one tube of SG-II concentrate at  **Room temperature**  **13000 rpm, 00:05:00** <sup>5m</sup> to deposit DMSO.

60.3 Wrap SG-II WS tube with foil, transfer  **17.5 µL** *supernatant* of SYBR Green II 10,000X concentrate to SG-II WS tube in biosafety cabinet.

Any step involving SYBR Green II should be operated in dark room or at least dim light.


60.4 Adhere black film on the top of a microplate lid.

Black Vinyl Films for Fluorescence and  
Photoprotection

VWR 89087-692

Microplate Lids  
Polystyrene

Greiner Bio-One 07000288

60.5 Load  **10 µL** SG-II WS to each well in the microplate with 0.5 mL tip of stepper, and cover the plate with the black-film lid.

Finntip™ Stepper Pipette Tips  
500 µL

Thermo Scientific™ 9404170

96-Well Black Microplates  
Polystyrene

Greiner Bio-One 655076


Wipe or dab the liquid drop on the outside of the tip, avoid wiping the tip open before dispensing the liquid.

61 After incubation, vortex each tube for  **00:00:02** and then place into the fridge to stop the reaction.

62 Allow samples to reach  **Room temperature** for  **00:02:00** before loading the microplate.

63 Label the microplate.

64 Organize tubes in 96-well microtube rack in the same order as how microplates are loaded.

65 Load  190 µL working sample to the microplate by reverse pipetting.

45m

Blank must be included in each plate.

Pink area: RNA standard solutions for RNA standard curves



Blue area: DNA standard solutions for DNA standard curves

Yellow area: Samples and blanks

	1	2	3	4	5	6	7	8	9	10	11	12
A	R1A	R2A	R3A	R4A	R5A	R6A	1A	1A	1B	1B	1C	1C
B	R1B	R2B	R3B	R4B	R5B	R6B	2A	2A	2B	2B	2C	2C
C	R1C	R2C	R3C	R4C	R5C	R6C	3A	3A	3B	3B	3C	3C
D	D1A	D2A	D3A	D4A	D5A		4A	4A	4B	4B	4C	4C
E	D1B	D2B	D3B	D4B	D5B		5A	5A	5B	5B	5C	5C
F	D1C	D2C	D3C	D4C	D5C		6A	6A	6B	6B	6C	6C
G	9A	9A	9B	9B	9C	9C	7A	7A	7B	7B	7C	7C
H	Blank A	Blank A	Blank B	Blank B	Blank C	Blank C	8A	8A	8B	8B	8C	8C

	1	2	3	4	5	6	7	8	9	10	11	12
A	Blank A	Blank A	Blank B	Blank B	Blank C	Blank C	17A	17A	17B	17B	17C	17C
B	10A	10A	10B	10B	10C	10C	18A	18A	18B	18B	18C	18C
C	11A	11A	11B	11B	11C	11C	19A	19A	19B	19B	19C	19C
D	12A	12A	12B	12B	12C	12C	20A	20A	20B	20B	20C	20C
E	13A	13A	13B	13B	13C	13C	21A	21A	21B	21B	21C	21C
F	14A	14A	14B	14B	14C	1C	22A	22A	22B	22B	22C	22C
G	15A	15A	15B	15B	15C	15C	23A	23A	23B	23B	23C	23C
H	16A	16A	16B	16B	16C	16C	24A	24A	24B	24B	24C	24C

Wipe or dab the liquid drop on the outside of the tip, avoid wiping the tip open before dispensing the liquid.

66 Shake black film covered microplate at  Room temperature for  00:10:00

10m

67 Setup microplate reader:

Plate: Greiner F bottom chimney well PP 96 well;

Shake: Continuous 5s at 600 rpm

Endpoint reading: Ex 490 nm/Em 520 nm;

Varioskan LUX Multimode Microplate  
Reader

Thermo Fisher VL0L00D0

68 Read fluorescence and export data to excel sheet.



## 69 RNA standard curve

## 69.1 Concentrations of RNA standards in the microplate

Standard	2 ug/mL (uL)	Tris+WS (uL)	0.05% STEB (uL)	SG II WS (uL)	Final (ng/mL)
R1	0	700	250	50	0
R2	10	690	250	50	20
R3	25	675	250	50	50
R4	50	650	250	50	100
R5	100	600	250	50	200
R6	150	550	250	50	300

69.2 Slope of fluorescence in Tube A vs concentration of RNA standard gives  $m_{\text{RNA}+\text{DNase}}$  ( $\approx 0.03$ )

Slope of fluorescence in Tube B vs concentration of RNA standard gives  $m_{\text{RNA}+\text{RNase}}$

69.3 Calculate  $\rho$ 

$$\rho = \frac{m_{\text{RNA}+\text{RNase}}}{m_{\text{RNA}+\text{DNase}}}$$

## 70 Total RNA of the samples

$$\mu\text{g Total RNA (ml assay)}^{-1}$$

$$= 0.001 \times \frac{(RFU_A - RFU_{A\text{Blank}}) - (RFU_C - RFU_{C\text{Blank}})}{(1 - \rho)} / m_{\text{RNA}+\text{DNase}}$$

Where,

$RFU_A$  and  $RFU_C$  are the fluorescence in Tube A and Tube C of the same sample.

$RFU_{A\text{Blank}}$  and  $RFU_{C\text{Blank}}$  are the fluorescence in Tube A and Tube C of the blank.

## 71 DNA standard curve

## 71.1 Concentrations of DNA standards in the microplate

Standard	1 ug/mL (uL)	Tris+WS (uL)	0.05% STEB (uL)	SG II WS (uL)	Final (ng/mL)
D1	10	690	250	50	10
D2	20	680	250	50	20
D3	40	660	250	50	40
D4	70	630	250	50	70
D5	100	580	250	50	100

71.2 Slope of fluorescence in Tube A vs concentration of DNA standard gives  $m_{\text{DNA}+\text{DNase}}$

Slope of fluorescence in Tube B vs concentration of DNA standard gives  $m_{\text{DNA}+\text{RNase}}$  ( $\approx 0.12$ )

71.3 Calculate  $\delta$

$$\delta = \frac{m_{DNA+DNase}}{m_{DNA+RNase}}$$

72 Total DNA of the samples

$$\begin{aligned} & \mu g \text{ Total DNA (ml assay)}^{-1} \\ &= 0.001 \times \frac{(RFU_B - RFU_{BBlank}) - (RFU_C - RFU_{CBlank})}{(1 - \delta)} / m_{DNA+RNase} \end{aligned}$$

Where,

$RFU_B$  and  $RFU_C$  are the fluorescence in Tube B and Tube C of the same sample

$RFU_{BBlank}$  and  $RFU_{CBlank}$  are the fluorescence in Tube B and Tube C of the blank.

73 Dilution factor=40

If,

- Sample is extracted by 1 mL extraction reagent
- In Set 1, sample is diluted to 1/10
- In Set 3, diluted by Tris and all working solutions to 250/950
- In microplate, diluted by SG-II WS to 190/200