



Sep 14, 2022

Midbrain organoid generation from mfNPC

Rachel Bates¹

¹UCL

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.6qpvr4x1pgmk/v1](https://doi.org/10.17504/protocols.io.6qpvr4x1pgmk/v1)

[rachel.bates](#)

ABSTRACT

This protocol describes our method for the differentiation of human floor plate neural progenitor cells into human midbrain-like organoids (hMLOs). This protocol has been developed using a combination of several published protocols.

Adapted from

Jo J, Xiao Y, Sun AX, Cukuroglu E, Tran HD, Göke J, Tan ZY, Saw TY, Tan CP, Lokman H, Lee Y, Kim D, Ko HS, Kim SO, Park JH, Cho NJ, Hyde TM, Kleinman JE, Shin JH, Weinberger DR, Tan EK, Je HS, Ng HH (2016). Midbrain-like Organoids from Human Pluripotent Stem Cells Contain Functional Dopaminergic and Neuromelanin-Producing Neurons.. Cell stem cell.
<https://doi.org/10.1016/j.stem.2016.07.005>

Mohamed NV, Sirois J, Ramamurthy J, Mathur M, Lépine P, Deneault E, Maussion G, Nicoulet M, Chen CX, Abdian N, Soubannier V, Cai E, Nami H, Thomas RA, Wen D, Tabatabaei M, Beitel LK, Singh Dolt K, Karamchandani J, Stratton JA, Kunath T, Fon EA, Durcan TM (2021). Midbrain organoids with an SNCA gene triplication model key features of synucleinopathy.. Brain communications.
<https://doi.org/10.1093/braincomms/fcab223>

DOI

[dx.doi.org/10.17504/protocols.io.6qpvr4x1pgmk/v1](https://doi.org/10.17504/protocols.io.6qpvr4x1pgmk/v1)

PROTOCOL CITATION

Rachel Bates 2022. Midbrain organoid generation from mfNPC. **protocols.io**
<https://protocols.io/view/midbrain-organoid-generation-from-mfnpc-cgastsee>

FUNDERS ACKNOWLEDGEMENT

ASAP
Grant ID: ASAP2022

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 07, 2022

LAST MODIFIED

Sep 14, 2022

PROTOCOL INTEGER ID

69682

Day 0

- 1 Human floor plate neuronal progenitor cells (mfNPC) were derived using Smits 2019 protocol. They were maintained for a minimum of 5 passages before being used to generate organoids in maintenance medium.

Fedele S, Collo G, Behr K, Bischofberger J, Müller S, Kunath T, Christensen K, Gündner AL, Graf M, Jagasia R, Taylor V (2017). Expansion of human midbrain floor plate progenitors from induced pluripotent stem cells increases dopaminergic neuron differentiation potential.. Scientific reports.

<https://doi.org/10.1038/s41598-017-05633-1>

Smits LM, Reinhardt L, Reinhardt P, Glatza M, Monzel AS, Stanslowsky N, Rosato-Siri MD, Zanon A, Antony PM, Bellmann J, Nicklas SM, Hemmer K, Qing X, Berger E, Kalmbach N, Ehrlich M, Bolognin S, Hicks AA, Wegner F, Sternecker JL, Schwamborn JC (2019). Modeling Parkinson's disease in midbrain-like organoids.. NPJ Parkinson's disease.

<https://doi.org/10.1038/s41531-019-0078-4>

1.1 mfNPC maintenance medium

☒ DMEM/F-12 Thermo

[M] 50 % volume

Fisher Catalog #11320033

[M] 50 % volume

☒ Neurobasal Plus Medium Gibco,

ThermoFisher Catalog #A3582901

1:50

☒ N21-MAX Media Supplement (50X) R&D

Systems Catalog #AR008

☒ N-2 max supplement R&D

1:100

Systems Catalog #AR009

☒ Glutamax (100x) Gibco - Thermo

[M] 1 % volume

Fischer Catalog #35050-061

☒ MEM, NEAA, no glutamine Thermo

[M] 1 % volume

Fisher Catalog #10370021

☒ SB431542 Cell Signaling

[M] 10 micromolar (µM)

Technology Catalog #14775

☒ SAG 1 mg Stemcell

[M]0.5 micromolar (µM) Technologies Catalog #73412

[M]250 millimolar (mM)

☒ Stemolecule LDN-193189 Stemgent - Bio-

connect Catalog #04-0074

☒ L-Ascorbic acid Sigma

[M]200 micromolar (µM) Aldrich Catalog #A4403

☒ CHIR99021 R&D

[M]3 micromolar (µM) Systems Catalog #4423

2 mfNPCs were detached using accutase at 37 °C for 00:03:00 .

3m

3 Re-suspend cells in d0 induction medium and plate 9,000 cell/well in ultra-low attachment U-bottomed 96 well plates.

BIOFLOAT plate from Facellitate work best for us compared to Corning or Nunc at producing uniform EBs.

Day 2

4 Change the medium to mfNPC medium supplemented with [M]0.0001 mg/mL

☒ FGF8 Novus

Biologicals Catalog #423-F8-025

[M]0.0001 mg/mL

☒ recombinant human SHH protein Qkine being careful not to touch the organoid.

Day 4

5 Change the medium to patterning I medium adding 300 µL per well.

5.1 patterning I medium

☒ DMEM/F-12 Thermo

[M]50 % volume

Fisher Catalog #11320033

[M]50 % volume

☒ Neurobasal Plus Medium Gibco,

ThermoFisher Catalog #A3582901

1:50

☒ N21-MAX Media Supplement (50X) R&D

Systems Catalog #AR008

☒ N-2 max supplement R&D

1:100

Systems Catalog #AR009

	Glutamax (100x) Gibco - Thermo
[M] 1 % volume	Fischer Catalog #35050-061
	MEM, NEAA, no glutamine Thermo
[M] 1 % volume	Fisher Catalog #10370021
	SAG 1 mg Stemcell
[M] 0.5 micromolar (μM)	Technologies Catalog #73412
	L-Ascorbic acid Sigma
[M] 200 micromolar (μM)	Aldrich Catalog #A4403
	CHIR99021 R&D
[M] 3 micromolar (μM)	Systems Catalog #4423
	FGF8 Novus
[M] 0.0001 mg/mL	Biologicals Catalog #423-F8-025
[M] 0.0001 mg/mL	recombinant human SHH protein Qkine

Day 6

- Change medium to patterning II medium with reduced CHIR adding [300 μL](#) per well.

6.1 Patterning II medium

	DMEM/F-12 Thermo
[M] 50 % volume	Fisher Catalog #11320033
[M] 50 % volume	
	Neurobasal Plus Medium Gibco,
	ThermoFisher Catalog #A3582901
1:50	
	N21-MAX Media Supplement (50X) R&D
	Systems Catalog #AR008
	N-2 max supplement R&D
1:100	Systems Catalog #AR009
	Glutamax (100x) Gibco - Thermo
[M] 1 % volume	Fischer Catalog #35050-061
	MEM, NEAA, no glutamine Thermo
[M] 1 % volume	Fisher Catalog #10370021
	L-Ascorbic acid Sigma
[M] 200 micromolar (μM)	Aldrich Catalog #A4403
	CHIR99021 R&D
[M] 0.7 micromolar (μM)	Systems Catalog #4423
	FGF8 Novus
[M] 0.0001 mg/mL	Biologicals Catalog #423-F8-025
[M] 0.0001 mg/mL	recombinant human SHH protein Qkine

- [Cultrex Reduced Growth Factor Basement Membrane Extract R&D](#)

Thaw [Systems Catalog #343301001](#)

[Overnight](#) at [4 °C](#) ready for next day.

Day 8 30m

- 8 Carefully remove as much medium as possible from each well being careful not to touch the organoid.
- 9 Add **15 μ L** of **Cultrex Reduced Growth Factor Basement Membrane Extract R&D Systems Catalog #343301001** to each well and return to **37 °C** incubator for **00:30:00**.
- 10 Carefully add **300 μ L** of tissue induction medium to each well.

Some organoids may float others with remain attached to plate, this does not affect the organoid.

10.1 Tissue induction medium

100 % volume

Neurobasal Plus Medium Gibco, ThermoFisher Catalog #A3582901

1:50

N21-MAX Media Supplement (50X) R&D Systems Catalog #AR008

1:100

N-2 max supplement R&D Systems Catalog #AR009

1 % volume

Glutamax (100x) Gibco - Thermo Fischer Catalog #35050-061

1 % volume

MEM, NEAA, no glutamine Thermo Fisher Catalog #10370021

200 micromolar (μ M)

L-Ascorbic acid Sigma

Aldrich Catalog #A4403

0.0001 mg/mL

FGF8 Novus

Biologicals Catalog #423-F8-025

0.0001 mg/mL

recombinant human SHH protein Qkine

0.00025 mg/mL

Laminin from Engelbreth-Holm-Swarm murine sarcoma basement membrane Sigma – Aldrich Catalog #L2020





Insulin Sigma

0.025 mg/mL Aldrich Catalog #I2643-25MG

Day 9 30m

- 11 Using sterilised scissors or a scalpel blade cut a pasteur pipette at the widest point.

this will provide a sufficient bore to pick up the embedded organoid without damaging the cultrex.


- 12 Either add  **2 mL** of Tissue induction medium to [Ultra-low attachment 6-well plate Corning Catalog #CLS3471](#) or add  **500 µL** to Ultra-low attachment 48 well plate.
- 13 using the cut pasteur pipette carefully pick up each organoid and move to new plate. For 6 well plate add up to 6 organoids per well or 1 per well in the 48-well plate.
- 14 Incubate plates at  **37 °C** on a orbital shaker set to  **70 rpm**.

Day 10

- 15 Change medium to Differentiation medium.

15.1 Differentiation medium


[M] 100 % volume


 [Neurobasal Plus Medium Gibco,](#)
ThermoFisher Catalog #A3582901

1:50

 [N21-MAX Media Supplement \(50X\) R&D](#)
Systems Catalog #AR008

1:100

 [N-2 max supplement R&D](#)
Systems Catalog #AR009

 [Glutamax \(100x\) Gibco - Thermo](#)

[M] 1 % volume

Fischer Catalog #35050-061

 [MEM, NEAA, no glutamine Thermo](#)

[M] 1 % volume

Fisher Catalog #10370021

 [L-Ascorbic acid Sigma](#)

[M] 200 micromolar (µM)

Aldrich Catalog #A4403

 [Db-cAMP \(dibutyl-cyclic](#)

[M] 125 micromolar (µM)

AMP) Sigma Catalog #D0627

10ng/mL  [BDNF Qkine Catalog #QK050](#)

10ng/mL  [GDNF Qkine Catalog #QK051](#)

- 16 After day 10 perform a 75% medium change every 2-3 days.

17